



# **Version 4**

## **User Manual**

HiTrack was conceptualized and designed with invaluable assistance from a wide variety of dedicated professionals, including hospital-based newborn hearing screeners, audiologists, and Early Hearing, Detection and Intervention program coordinators. Their knowledge and experience is evident throughout the software and is acknowledged with respect and appreciation.

™ The NCHAM logo and HiTrack are trademarks of NCHAM. ®Natus ALGO 2 is a registered trademark of Natus Medical, Inc. ®Bio-logic and AuDX are registered trademarks of Bio-logic Systems Corp. ™ABaer and HATS are trademarks of Bio-logic Systems Corp. All other marks belong to their respective owners.

For support: Phone our Help Desk at 435.797.3584 or E-mail us at [ncham.helpdesk@usu.edu](mailto:ncham.helpdesk@usu.edu)

National Center for Hearing Assessment & Management (NCHAM)

Utah State University 2615 Old Main Hill -Logan, UT 84322 Tel: 435.797.3584

[www.hitrack.org](http://www.hitrack.org)

Revised 6/2018

**Table of Contents**

..... 1

Version 4..... 1

User Manual..... 1

    Table of Contents..... 3

        Introduction..... 11

        Hearing Screening Protocols..... 11

        Data Management..... 11

        State or Regional Reporting..... 12

        This User Manual..... 12

Opening and Configuring HiTrack..... 13

User Interface Layout..... 14

Settings and Administration - System Settings..... 15

    Users..... 15

    User Groups..... 17

    Access Rights..... 17

    Facilities..... 19

    Merge Sources..... 20

    Facility Contact Info..... 20

    Facility Groups..... 21

    Nurseries & Pick Lists..... 21

    Advanced Settings..... 21

    Centralized HiTrack Setting..... 21

Facility Specific Settings..... 23

    Medical People..... 23

    Letter Settings..... 23

    Auto-Fills..... 24

    Milestones..... 24

    Other Settings..... 24

User Specific Settings..... 24

    Automatic Email..... 24

    Password..... 25

Sharing Program Information with Screening Software.....	26
Tracking Features Overview.....	27
The Filing Cabinet .....	27
Action Needed Drawer .....	27
No Action Drawer.....	29
All Folders .....	30
Baby List.....	30
Sorting.....	31
Filters .....	31
Diagnostic Result Filters.....	31
Returning to the List .....	32
Facility Group Baby List .....	32
Combining Duplicate Records.....	32
The Baby Basket.....	33
Collecting Screening Information using HiTrack.....	35
Facility Protocol Setting .....	35
Other Settings.....	35
Procedure Overview: .....	35
Step 1: Creating baby records .....	36
Step 2: Printing a list of infants to be screened.....	40
Step 3: Entering screening results and finishing screening records.....	41
Step 4: Entering Outpatient screening results in a Two-Stage Protocol.....	42
Collecting Screening Information using Equipment Software .....	43
Facility Protocol Setting .....	43
Facility Merge Sources .....	43
Matching Rules .....	43
Procedure Overview: .....	43
Step 1: Collecting Screening Information using Equipment Software.....	44
Step 2: Starting the Merge process .....	44
Step 3: Resolving any pending Merge batches .....	46
Collecting Diagnostic Hearing Information.....	48
Set up and IT Configuration .....	48
Matching Rules and Required Fields.....	48

Audiologist Onboarding Process .....	49
Suggested Audiologist Welcome Message .....	50
Configuration .....	50
Managing Screening and Tracking Data.....	51
Toolbar .....	51
Tracking Tasks.....	52
Procedure Overview: .....	52
Tracking Tasks - Step 1: Updating baby records.....	53
Tracking Tasks - Step 2: Generating Letters .....	73
Tracking Tasks - Step 3: Generating Reports .....	75
General Filtering Options .....	75
Available Reports.....	76
Report Advanced Criteria .....	78
Printing and Exporting Reports .....	78
Milestones .....	79
Tracking Tasks - Step 4: Transferring Data to a Central Database .....	80
Maintaining the HiTrack Database .....	84
Program Information.....	84
Change Log and Audits .....	84
State/Regional Transfer .....	84
Back-ups and Conversion .....	85
Key Maintenance .....	86
Advanced Tools .....	86
Classification Browser Tool .....	86
Data Integration .....	88
Helpful Hints for Data Entry in HiTrack .....	89
Navigating screens and fields .....	89
Opening Records.....	90
Entering Dates.....	90
Using drop-down menus .....	90
Entering File Locations.....	90
Data Grids .....	90
HiTrack's Web Companion.....	91

Web Browser Settings.....	91
Differences in Web Companion.....	92
Manually Uploading Equipment Files to the Web Companion.....	93
Automatic Uploading Equipment Files using the Merge Companion.....	94
Appendix A: Linking HiTrack with.....	97
Equipment Manufacturer Software.....	97
Manufacturer Software Instructions .....	97
About File Locations.....	97
File Processing.....	97
Program Information Sharing .....	97
Free Text Program Information Support.....	98
Automating the HiTrack Merge.....	98
Bio-logic System Corp ABaer and AOAE .....	99
(Version 3.17.0.i.020).....	99
Link Setup.....	99
Data Entry .....	100
Export tips.....	102
Natus ALGO 2e Color or 3 with DataBook.....	103
Link Setup.....	103
Data entry and export tips .....	104
Natus ALGO 2 or 2e with DataBook.....	105
Link Setup.....	105
Data entry and export tips .....	105
Otodynamics EZ-Screen2 .....	108
(Version 6.25.0.0).....	108
Link Setup.....	108
Data Entry .....	109
Export tips.....	111
Appendix B: Report Criteria Documentation .....	112
Classifications .....	112
Inpatient Classifications.....	112
Outpatient Classifications .....	113
Evaluation Classifications.....	114

Intervention Classifications .....	115
Amplification Report .....	116
Report Specific Filtering Options .....	116
JCIH 2007 Benchmark.....	116
CDC Survey Report .....	117
Report Specific Filtering Options .....	117
Record Browse List Feature.....	117
Lost to Follow-Up % and Number (Diagnosis).....	117
Age Calculations and Types of Loss Considered .....	117
Non-Resident CDC Counts.....	117
Adoption Counts .....	118
Hearing Loss Type and Severity Summary .....	118
Hearing Status for CDC Reporting.....	119
Missed ears.....	120
Warning List Feature .....	120
User Defined Report Columns.....	121
CDC Survey Report Counts.....	122
Survey Counts.....	122
Early Intervention Counts .....	124
Additional Cases Not Reported .....	127
Demographics Report .....	127
Data Cleaning Reports .....	128
Duplicate Medical ID.....	128
Probable Duplicates .....	128
Custom Duplicates .....	128
Report Specific Filtering Options .....	128
Combining Duplicates .....	128
Early Intervention Report .....	128
Report Specific Filtering Options .....	128
Hearing Status Report.....	129
Report Specific Filtering Options .....	129
Needing Evaluation Report.....	129
Report Specific Filtering Options .....	129

Needing Outpatient Report .....	130
Report Specific Filtering Options .....	130
Outpatient Screening Report.....	130
Report Specific Filtering Options .....	130
Physicians Report.....	131
Report Specific Filtering Options .....	131
Risk Indicator Report.....	131
Report Specific Filtering Options .....	131
Screeener Profile Report.....	132
Report Specific Filtering Options .....	132
Screening Results Report.....	132
Report Specific Filtering Options .....	132
Flow Chart .....	133
Inpatient Section .....	133
Outpatient Section.....	133
Hearing Status .....	133
Confirmed Loss - Other Category.....	134
Evaluation in Process .....	134
Unmonitored Infants with Risk Indicators.....	134
Pass Infants monitored for Risk Indicators.....	134
Report Specific Filtering Options .....	135
Grouping Options.....	135
Hospital Summary .....	136
Report Specific Filtering Options .....	136
Grouping Options.....	136
Quarterly Report .....	137
Percent Screened as Inpatient .....	137
Percent Passing Inpatient Screening .....	137
No Info .....	137
Deceased .....	137
Percent Completing Outpatient Screening.....	137
Referred For Diagnostics .....	137
Other Dx.....	138

No Follow-Up .....	138
Report Specific Filtering Options .....	138
Grouping Options.....	138
Technology Summary.....	139
Report Specific Filtering Options .....	139
Grouping Options.....	139
Milestone Summary Report .....	140
"Screened by" Milestone: .....	140
"Hearing Identified by" Milestone:.....	140
"Enrolled by" Milestone: .....	140
Report Specific Filtering Options .....	141
Grouping Options.....	141
User Defined Reports.....	142
Advanced Criteria .....	143
Adding criteria .....	143
Modifying criteria .....	143
Boolean "OR" Operation .....	143
Appendix C: Classification Combinations.....	144
Appendix D: Merge Matching Criteria and Overwrite rules .....	150
Step 1: Filters are used to find a list of possible matches .....	150
Equipment Merges .....	150
HiTrack 3.5 ETF Merges .....	150
HiTrack 4 ETF4 Merges .....	150
Step 2: Calculating match score for possible matches.....	151
Baby Field Matching Scores.....	151
Other Scores.....	151
Match Score Examples .....	152
Field Overwriting.....	153
Overwriting Notes.....	153
Appendix E: Most Conclusive Result Calculation .....	154
Overview.....	154
Most Conclusive Result Ranking Rules .....	154
Most Conclusive Result Window .....	155

Example #1 .....	155
Example #2 .....	155
“Same Day” Most Conclusive Result Option.....	155
Appendix F: Automatic Letter Generation Criteria.....	156
Criteria .....	156

## Introduction

One of the greatest challenges for any Early Hearing Detection and Intervention (EHDI) program is to ensure that all infants are screened for hearing loss and that those who do not pass receive appropriate and timely follow-up care. In response to hospital, regional, and national needs, NCHAM developed the HiTrack data management system which is now being used in hundreds of hospitals across the United States. Program coordinators operating or monitoring hospital-based screening programs have discovered that a computer-based patient information system contributes significantly to program success. HiTrack can support any screening program using otoacoustic emissions (OAEs) or automated auditory brainstem responses (ABRs) to screen infants for hearing loss.

## Hearing Screening Protocols

Most programs use either a single or two-stage protocol to conduct their hearing screening.

A Two Stage screening protocol is designed to screen the baby with cheap and efficient technology (for example, OAE) while the baby is at the hospital. This stage of screening is called Inpatient screening because the mother has not left the hospital. If the baby is not screened or the outcome of the Inpatient screening is not satisfactory the hospital then makes appointments to have the mom and baby come back for Outpatient screening with a more advanced technology (Automated ABR).



Single Stage screening protocols are also becoming more popular as Automated ABR screening becomes more cost effective. In this protocol, infants who do not pass Inpatient screening are referred directly for diagnostic evaluation.

## Data Management

The HiTrack 4 Series represents the cutting-edge of EHDI Data Management. The HiTrack 4 database is built around a dynamic framework to provide the logic needed for managing the complex requirements of EHDI tracking and follow-up. The HiTrack client offers a task-oriented user interface and is designed for maximum flexibility.

The Windows user interface offers all of the HiTrack functionality in a quick, responsive Windows application that requires minimal configuration. The Windows Client allows HiTrack users to run their own individual database, or connect to a remotely hosted database to access shared data.

The [Web Companion](#) user interface offers virtually all of the same interface and features of the Windows user interface while providing the versatility of being able to securely access HiTrack from any computer connected to the internet.

## State or Regional Reporting

If you will be reporting your data to a central EHDI coordinator please contact them for instructions on how you will be using HiTrack: There are two common Database Layouts that can be used with HiTrack:

<p>Distributed</p> 	<p>In this layout each facility in your region or state has their own fully functional database dedicated to their own HiTrack data.</p> <p>A separate Central database is maintained at the state or region. The EHDI coordinator maintaining this database will populate their database by collecting Data Transfer files from each facility in their area of responsibility. They will typically set a schedule for how often (e.g. once a month) the facilities will sent their files.</p> <p>Statewide or regional <a href="#">Program Information</a> lists, such as physicians, facilities and audiologists are often provided for import into the facility databases. This ensures that records in the Central database contain unique and synchronized Program Information entries.</p> <p>This is referred to as a Distributed layout because the Central database is distributed across all of the separate facility databases.</p>
<p>Centralized</p> 	<p>In this layout a single Central database is maintained at the state or region and it is shared with each of the facilities in the area of responsibility. This sharing is typically done with HiTrack's <a href="#">Web Companion</a>.</p> <p>HiTrack provides an <a href="#">advanced setting</a> to enable data management features for a Centralized layout. Please review this setting with the <a href="#">HiTrack Help Desk</a>.</p>

Some states or regions may decide to combine these layouts. For example the Central coordinator may collect Data Transfer files from each hospital and also share their database with clinic offices using the Web Companion.

There are a few additional tasks required if your state or region has arranged for you and your Facility to maintain your own database in a Distributed database layout.

These tasks will be noted in this manual with a "Distributed" icon.



If you do not report your data to a central EHDI coordinator or if you use HiTrack's Web Companion the Data Transfer tasks can be omitted.

## This User Manual

This hyperlinked User Manual will help you configure HiTrack and will help you to become familiar with its functions. Clicking on the [blue underlined links](#) inside this document will scroll the page down to related sections.



This manual does not cover software installation. Please refer to the HiTrack 4 Installation Guide for system requirements and installation instructions.

## Opening and Configuring HiTrack


To configure HiTrack you will need to log in. During the installation process you will be asked to create a username and password to use in HiTrack. If you are starting a new database you were also asked to enter in your Facility name. If you already have an existing HiTrack 3.5 database you should contact the [HiTrack Help Desk](#) for conversion steps.

When you are ready to configure HiTrack, log in using the following steps:

1. Enter the Username and Password (passwords are case sensitive)



2. Click on the Start button or press the Enter key and the following box will appear:



3. Click on the drop-down box and select a Facility or Facility Group. Names of Facility Groups are preceded by an asterisk.
4. Click on the Continue button or press the Enter key.

## User Interface Layout

After successfully logging in the HiTrack user interface opens. The user interface is made up of the following sections:

The screenshot shows the HiTrack user interface for 'All Data for Green River Regional (15)'. The interface is divided into several sections:

- File Cabinet & File Card:** Located on the left side, it displays patient information for 'Hanson, James', including Medical ID (1457896), Birth Date (10/8/2012), and Action (Need Inpatient Screening). It also shows a 'Record Overview' with various service categories like Tracking Options, Demographics, Screening, Risk Indicators, and Diagnostics.
- Facility Selector:** A dropdown menu at the top center showing 'Green River Regional'.
- Tool Bar:** A horizontal bar at the top right containing icons for Letters, Reports, Merge, Tools, Settings, and Log Out.
- Action Buttons:** A row of buttons at the bottom including 'Add Baby', 'Delete Baby', and 'Print List'.
- Work Area:** The central area containing a table of patient data and a 'Quick Add' form for 'Hanson, James'.

The table in the Work Area displays the following data:

Baby Name	Medical ID	Birth Date	Inpatient	IP Date	Outpatient	OP Date
Hanson, James	1457896	10/8/2012			No Out. Scrn->Eval	10/9/2012
Mortonson, Emily	9568421	10/8/2012				
Douglas, Sally	8974652	10/8/2012	Referred	10/8/2012	Referred	10/8/2012
Hatch, Felicia	8723691	10/6/2012	Inconclusive	10/8/2012	No Out. Scrn->Eval	10/8/2012 8:
Rose, Kevin	8971246	10/5/2012	Referred	10/8/2012	Inconclusive	10/8/2012
Stevens, Rebecca	49876324	10/3/2012				
Willard, Kevin	7984632	10/1/2012			No Out. Scrn->Eval	10/8/2012
Smith, James	9451278	10/1/2012	Referred	10/8/2012		
Douglas, Cammie	8974521	10/7/2012	Passed	10/8/2012		
Fredricks, Cloey	8974563	10/7/2012	Passed	10/8/2012		
Morgonson, Alice	8679124	10/6/2012	Passed	10/8/2012		
Rose, Brittny	5486128	10/5/2012	Inconclusive	10/8/2012	Passed	10/8/2012

The 'Quick Add' form for 'Hanson, James' includes the following fields:

- Stage: Inpatient
- Type: TEOAE
- Right: -Blank-
- Left: -Blank-
- Date: 10/ 9/2012
- Time: 00 : 00
- Source: -Blank-
- Facility: Green River Regional

## Settings and Administration - System Settings

On a new installation there are specific items that need to be setup before you can begin entering baby data into the system. Click on the Settings dial icon in the toolbar in the upper right corner of the screen. This opens the first Settings screen shown below.

Username	Group	Expiration date	Active
demo	System Administrators		Y
Dora	System Administrators	08/29/2013	Y
Kevin	Screeners		Y
Michelle	Facility Administrators		Y

Settings are grouped into three categories listed on the left, System-wide settings, Facility specific settings and User specific settings.

### Users

The Users list is the first system setting and allows you to maintain a list of users, their passwords, expiration dates and facility assignments.

To add a new user login:

1. Click on the "Add" button below the list.

The Add User screen is shown

The screenshot shows the 'Add User' form with the following fields and options:

- Profile:**
  - Username: [text input]
  - Password: [text input]
  - Confirm Password: [text input]
  - Password Strength: Weak (with a help icon)
  - Account Expires: [checkbox] [date dropdown]
  - Account Active: [checked checkbox]
  - User Group: [Data Entry Staff] (dropdown)
  - Facility Group: [no group access] (dropdown)
- Facility Access:**
  - Alton Mem Hospital [checkbox]
  - Creekside Hospital [checkbox]
  - Green River Medical [checkbox]
  - JJA-Alton Mem Hospital [checkbox]
  - Logan Regional Hospital [checkbox]
  - Mountain View Medical Center [checkbox]
  - UNKNOWN [checkbox]
  - Z HT35 Conv. [checkbox]
- Demographics:**
  - Name-Last: [text input]
  - First: [text input]
  - Middle: [text input]
  - Address: [text input]
  - Phone: [text input]
  - Ph 2: [text input]
  - Email: [text input]
  - City/State/Zip: [text input]
  - County: [text input]
  - Birth Date: [checkbox] [date dropdown]
  - Language: [blank] (dropdown)
  - Title/Suffix: [text input]
  - Race: [blank] (dropdown)
  - Gender: [-Blank-] (dropdown)
  - Education: [blank] (dropdown)
  - Ethnicity: [blank] (dropdown)
  - Deceased: [checkbox]

Buttons: Find, Save, Cancel. Checkbox: Require User to Change Password at Next Login.

2. Enter in a user name and password and select an expiration date in the future using the dropdown to the right of the date field.

Select a User Group that provides an appropriate level of access for this person. (described in the next section)

Using the checkboxes in the Facility Access list identify the Facility to be used by this person. Select more than one Facility if this person will be signing into this HiTrack system from multiple facilities.

HiTrack requires strong passwords and will automatically indicate when strong password is entered. Click the "?" for information about strong password criteria.

3. Use the "Find" button to locate a person record already in HiTrack (i.e. a Screener, Doctor or Audiologist.) or enter in values for the appropriate demographic fields.


Colored fields are mandatory.

4. Click the "Save" button below the demographic fields to add this user.

## User Groups

The User Groups section allows you to assign access rights to the users within a group. Access rights allow administrators to control each user's ability to make changes to data and grant access to specific screens.

The built-in Administrators group allows all access rights. A few other stock User Groups have been built-in to provide a starting point for program coordinators that are new to HiTrack 4. It is best to create your own User Groups rather than modifying stock HiTrack User Groups.



## Access Rights

<i>Access Right</i>	<i>Explanation</i>
System Administrator	Allows access to System-wide Settings.
Facility Administrator	Allows access to Facility related Settings and features. Including: User Maintenance (for Facilities already assigned to this User), Merge Sources, Letter Settings, Required Fields and other settings.   This Access Right is especially useful in <a href="#">Centralized Layouts</a> because it allows local Facility staff to manage everything related to their own settings, including new

	user accounts.
Edit Existing Records	This access right grants the ability to modify existing data in HiTrack.
Add New Records	This access right grants the ability to add additional data to HiTrack.
Baby	Allows access to view the following baby related areas: Demographics, Tracking History, Risk Indicators, Letters, Recommended Actions, Transfer History, EI Services and Record Overview.
Diagnostics	Allows access to Diagnostics, Hearing Disposition, and Amplification areas. Requires -Baby- Access Right.
Letters	Allows access to Letter features. Requires -Baby- Access Right to work with individual baby letters.
Medical People Settings	Allows access to Medical People from the Settings Screen.
Merge Data	Allows basic access to Merge tool. Users can merge new data into the system. Requires –Add New Records- Access Right.
Resolve Merge Data	Allows access to features used to make corrections to incoming Merge data.
Clear Merge Data	Allows access to features used to delete pending Merge data.
Record Combination	Allows advanced users to permanently combine baby records from the Baby List.
Reports	Allows access to Reports.
Screening	Allows access to view the Screening area. Requires -Baby- Access Right.
Tools	Allows access to basic Tools, including Key Maintenance, Change Logs, Exporting Program Information, the Database Manager and Conversion tools.
State / Regional Transfer	Allows access to Tools that create data transfer files for State or Regional data collection agencies.
Advanced Tools	Allows access to advanced features in the Tools section.
External Links	Allows access to external database links on the patient file folder.
Lab Testing	Allows access to Lab Testing data on the patient file folder.

## Facilities

You may add or edit Facilities in HiTrack from the Facilities category. Clicking the Add button at the bottom of the Facilities tab will show the following screen.



If you [will be reporting your HiTrack data](#) to a Central coordinator you must contact them and decide on a mutually agreeable way of representing your Facility name. The spelling of your Facility name must match what your Central coordinator is using.

This screen also requires you to set the [EHDI Protocol](#) for your Facility. Based on the screening protocol baby records are automatically filed into folders to help you track their hearing screening and evaluation.

Listed below are the available protocols and typical usage information.

<i>Protocol</i>	<i>Use if...</i>
Two-Stage Automated Pass/Manual Non-Pass	Use this Protocol if you want to manually review Inpatient and Outpatient results for infants who have not Passed.
Two-Stage Automated	Use this Protocol if you merge completed Inpatient and Outpatient screening results from screening equipment or other HiTrack systems.
Two-Stage Manual	Use this Protocol if you want to manually review Inpatient and Outpatient results for all infants.
Two-Stage Inpatient Automated	Use this Protocol if you merge completed Inpatient screening results from screening equipment or other HiTrack systems, but want to review Outpatient screening results by hand.
Two-Stage Dual Equipment Auto-Pass Only	Use this Protocol if you want to require A-ABR tests for infants who did not pass OAE screening and keep an active list of infants who have not Passed. (Infants who do not Pass Outpatient must

	be finished by hand.)
Two-Stage Dual Equipment Automated.	Use this Protocol if you merge completed OAE and A-ABR Inpatient screening results from screening equipment. (Infants in Outpatient screening must be finished by hand.)
Single-Stage Dual Equipment* Automated	Use this Protocol if you merge completed OAE and A-ABR Inpatient screening results from third-party screening equipment. Infants who do not pass screening will be recommended for evaluation.
Single-Stage Dual Equipment Auto-Pass Only*	Use this Protocol if you want to require A-ABR tests for infants who did not pass OAE screening and keep an active list of infants who have not Passed. Infants who do not pass A-ABR screening will be recommended for evaluation.

## Merge Sources

After saving changes to a facility click on the "Merge Sources" button to define the setup information needed to merge data into the HiTrack program. Merge sources included data from screening equipment software, hospital patient information system exports and data from other HiTrack databases.

Each merge source is entered by selecting a type of data file, specifying a name and browsing for a folder on your computer. The folder you select here is where HiTrack will look to [find the file](#) to be merged. Because some facilities use one device in the Well Baby nursery and one device in the NICU, HiTrack allows you to enter additional merge sources for each device.


Name/Type	Source Path
HITRACK ETF4 File - (HT4 Transfer File (ETF4))	C:\Merge Files\
Hosp. Admissions File - (Demographic Data (Expanded)	H:\Fileshare\
Physician List - (HT3.5 PI pednames.txt (pednames.txt))	H:\Fileshare\exports\
XML - (HT4 XML Data File)	W:\

## Facility Contact Info

The Contact Info button at the bottom of the Facilities tab is used to enter contact information for parent and physician letters. This information, including hospital Phone Number, should be entered prior to using HiTrack to generate letters.

## Facility Groups

From the Facility Group list, you may create groupings of facilities. This allows Central coordinators to easily manage and report on data from more than one place. Facility Groups also allow Central coordinators generate Letters as a state or regional entity. Contact Information should be entered for your Facility Group if you will be generating Letters as a Facility Group.



The Facility Group screen also allows you to assign access rights for users who need access to more than one facility in one easy step.

## Nurseries & Pick Lists

The Nurseries, Nursery Groups and Pick Lists settings allow you to customize the data entry screens for your facility. Under Pick Lists you can customize the following categories: Race, Insurance Type, Languages and Custom Fields. For example, if you added "Free Healthcare" as an insurance type, users could select it from a list of insurance types when entering information. The Code field allows you to synchronize data coming from screening equipment data files.

The Custom Fields feature allows you to setup new pick list categories and specify allowable entries. For example a Custom field category for "Eye Color" could be added with possible entries such as Blue, Green, and Brown. After being added this user defined custom field could then be used later in baby records and on Reports as Advanced Filter Criteria. Free-form text entry Custom fields are also allowed.



## Advanced Settings

The last category listed in the System settings is Advanced Settings. The Advanced Settings screen provides access to four setting categories. In most cases these settings will not need to be reconfigured, however experienced program coordinators may want to use these settings to customize HiTrack.

Features in Advanced Settings determine how Most Conclusive Results are calculated. Please refer to [Appendix E for more information about Most Conclusive Results](#) or consult the [HiTrack Help Desk](#).

## Centralized HiTrack Setting

As described above, HiTrack can be used in a [Centralized Database Layout](#). Use the Centralized Database checkbox in Advanced Settings to turn on this feature.

The Centralized setting changes how records are managed when baby transfers are entered. Newborns that become sick or need special care will often be transferred to a different facility

better suited to handle their conditions. (Note this type of Transfer is **not** the same as the [Data Transfer](#) feature.) In these cases the receiving facility may need access to incoming baby records so that they can conduct hearing screening.

A second Advanced Setting called “Sending Facilities retain access to Transferred out records” also determines how transferred baby records can be accessed. When enabled, sending facilities keep their access to the baby record in addition to sharing access with the receiving Facility. Although both facilities report on the record the record shows singly when accessed from a Facility Group.

The following table summarizes how records are handled by the Centralized setting:

	<b>Centralized</b>		<b>Distributed (Centralized Off)</b>
<b>Feature:</b>	<b>In System Transfer</b>	<b>Out of System Transfer</b>	<b>In system and Out of System behave the same.</b>
Baby's Tracking Status	"In Progress"	"Transferred Out"	"Transferred Out"
Baby shown in Sending Facility Baby List and Reports?	Yes if “Sending Facilities retain access...” is enabled. Otherwise No.	Yes	Yes
Baby shown in Receiving Facility Baby List and Reports?	Yes	No (N/A)	No
Baby shown in Sending Facility's "Transferred Out" Folder?	Yes, but with view only access, unless “Sending Facilities retain access...” is enabled.  The Receiving Facility name is available from the baby's Transfer History.	Yes	Yes

## Facility Specific Settings

The setting categories listed below your Facility or Facility Group name allow you to customize features for your Facility or Facility Group. Access to this area in HiTrack can be granted to groups of users with the [Facility Administrator Access Right](#).

The "Save" button at the bottom of each setting screen will save your changes. Some changes, especially User access related settings, will not take effect until you leave HiTrack and log back in.

## Medical People

The Medical People list allows you to add physicians, audiologists, screeners, and Early Intervention personnel for a specific Facility or Facility Group. Each data entry screen allows for fully detailed contact information.

If you will be sending letters to Physicians and Early Intervention contacts you should enter detailed address information so it can be included on the letters.

The Facility Assignment list at the bottom of the Medical Person data entry screen allows you to decide which facilities will see this person on their data entry lists. When adding a medical person to a Facility Group they are automatically assigned to all Facilities within the Group.

Older non-current personnel can be suppressed from new data entry screens by using a "Show?" option checkbox. This allows you to keep your historical data intact while maintaining streamlined data entry for new records.

The "Review Usage" feature (controlled by the Advanced Tools Access Right) allows you to delete unused entries or combine and replace duplicates by first selecting a Medical Person from the list.

## Letter Settings

Under the Letter Settings tab you can select which letters should be generated for the currently selected facility and indicate if you will be using HiTrack's built in RTF Letter formatting or Microsoft Word's Mail Merge feature to generate Letters.

Letter Templates can also be customized from this screen by clicking the Letter Templates button. After selecting the type of letter and the target language you can use the Rich Text editor or Microsoft Word to format the text, font and placement of patient data fields in the each letter. These fields contain the variable information that can be included in a letter such as a baby's name, contact name, mailing address, screening results, and so forth.

Letter Templates for Green River Regional

Letter Type: Inpatient Refer for Outp. Scn (PC) Language: ENGLISH

System Default Inpatient Refer for Outp. Scn (PC) Template

Georgia 10 [Rich Text Editor] Undo Import RTF Help Prev

<header>  
Dear Ms./Mr. <contact last name>,  
  
*Congratulations on the birth of your child!* As you may be aware, we provide hearing screening to hospital. The results for your baby are shown below:  
  
Inpatient Screen  
Right Ear: <right ear inpatient results>  
Left Ear: <left ear inpatient results>  
  
The results indicate that the hearing screening should be repeated. This does not necessarily mean loss. The most common reason for needing this second test is that some babies still have birth debris blocks sound from entering the ear at the time that the initial screening was performed.  
  
However, if a child does have a hearing loss, it is very important to find out as early as possible. That you bring your baby in for this second screening.  
  
If we have not already contacted you to set up an appointment, please call our office at <hospital phone number>  
  
Sincerely,  
<footer>



When using Microsoft Word to generate letters it is strongly recommended that each template does not exceed an uncompressed file size of 900 KB. The Microsoft Word “docx” format used by HiTrack is automatically compressed so the file size of the letter template itself cannot be used to determine if it is too large. For more information and support about letter template size please call the Help Desk.

Older versions of Office require Microsoft Office 2003 Compatibility Patch to process HiTrack's Letters. Please go to the Microsoft Office website and download any service packs for your Office program before installing the patch. Microsoft provides the patch at: [Patch File](#)

### **Auto-Fills**

If your site manually enters baby demographics and screening data the Auto-Fill feature can simplify data entry. For example, if your hospital uses a hand-held A-ABR device you may want to use A-ABR as an auto-fill for Test Type. This feature reduces the number of fields that you and your staff need to fill in by hand.

### **Milestones**

These three settings allow program coordinators to customize at what age babies should qualify for screening, loss identification and intervention enrollment milestones.

### **Other Settings**

These are other settings that can be changed and customized such as Required Fields, starting date defaults, birth weight entry format and advanced feature settings. The question mark icon to the left of each setting will explain the purpose of the setting. The "Reset" button to the right of each setting value allows you to reset the setting to the default value.

Commonly Configured Settings include:

Show Quick Add
Birth Weight Format
Auto Export PI (Program Information) (Windows user interface only)
Baby List Starting DOB range

### ***User Specific Settings***

The setting categories below your login name allow you to customize features for your login. The Auto Fills, Default Paths and Other sections listed here are similar to the settings found in the Facility specific category however, Automatic Email and Password are unique to the user specific settings. Remember to use "Save" after making changes.

### **Automatic Email**

Some features in HiTrack can generate an email to be sent out to interested parties. For example, when a user does a State Data Transfer this feature can automatically send an email to a

predefined state contact. Users can set up their email in HiTrack by SMTP or MAPI. For SMTP configurations, please contact your email provider. If MAPI is used, the computer's default MAPI-compliant e-mail program such as Microsoft Outlook, Outlook Express, Mozilla Thunderbird, etc... will be used. To enable Automatic E-Mail Generation for state data transfers, please check the corresponding box, and supply a valid contact email address. This address is where the automatic emails will be sent.

### **Password**

This is where a user can choose to change their password. Passwords cannot be reused within 12 months.

## ***Sharing Program Information with Screening Software***

Your lists of Doctors, Screeners, Audiologists and other pick lists are collectively called "Program Information". Your Program Information is always entered and updated in HiTrack as described in the previous section. Although it is easy to create infant records directly in HiTrack, most screening programs find it more convenient to collect demographic information and screening results using a screening software.


To reduce data entry time and improve data entry accuracy, several manufacturers of screening software have designed their systems to import the Program Information from HiTrack. This Program Information is made available to screeners through pop-up prompts or drop-down menus as they create infant records and document screening results. HiTrack can often be installed on the same computer as the screening software or it can reside on another computer if that is more convenient for data management purposes.

To create the necessary link between HiTrack and your screening equipment software, refer to the [Linking HiTrack with Equipment Manufacturer Software](#) section. Find the correct instructions based on the type(s) of data collection screening software you are using and whether the screening software and HiTrack are installed on the same or different computers.

## Tracking Features Overview


### The Filing Cabinet

HiTrack organizes your EHDI data to help you track the progress of babies needing attention in your database. This is done by categorizing babies into folders representing what they *need next* in the flow of follow-up. At the broadest level baby records are considered "Complete" with no further action needed or considered "In Progress" with on-going recommendations.



The File Cabinet organizes babies into drawers and folders representing these categories and provides powerful "To Do List" filtering based on the protocol setting for your facility. After a drawer and folder are selected the HiTrack Baby List is updated to show babies for that folder.

### Action Needed Drawer



The folders within this drawer represent actions for babies with tracking in progress based on your protocol. For example: The [Two-Stage Automated Pass/Manual Non-Pass](#) protocol requires manual review for all infants who do not pass screening. In other words you must "Finish" each record by hand to complete the screening step. So as baby records

The counts shown in the Action Needed folders show how many babies are in each action category based on the current Baby List DOB filter.

The [Two-Stage Automated Pass/Manual Non-Pass](#) protocol includes the following folders under Action Needed:

Inpatient Screening	<i>(Shows all babies in Inpatient subcategories)</i>
Need Inpatient Screening	Babies yet to be screened
Incomplete Inpatient Screening	Babies in process of being screening waiting to be finished.
Outpatient Screening	<i>(Shows all babies in Outpatient subcategories)</i>
Need Outpatient Appt.	Babies who have finished Inpatient screening without passing
Need Outpatient Screening	Babies with appointments for Outpatient Screening and babies who have not passed Outpatient screening that need to be finished.

Evaluation	<i>(Shows all babies in Evaluation subcategories)</i>
Need Evaluation Appointment	Babies needing appointments for diagnostic or other evaluation
Need Evaluation	Babies with appointments for evaluation
Risk Monitoring	Babies who've passed screening but are being monitored because of Risk Indicators
Needs Intervention	Babies that need Early Intervention
Need to Locate	Babies with invalid contact information who need follow-up

## No Action Drawer

After screening, most babies will need no further action. Infants who do not need further action are categorized into the following folders:

Complete	<i>(Shows all babies in Complete subcategories)</i>
Passed Screening	Babies that passed Inpatient or Outpatient Screening
Evaluation Stopped	Babies that have completed evaluation and need no further follow-up. e.g. Normal Hearing records.
Not Evaluated	<i>(Shows all babies in No Evaluation subcategories)</i>
Contacted but Unresponsive	This matches the CDC category for documenting why a baby did not receive a hearing evaluation.
Unable to Contact	(see above)
No Eval. – Medical Reasons	(see above)
PCP did not Refer for Evaluation	(see above)
No Services	<i>(Shows all babies in No Services subcategories)</i>
Contacted but Unresponsive	This matches the CDC category for documenting why a baby did not receive Early Intervention services.
Unable to Contact	(see above)
Transferred Out	<p>Babies that have been transferred out of your facility</p> <p>The Transferred Out folder works in one of two ways depending on your <a href="#">database layout</a>. When using a Centralized database (i.e. <a href="#">Web Companion</a>) the "Transferred" Baby List folder includes a read only view of babies that have transferred out of your Facility that are being listed under a different facility. (The Centralized Database setting can be changed in <a href="#">Advanced Settings</a>.) In Centralized databases the total number of infants shown in this folder may be more than the number of infants at your facility</p> <p>If you access HiTrack using the Windows Client you typically are participating in a <a href="#">Distributed database layout</a>. In this case the records marked as being transferred are included in your periodic State/Regional transfer. These babies are shown in the "Transferred" Baby List folder and are marked as "Transferred Out" because they are not being tracked in your portion of the Distributed database.</p> 

Refused	Babies whose parents have refused screening or evaluation
Discontinued	Babies whose tracking has been discontinued
Deceased	Babies who died

The File Cabinet folders are based on the [Protocol](#) selected for your facility. For example, automated Protocols, such as "Two-Stage Automated", omit the "Incomplete Inpatient" folder because this protocol anticipates completed screening data for babies merged from screening equipment software. This allows HiTrack to give you the most efficient data management support and keeps your data in convenient "To Do" lists.

**All Folders**

The "All Folders" drawer at the top left of the File Cabinet gives you ready access to a complete, unfiltered view of the Baby List.

**Baby List**

The Baby List is the heart of HiTrack and is a launching point for baby tracking. The Baby List is loaded after selecting a folder from the File Cabinet:

**All Data for Green River Regional (14 of 14)**

DOB From: 4/1/2010 To: 10/10/2012  Advanced Find  Milestones

Baby Name	Medical ID	Birth Date	Inpatient	IP Date	Outpatient
Mortonson, Emily	9568421	10/8/2012			
Douglas, Sally	8974652	10/8/2012	Referred	10/8/2012	Referred
Hatch, Felicia	8723691	10/6/2012	Inconclusive	10/8/2012	Referred
Rose, Kevin	8971246	10/5/2012	Referred	10/8/2012	
Stevens, Rebecca	49876324	10/3/2012			
Willard, Kevin	7984632	10/1/2012			No Out. Scrn->Eval
Smith, James	9451278	10/1/2012	Referred	10/8/2012	
Douglas, Cammie	8974521	10/7/2012	Passed	10/8/2012	
Fredricks, Cloey	8974563	10/7/2012	Passed	10/8/2012	

This list of babies is generated based on a specified range of birthdates and the current facility or facility group. Additional columns may be viewed using the scroll bar (not shown) at the bottom of the list. The list can be refreshed at any time by clicking the "Refresh" button.

## Sorting

The baby list may be sorted by clicking on the desired column name at the top of the list. Repeated clicks to the same column header will toggle between ascending and descending order.

## Filters

By selecting the Advanced Find box you can search using Medical ID, Last Name, First Name, Result, Birth Facility, Probable Duplicates, Notes, and even by appointments. Deselect the box to return to basic filtering. Auto-Filter provides additional flexibility for arranging your data in the Baby List by allowing column specific filtering.

## Diagnostic Result Filters

HiTrack's Baby List provides four diagnostic filters to help ensure accurate follow-up and documentation. Three of them are mutually exclusive and one of them is not.

The first three filter options are based on Hearing Disposition entries. If any Hearing Disposition entries have been made a baby will be included in one of these three baby list result filters:

Loss Confirmed	Hearing Disposition with a type and degree of hearing loss specified. Confirmed Date must be entered.
Normal Confirmed	Hearing Disposition with "N/A No Loss" entered. Confirmed Date must be entered for both ears.
Unconfirmed	Hearing Disposition entered with a blank Confirmed Date

These options are mutually exclusive; a baby will only show once for one of these three filter categories if disposition is entered.

The fourth diagnostic-related filter is "Has Diagnostics". This filter is not mutually exclusive. A baby will be included in this filter if they have any of the Diagnostics/"Manage" entries. This means that babies shown in this filter could have any possible screening or evaluation outcome, including Pass Screening, Unconfirmed Confirmed hearing and others.

Although in some jurisdictions it may be preferable to have diagnostics entered in conjunction with a baby's hearing disposition, it is not a requirement in HiTrack.



"How can I find records that are missing documentation of diagnosis?" While using the Has Diagnostic filter refer to the Messages area in the yellow File Card area to see what their Hearing Disposition outcome is. This will help you see records that may need additional documentation.



The "Refresh" button will blink red when filter criteria has changed to remind you to reload the list.

## Returning to the List

After reviewing a baby record or working in another part of HiTrack you can return to the last Baby List by clicking on the Return to Baby List icon at the top of the screen.



## Facility Group Baby List

When a Facility Group is selected in the Toolbar the baby list is not automatically loaded. This allows Central coordinators, who may work with thousands of records, to first specify their Baby List filter criteria and then load the list.

The Transferred folder Baby List also works in a different way when using a Facility Group. This list includes babies transferred from any facility included in the Facility Group.

## Combining Duplicate Records

De-duplication is an important Tracking task simplified for you by HiTrack with the Combine feature. By using this feature, all related record details from the two records will be permanently combined, including contact info, screening results and notes.

To combine two records from the Baby List follow these steps:

1. From the Baby List, select the baby record you would like to discard after combining.
2. Fill in the checkbox for "Combine" located at the bottom of the Baby List Screen.

This will show two new fields to the right. The first field shows the first baby that was selected.

3. Select the duplicate record from the list. Note that this fills in the second field.
4. Click on the "Preview" button after both records have been selected.

(Uncheck the Combine box to restart the process.)

5. HiTrack will then provide a side-by-side view of the records that are going to be combined.

Use the checkboxes in the list on the left to identify the field values you want to keep in the combined record. Leaving any field unchecked will discard the value.

Use the "Combine" button to confirm the record combination.

This powerful feature can be granted or suppressed as a specific access right while customizing [User Groups](#).

HiTrack also provides Data Cleaning Reports to help identify records that are probable duplicates. Please refer to information on [these reports](#) in Appendix A.

## The Baby Basket

The Baby Basket is a feature that allows you to work with more than one baby at a time. It can be shown by clicking on the blue bassinet icon to the left of the tool bar. The Baby Basket will appear on the right side of the screen.

Babies can be added to the basket by dragging the desired baby row from the Baby List onto the Baby Basket list. Multiple babies can be added by holding CTRL or Shift on the keyboard as babies are selected to drag to the Baby Basket.

The arrows located beneath the Basket list also allow you to work with the list. The button with double arrows pointing to the left clears the Baby Basket. The button with double arrows pointing to the right adds all babies from the current Baby List to the Baby Basket. The buttons with single arrows will add or remove the currently selected baby.

Once there are babies in the basket HiTrack provides a variety of actions that can be activated by clicking "Proceed".

**"Individual Diagnostic Report"** will generate a summary report for each baby in the Basket.

**"Create Data Transfer File"** will create an encrypted data file for a State/Regional transfer based on the babies in the Basket. This action is only available in the Windows version of HiTrack.

**"Create Unencrypted XML File"** will create an unencrypted XML file based on the information of the babies in the Baby Basket.

**"Email Medical IDs"** option will create an email using the values found in the Automatic Email link inside of settings. The body of this email will be populated with the medical IDs of the infants that are in the Baby Basket. This feature is for Windows only. For the web version of HiTrack there is the Get Baby IDs option that will populate a text box with the IDs of the Baby Basket babies. Then the user can copy and paste these values into their email client of choice.


**"Add Note/Custom Field"** option allows a common note to be added for each baby in the Basket using just one step.

**"Assign Physician"** assigns the specified physician to all babies in the Baby Basket.

**"Request Letter"** requests the specified letter for all babies in the Baby Basket.

**"Add Recommended Action"** will add the specified recommended action to all babies in the Baby Basket.

**"Change Tracking Folder"** will move all babies in the Baby Basket to a chosen folder. In order to use this option, all babies in the Baby Basket must currently be in the same tracking folder.



**“Delete All”** will delete all babies in the Baby Basket permanently from the system. This action cannot be undone.

The Baby Basket can be printed by selecting "Print" from the Right-Click menu inside the list itself.

The Baby Basket is also a good way to "park" baby records that need additional data entry. Later, after their data entry is done, the baby can be removed from the Basket.

## Collecting Screening Information using HiTrack

This section describes how the Inpatient and Outpatient folders are used when your screening equipment has little or no data collection capabilities (such as handheld devices). If you use a different type of screening equipment software, please refer to the section that describes [collecting screening information using screening equipment software](#).

The steps below will cover how to enter basic demographic data on babies who need to be screened, add screening results and finish completed records. Baby records should generally remain in the Inpatient folders until screening is completed or the infant is discharged.

### Facility Protocol Setting

When using HiTrack to enter screening results from a hand-held device you will need to use one of the following Protocols: "Two-Stage Automated Pass/Manual Non-Pass", "Two-Stage Manual", "Two-Stage Dual Equipment Auto-Pass Only" or "Single-Stage Dual Equipment Auto-Pass Only" [See Facility set up](#) and [protocol descriptions](#) for more information.



### Other Settings

Quick Add should be turned on for your Facility ([see other settings](#).) because it allows screeners to easily add screening results and finish baby records directly from the Baby List in the Windows user interface.

### Procedure Overview:

For Inpatient Screening you will use the "Inpatient Screening List". (The name of your current list is shown above the filter criteria on the Baby List.) If you are on a different list use the File Cabinet to select "Action Needed" and then the "Inpatient Screening" folder before starting your screening procedures.

**Step 1: Creating baby records.** The screener or data entry clerk creates a new record for each infant born or transferred into the facility by entering basic infant demographic information.

**Step 2: Printing a list of babies to be screened.** The screener prints a list showing all babies to be screened. The screener then manually records screening results on the printed list. If screening is performed where the screening software is located, the screener does not need to print a list; results can be entered directly into HiTrack as explained in Step 3.

**Step 3: Entering screening results and finishing screening records.** The screener adds screening results to an individual infant's record. The screener "finishes" an infant record indicating that the screening has been completed successfully or that the infant has been discharged from the hospital by using the Tracking Options link. Infants should remain on the browse list until they pass the screening or are discharged.

**Step 4: Entering outpatient screening results in a Two-Stage Protocol.** When an infant returns for Outpatient screening, the screener locates the baby record on the Outpatient Screening list and updates it with Outpatient results. The updated record is then "finished".

**Note:** Screening records may be entered manually or merged in from your screening equipment. See Tracking Tasks, Step 1, for details.

## Step 1: Creating baby records

For every infant born or transferred into your program, you must create a new child record. Click on a field or press the Tab key (to move forward) or press Shift-Tab (to move backward) to select a field.

To add a new baby record:

1. Click on the "Add" button at the bottom of the [Baby List](#).

A blank Demographics screen is shown.

The screenshot shows a software window titled "Baby Demographics". It contains two main sections: "Baby" and "Primary Contact". The "Baby" section has fields for Medical ID, Name-Last, First, Middle, Birth-Date (01/29/2013), Time (00:00), Order (Single), Weight, Gender (-Blank-), Birth Facility (Green River Medical), Nursery ([blank]), Location, Insurance ([blank]), Gest. Age, Physician ([blank]), Race ([blank]), ICU Days, Alt. ID, and Deceased. The "Primary Contact" section has fields for Name-Last, First, Middle, Medical ID, Address, City/State/Zip, Phone, Ph 2, Email, County, Birth Mother (checked), Birth Date, Language ([blank]), Title/Suffix, Race ([blank]), Gender (-Blank-), Education ([blank]), Ethnicity ([blank]), and Deceased. A "Save" button is at the bottom right.

2. Enter the baby's information in the appropriate fields. Colored fields indicate mandatory data.

The data field drop-down menus offer preset options that standardize data entry. Some fields might be auto-filled. If the auto-filled information is incorrect, simply enter the information for the blank fields.

3. Enter the Primary Contact's information in the appropriate fields

If the Primary Contact is someone other than the birth mother (e.g., foster or adoptive parent), uncheck the "Birth Mother" box. (To add more contacts return to "Demographics" after saving.)

4. Click the "Save" button at the bottom of the screen

The Baby List is shown with your new baby highlighted. This baby's yellow File Card will be shown on the left.

*To enter additional notes:*

1. If the baby File Card is not shown on the left select the baby from the Baby List.
2. Click on "Demographics" from the File Card  
This baby's Demographics screen is shown.
3. Click on the "Notes" tab above the demographics screen or click on the yellow and blue Note preview box.  
This baby's Note history is shown. If no notes are present a "New Entry" field is selected and the cursor is set for you to type.
4. Type in your note. The date and time are automatically stored with each note entry. In HiTrack 4 each note should be entered as a separate new entry unless you are making revisions to a prior note.
5. Click the "Save" button at the bottom of the screen

To enter additional contacts:

1. If the baby File Card is not shown on the left select the baby from the Baby List.
- 2a. Click on "Demographics" from the File Card on the Windows user interface.  
This baby's Demographics screen is shown.  
Click on the "Contacts" tab above the demographics screen.
- 2b. [Web Companion](#) users should click on the Contacts link on the File Card.
3. A list of all contact roles is shown. Persons may be shown more than once in this list because some contacts play more than one role. For example "Birth Mother" and "Primary Contact."
4. Click on "Add"  
A blank Person Demographics screen is shown.

The screenshot shows a software window titled "Demographics" with three tabs: "Demographics", "Contacts", and "Notes". The "Contacts" tab is active. The window is divided into two main sections: "Contact Info" and "Contact Role".

**Contact Info:**

- Name-Last: [text input]
- First: [text input]
- Middle: [text input]
- Medical ID: [text input]
- Address: [text input]
- City/State/Zip: [text input]
- County: [text input]
- Phone: [text input]
- Ph 2: [text input]
- Email: [text input]
- Birth Date: [checkbox] [text input]
- Language: [blank] [dropdown]
- Title/Suffix: [text input]
- Race: [blank] [dropdown]
- Gender: [-Blank-] [dropdown]
- Education: [blank] [dropdown]
- Ethnicity: [blank] [dropdown]
- Deceased: [checkbox]

**Contact Role:**

- Primary Contact
- Other Contact
- Birth Mother
- Father
- Legal Guardian
- Kin
- Foster Parent
- Parent
- Primary Care Provider
- Attending Physician
- Attending Nurse
- EI Service Provider
- Other Health Care Provider
- Medical Home

At the bottom of the window are "Save" and "Cancel" buttons.

5. Enter the contacts information in the appropriate fields and select appropriate contact role types from the checkboxes below the demographics fields. Colored fields indicate mandatory data.  
Select "Primary Contact" (in addition to any other roles you select) if this is the person that should receive mail and contact for this baby.
6. Click the "Save" button at the bottom of the screen

*Selecting Risk Indicators:*

1. If the baby File Card is not shown on the left select the baby from the Baby List.
2. Click on "Risk Indicators" from the File Card

A blank Risk Indicator list is shown.

JCIH 2007		Other Indicators
Result	Indicator	Detail
Unknown	Caregiver concern regarding hearing, speech, language, or developmental delay.	
Unknown	Family history of permanent childhood hearing loss.	
Unknown	Neonatal intensive care (All infants with or without risk factors) for greater than 5 days, including any of the following: ECMO,* assisted ventilation, exposure to ototoxic medications (gentamycin and tobramycin) or loop diuretics (furosemide/lasix) In addition, regardless of length of stay: hyperbilirubinemia requiring exchange transfusion.	
Unknown	In utero infections such as CMV, herpes, rubella, syphilis, and toxoplasmosis.	
Unknown	Cranio-facial anomalies including those that involve the pinna, ear canal, ear tags, ear pits, and temporal bone anomalies.	
Unknown	Physical findings such as white forelock, that are associated with a syndrome known to include a sensorineural or permanent conductive hearing loss.	
Unknown	Syndromes associated with hearing loss or progressive or late-onset hearing loss, such as neurofibromatosis, osteopetrosis, and Ushersyndrome; other frequently identified syndromes include Waardenburg, Alport, Pendred, and Jervell and Lange-Nielson.	
Unknown	Neurodegenerative disorders such as Hunter syndrome, or sensory motor neuropathies, such as Friedreich's ataxia and Charcot-Marie-Tooth syndrome	
Unknown	Culture-positive postnatal infections associated with sensorineural hearing loss, including confirmed bacterial and viral (especially herpes viruses and varicella) meningitis.	
Unknown	Head trauma especially basal skull/temporal bone fracture that requires hospitalization.	
Unknown	Chemotherapy	

3. Using the drop-down menus, select either Yes, No or Unknown to indicate if the specific risk criteria are present.
4. To enter a Risk Indicator not listed in the JCIH 2007 list, select the "Other Indicators" tab.
5. Click the "Save" button below.

## **Step 2: Printing a list of infants to be screened**

After you create infant records in Screening, you can print a list of all infants and the respective screening actions required. This list includes the infants whose records you entered today, as well as those entered previously who have not yet passed the screening. Carry this printout with you to record screening outcomes.

If you have direct access to the screening software while conducting screening, you can conduct the screening and enter results directly in the screening software as explained in Step 3.

*To print a list of babies needing Inpatient screening:*

1. Verify that you are on the "Inpatient Screening List" (shown above the baby list filter criteria.)
2. Sort the list in a way that is helpful for your screening staff by clicking on the header column for "Baby Name", "Medical ID" or "Birth Date". A second click will reverse the sort order.
3. Click on "Print List" at the bottom of the Baby List to see a preview and then "Print" when your printer is ready.

*To print a list of babies needing Outpatient screening:*

1. Verify that you are on the "Outpatient Screening List" (shown above the baby list filter criteria.) If you are on a different list use the File Cabinet to select "Action Needed" and then the "Outpatient Screen" folder.
2. Sort the list in a way that is helpful for your screening staff by clicking on the header column for "Baby Name", "Medical ID" or "Birth Date". A second click will reverse the sort order.
3. Click on "Print List" at the bottom of the Baby List to see a preview and then "Print" when your printer is ready.

### Step 3: Entering screening results and finishing screening records

At the end of a screening session you should enter the screening results for all the infants you have screened.

To enter screening results:

1. If the baby File Card is not shown on the left select the baby from the Baby List.

[Quick Add](#) will be activated below the list.

Method	Result	Type	Date
IP Right (A-ABR 1st)	Refer	TEOAE	5/30/2002
IP Left (A-ABR 1st)	Could Not Test	TEOAE	5/30/2002
IP Head (A-ABR 1st)	Refer	TEOAE	5/30/2002

Note that "Inpatient" is pre-selected as the Stage for this baby.

2. Enter screening results for the selected baby in the fields on the "Results" tab. If your hand-held device provides a test ID or number you may enter it into the "ID" column. (Go directly to the "Finish" for infants who were not screened.)
3. If screening has been completed for this baby click on the "Finish" tab. ([Web Companion](#) users should click on the [Tracking Options](#) link.)

Finish Options will be shown on the Windows user interface.

Method	Result	Type	Date
IP Right (A-ABR 1st)	Refer	TEOAE	5/30/2002
IP Left (A-ABR 1st)	Could Not Test	TEOAE	5/30/2002
IP Head (A-ABR 1st)	Refer	TEOAE	5/30/2002

4. Use the Outcome drop-down list to pick the appropriate outcome for this baby. The available options in this list will change depending on what has been entered for the record. For example if no results are entered the "Missed" option will be shown.
5. If desired add a short note explaining reason of the outcome. This is optional.
6. Click on the "Save" button below.

This baby will be moved out of the Inpatient Screen folder.

#### Step 4: Entering Outpatient screening results in a Two-Stage Protocol

To enter in Outpatient screening results use the File Cabinet to select "Action Needed" and then the "Outpatient Screen" folder.

*To enter outpatient screening results:*

1. If the baby File Card is not shown on the left select the baby from the Baby List.  
[Quick Add](#) will be activated below the list.  
"Outpatient" is pre-selected as the Stage for this baby.
2. Enter screening results for the selected baby in the fields on the "Results" tab. If your hand-held device provides a test ID or number you may enter it into the "ID" column.  
(Go directly to the "Finish" for infants who were not screened.)
3. If screening has been completed for this baby or the baby was not screened click on the "Finish" tab.

([Web Companion](#) users should click on the [Tracking Options](#) link.)

Finish Options will be shown.

Method	Result	Type	Date
OP Right (A-ABR 1st)	Refer	A-ABR	6/1/2002
OP Left (A-ABR 1st)	Pass	A-ABR	6/1/2002
OP Head (A-ABR 1st)	Refer	A-ABR	6/1/2002

Outcome: Finish OP with Evaluation Needed

Reason:

Save

4. Use the Outcome drop-down list to pick the appropriate outcome for this baby. A "Finish OP with Evaluation Needed" option will be present if this baby referred.
5. If desired add a short note explaining reason of the outcome. This is optional.
6. Click on the "Save" button below.

This baby will be moved out of the Outpatient Screen folder.

## **Collecting Screening Information using Equipment Software**

You can merge screening information into HiTrack electronically if your hearing screening program uses a screening system with its own compatible screening software. Other screening systems software changes frequently so if you have additional questions beyond what is covered here in the manual, please contact the screening manufacturer.

### **Facility Protocol Setting**

If you plan to manually review non-pass babies imported from an equipment file then you should use one of the following protocols: "Two-Stage Automated Pass/Manual Non-Pass", "Two-Stage Dual Equipment Auto-Pass Only" or "Single-Stage Dual Equipment Auto-Pass Only"

If you merge completed screening results from an equipment file then you should use one of the following protocols: "Two-Stage Automated", "Two-Stage Inpatient Automated", "Two-Stage Dual Equipment Automated" or "Single-Stage Dual Equipment Automated"



[See Facility set up](#) and [protocol descriptions](#) for more information.

### **Facility Merge Sources**

While configuring HiTrack for your facility you should have added [Merge Sources](#) for the screening equipment you use. Each file imported into HiTrack is referred to as a "Batch" of screening equipment information. After merging the file will be renamed with a "Time Stamp" to indicate that this file has been imported.

### **Matching Rules**

In most cases Equipment merge files contain the newly collected screening data. However, in some scenarios, such as Outpatient follow-up screening, additional screening data may be collected. To account for this HiTrack looks for matches on Medical ID, birth date and last name. If this information does not completely match an existing infant record, the information is placed in the Pending Batches list. This allows the program coordinator to decide which field needs to be corrected and then merged into an existing record or added as a new infant record, as the case may be. (See [Appendix D](#) for complete explanations of matching criteria.)

### **Procedure Overview:**

**Step 1: Collecting Screening Information using equipment software.** The screener uses the screening equipment software to enter baby demographic information and collect screening results. An "export" or "save" feature is used in the software to create an equipment file for HiTrack.

**Step 2: Starting the Merge process.** The screener or person in charge of tracking uses the Merge icon to open the Merge screen, selects the Merge Source and begins the merge for this batch.

**Step 3: Resolving any pending Merge batches.** The screener or person in charge of tracking uses the Merge screen to fix invalid data or resolve questionable matches.

## Step 1: Collecting Screening Information using Equipment Software

This step varies because each screening equipment software system works differently. Some equipment software provides complete demographic and screening fields while other equipment only allows for a handful of minimum required fields. Each equipment software system also writes the equipment file at different times during the screening. Please refer to the section describing data entry into [Equipment Manufacturer Software](#) for details about specific software packages that have been tested with HiTrack.

These steps are specific to HiTrack's Windows user interface. Refer to the [Web Companion section](#) for steps on how to use merge equipment files into HiTrack's Web Companion.

## Step 2: Starting the Merge process

The Merge can start after using your screening equipment software to create an equipment file.

*To merge an equipment file*

1. Click on the "Merge" icon from the Toolbar

The Merge screen is shown.

Facility	Start	Stop	Babies Ad	Babies Up	Total Add	Total Upd	Status
Green Rive	01/29/2013	01/29/2013	1	0	7	3	Merge Completed

2. Use the Merge Source drop-down to select your equipment. The Source Location for your merge source will be shown.

The "Merge" button is activated if a file from the selected equipment is found. If it is not found, the "Merge" button's text will become "Search" indicating that a search is required to find your specified file. Clicking on "Search" will automatically search the root directories of your system for the file. If it is not found, it will allow the user to search manually for the file's location. This search feature does not allow users to search


for XML, ETF, or ETF4 files.

If the "Merge" button is not activated after selecting your equipment, please refer to the [Equipment Manufacturer Software](#) section for instructions on creating the equipment file.

3. Click the "Merge" button.

The "Merge Results" screen is shown and the batch is processed.

After the merge process successfully completes a "Merge Complete" message is shown. The number of rows added and updated is also reported on this screen.



4. If the process completes but some Pending data must be resolved a different message will be shown; "*Some data must be manually resolved in this batch. Use Pending Data to finish work on this batch.*"

This batch will remain in Pending data to be resolved. See Step 3 for instructions on resolving batches.

5. Click the "Ok" button to return to the main Merge Page.


Merge History will be updated with the results of this batch.

### Step 3: Resolving any pending Merge batches

When a batch needs manual resolution HiTrack keeps the batch in the Pending Data list. These batches are temporary and can be deleted by using the "Delete Batch" button.

1. Select a pending batch from the "Pending Data" grid on the Merge screen.

The Browse Pending Data screen is shown.



2. Baby records shaded in purple either (a) have data entry problems or (b) are part of a questionable match. After you click on a baby the Merge Item Resolution screen is shown.

(continued)

- 3a. A single column of fields is shown when invalid entries need to be corrected. Invalid birth dates and birth weights must be resolved. Blank deceased dates also require manual resolution. Click on the field to make a correction and click on "Save" to retain the changes.
- 3b. Two records are presented side by side for questionable matches. The Existing infant on the left is the record already in HiTrack and the Incoming infant is the record loaded from the current batch. The Merged Record shows you the completed infant record based on the checked 'Use Incoming' column.

Browse Pending Records | Resolve Pending Record

Advanced Overwrites

Smit BB 12/15/2012 12:00:00 AM Merge with this Candidate

Not Matching: last name: SMIT, 92.86 %

Field	Existing	Incoming	Use Incoming	Merged Record
IDNumber	985743	985743	<input type="checkbox"/>	985743
last name	Smit	Smith	<input type="checkbox"/>	Smit
first name	BB	BB	<input type="checkbox"/>	BB
birth datetime	12/15/2012 12:00:00 AM	12/15/2012 12:00:00 AM	<input type="checkbox"/>	12/15/2012 12:00:00 AM
multiple	N		<input type="checkbox"/>	N
birth order	0		<input type="checkbox"/>	0
birthing facility id	Green River Medical	Green River Medical	<input type="checkbox"/>	Green River Medical
gender	F	F	<input type="checkbox"/>	F
brn			<input type="checkbox"/>	
deceased	N	N	<input type="checkbox"/>	N
weight			<input type="checkbox"/>	
total days in any nicu			<input type="checkbox"/>	
gestation period			<input type="checkbox"/>	
nursery level code			<input type="checkbox"/>	
current location			<input type="checkbox"/>	
middle names			<input type="checkbox"/>	
responsible facility	Green River Medical	Green River Medical	<input type="checkbox"/>	Green River Medical
tracking status	In Progress	In Progress	<input type="checkbox"/>	In Progress
change datetime	01/29/2013 11:51:43 AM	01/29/2013 11:51:48 AM	<input checked="" type="checkbox"/>	01/29/2013 11:51 AM
source facility	Green River Medical	Green River Medical	<input type="checkbox"/>	Green River Medical
deceased date dy			<input type="checkbox"/>	
deceased date mo			<input type="checkbox"/>	

There are three possible actions to take for possible matches: "Merge with this Candidate", "Add as New Record" or "Delete Imported Record".

Use "Merge with this Candidate" when the records are for the same baby.

Use "Add as a New Record" when the record from the current batch is not the same baby shown on the right.

Use "Delete Imported Record" to discard the record from the current batch.

After choosing an action the Merge will process your request and show the Merge Results screen.

4. Click the "Resume Merge" button to process invalid data corrections or to continue merging this batch again.

## Collecting Diagnostic Hearing Information

HiTrack provides a few different ways to collect diagnostic hearing evaluation information. Data can be collected and recorded within the HiTrack user interface using the [Diagnostic Quick Add](#) and [Diagnostic Management](#) links or collected by audiologists without HiTrack access using the Hearing Evaluation Submission Form.

HiTrack's Hearing Evaluation Submission Form is a web based form that allows diagnostic data to be submitted by audiologist clinics without HiTrack access in version HiTrack 4.6.0. This form submits diagnostic data to the HiTrack Pending Merge list by authenticating the audiologist's name, state code and facility name. As the submissions are collected on HiTrack's Pending Merge list they can be reviewed for completeness and then automatically merged with the permanent HiTrack data or discarded.

## Set up and IT Configuration

This form is available to users that have [Web Companion](#) installed. A new [Facility](#) and associated [Merge Source](#) entry must be added to HiTrack by the state HiTrack user to collect the data. The Facility should be named "Diagnostic Submissions Data". The new Merge Source should be given a name and use the type "Clinic Data Source".

Once these items are added to HiTrack's system, Web Server IT staff are required to update a configuration file on the web server with the IDs of this new Facility and Merge Source.

## Matching Rules and Required Fields

A new set of matching rules is used to compare diagnostic submissions to existing HiTrack records. See below:

<i>Field</i>	<i>Required and Note</i>
Patient ID/Medical ID	Required. A state-wide unique identifier should be used here when available. Otherwise, the audiologist's local patient ID should be used here.
Child Last and First Name	Required and used to match.
Child DOB	Required and used to match.
Contact Last and First Name	Required and used to match.
Contact DOB	Required and used to match.
Contact Phone	Required and used to match.
Contact Zip	Required and used to match.
Child's Birth Facility	Required. Must match a HiTrack Facility name.
Audiologist Name and Code	Required. Both must match HiTrack's name entered in Medical People.
Diagnostic Facility	Required. Must match a HiTrack Facility name.

As with HiTrack's [other match rules](#), incoming records that match a majority of fields are automatically combined. Partial matches are left in the Pending List and records with no matches are added as new records in HiTrack. These new unmatched records will show in the Baby List under the new "Diagnostic Submissions Data" facility as records without screening results and can be manually matched to existing records at any time.

## **Audiologist Onboarding Process**

Below is a suggested process of introducing new audiologists to the HiTrack Hearing Evaluation Submission Form.

1. Determine if the audiologist is capable and willing to submit relevant HiTrack diagnostic data.
2. Verify their related entries in HiTrack:
  - Check that their name is spelled and entered correctly in HiTrack's Medical Person Edit screen.
  - Check that their clinic/facility/site name is entered correctly in HiTrack's Facility List
3. Provide authentication details in an email and/or phone call (see suggested text below):
  - The HiTrack Web URL (similar to) "[https://yourstate.gov/HiTrack\\_4/HiTrack\\_4/HiTrack\\_4\\_Web\\_PL/SubmitHearingEvaluation.aspx](https://yourstate.gov/HiTrack_4/HiTrack_4/HiTrack_4_Web_PL/SubmitHearingEvaluation.aspx)" Suggest that they bookmark this site after opening it the first time.
  - Their State Audiologist Code from the Medical Person Edit screen. Create a new code for them if they do not already have one.
  - The list of HiTrack Facilities to use to enter in for Birth Facilities and Diagnostic Facility. This can be generated by saving the Facilities list to an Excel file.
4. Arrange with them to submit a test record that can be discarded from your Pending List after successful submission.
5. Arrange for periodic follow-up with them about their data quality.

## Suggested Audiologist Welcome Message

Below is a message that could be sent to welcome new Audiologists to HiTrack's Hearing Evaluation Submission form.

Dear Audiologist,

As you may know, our state uses an EDHI data management system called HiTrack to track and maintain early child hearing results and diagnosis. A new feature in HiTrack enables offsite audiologists such as yourself, to submit hearing evaluations for patients you see to the state department of health.

We hope that this form will allow us to replace the FAX forms that are currently in place.

Types of data collected by HiTrack's Hearing Evaluation Submission form:

- Degree, Type and Configuration of hearing loss
- Diagnostic ABR and OAE tests
- Tympanometry tests
- Behavioral tests
- A short note

To use this feature you will need:

- A State Issued HiTrack Audiologist Code
- The HiTrack Facility List
- A web browser configured to link to:  
[https://yourstate.gov/HiTrack\\_4/  
HiTrack\\_4/HiTrack\\_4\\_Web\\_PL/SubmitHearingEvaluation.aspx](https://yourstate.gov/HiTrack_4/HiTrack_4/HiTrack_4_Web_PL/SubmitHearingEvaluation.aspx)

Please feel free to contact me to learn more about setting up this form for you,

Thank You,

State EHDI Manager Staff

## Configuration

Please contact the [Help Desk](#) to enable this form.

## Managing Screening and Tracking Data

HiTrack provides the tools and features you need to manage your hearing screening and tracking data. This section describes common tasks that apply these tools and features effectively.

### Toolbar

The toolbar is located at the top of every screen in HiTrack and provides easy access to all of HiTrack's features.



On the left side of the toolbar is the Facility or Facility Group that is selected after logging in. This can be changed at any time while logged in. The currently selected Facility or Facility Group determines which babies you will see, which Letters will be generated and which [Settings](#) are being used.

On the right side of the toolbar are several icons:

	<a href="#">Letters</a> : HiTrack generates letters for parents and health care providers. This icon opens the list of pending letters waiting to be printed.
	<a href="#">Reports</a> : HiTrack provides reports to monitor hearing screening and evaluation performance. This icon provides a list of reports.
	<a href="#">Merge</a> : HiTrack can import data from screening equipment software, hospital data systems and other HiTrack databases. This icon opens the merge tools used to import data from these sources.
	<a href="#">Tools</a> : HiTrack features powerful tools to help manage your EHDI data. This toolbox icon opens a tool screen organized into helpful tabs.
	<a href="#">Settings</a> : This dial icon opens the Settings screen to allow customization and administration of HiTrack.
	<a href="#">Log out</a> : The exit door Log out icon will log a user out of HiTrack and return the Login screen.

## Tracking Tasks

As discussed in the [Tracking System Overview](#) the File Cabinet and Baby List allow you to manage your baby HiTrack data. The following sections will discuss common tracking tasks in detail.



Tracking tasks start at different times depending on your tracking role. State or regional users interact with HiTrack at different times than hospital users.

Some common tracking roles include:

Hospital users	Tracking Tasks begin after screeners have collected screening information and merged or entered it into HiTrack
State or regional users who collect data from multiple HiTrack databases in their state or region.	Tracking Tasks begin after the hospital Data Transfer files (ETF or ETF4 files) have been merged into their state or regional HiTrack system.
State or regional users who maintain data from multiple hospitals in a single, <a href="#">centralized database</a> .	Tracking Tasks begin after screeners have collected screening information and merged or entered it into HiTrack.  Encrypted Transfer Files do not need to be merged for state/regional users providing centralized data with HiTrack's Web Companion

Contact the [HiTrack Help Desk](#) for additional consultation about your tracking role.

### Procedure Overview:

**Step 1: Updating baby records.** Records for babies who do not pass screening are updated with additional screening results, demographics and recommendations. All in-progress records are finished and moved to the next action folder in the File Cabinet.

**Step 2: Generating Letters.** Letters are printed to communicate screening results to parents (or preferred contacts) and physicians so that follow-up appointments can be arranged.

**Step 3: Reporting Tracking Effectiveness.** Reports are generated as needed to evaluate follow-up success, Refer rates and other tracking measurements.

**Step 4: Transferred Data to a Central Database (optional).** This step creates an encrypted transfer file (ETF) for you to report your HiTrack data to a state or regional coordinator. This is only required for [Distributed](#) state or regional databases.



## Tracking Tasks - Step 1: Updating baby records

Updating baby records is an important daily tracking task. Babies who referred from screening that can be located by using the File Cabinet and selecting the "Need Evaluation" folder. The first subfolder "Need Evaluation Appointment" provides a quick list of babies that need appointments for diagnostic evaluation.

After a baby has been selected from the Baby List their File Card is shown to the left.

### Baby File Card

The File Card provides access to information specific to the selected baby. Relevant key information, such as contact information, risk indicator alert and pending letters notification also appear on the File Card providing you with on-demand assistance with your tasks.

The File Card is separated vertically into four main areas. The top of the card shows identifying information for the current baby. The middle section provides links to different parts of the baby's record:

The Recommended Actions section changes as the baby moves from folder to folder and suggests relevant tracking actions.

The Messages section displays helpful information about the selected baby. Hearing Disposition, Contact information and baby specific Milestone achievements are examples of helpful information provided in the Messages section

Use a single click to show the File Card for a baby from the Baby List.

The screenshot displays the HiTrack 4 software interface. At the top, it says "All Data for Green River Regional (14)". Below this are buttons for "All Folders", "No Action", and "Action Needed (7)". A "Facility:" dropdown is visible on the right. The main content area is divided into two columns. The left column, titled "Mortonson, Emily", shows her Medical ID (9568421), Birth Date (10/8/2012), Action (Need Inpatient Screening), and Status (In Progress). Below this are several sections: "Record Overview" with links for Tracking Options, Demographics, Screening, Risk Indicators, Diagnostics | Manage, Hearing Disposition, Amplification, and EI Services; "Letters" and "Transfer" sections; "Recommended Actions: Manage"; and "Messages: Milestones: -/-". The right column shows "DOB From:" and "To:" fields, a "Baby Name" list with "Mortonson, Emily" selected, and a "Quick Add for" section with "Results" and "Finish" buttons, a "Stage:" dropdown set to "Inpatient", a "Date:" field set to "10/11/2", and a "Source:" dropdown set to "-Blank-". An "Add Baby" button is at the bottom right.

## Record Overview

Record overview and Advanced Details provide a summary of data for a baby and provides links to obtain more detailed information.

Note that on this screen the File Card has been minimized to the top to maximize screen usage. Clicking on the yellow icon at the top will display the baby File Card again. Click on any "More" link to see detailed information on that part of the record

The screenshot displays the 'Overview for Winger, Robin' record page. At the top, there is a navigation bar with buttons for 'All Folders', 'No Action', 'Action Needed (15)', 'Facility: Great Basin Medical Center', 'Letters', 'Reports', 'Merge', 'Tools', 'Settings', and 'Log Out'. The user is identified as 'dora'.

The main content area is divided into several sections:

- Demographics:** Name: Winger, Robin; Medical ID: 5432-33471; Birth Date: 8/17/2006; Gender: F; Entered: 11/27/2007; Action: Need Outpatient Appt.; Nursery: Well Baby; Resp Facility: Great Basin Medical Center; Birth Facility: Great Basin Medical Center; Source: Primary Care; Provider: SPACH, BENNY; Primary Contact: Winger, Mika.
- Most Conclusive Result:** A table showing results for IP Head, IP Left, and IP Right, all with a 'Refer' status and 'A-ABR' test type on 8/17/2006.
- Diagnostic Tests:** A table showing 'Left' and 'Right' ear tests for 'TEOAE' on 8/17/2006.
- Recommended Action:** 31 Aug 2006 Broken Appt Needs Outpatie.
- Hearing Status:** None.
- Letter History:** None.
- Notes:** Mother requested follow 11/9/2011 4:18 PM.
- Tracking History:** Broken Appointment 11/27/2007 3:01 PM; OP Appointment made 8/31/2006.
- Audiogram:** None.
- Risk Indicators:** Apgar scores (i) Yes.
- Enrollment:** None.
- Amplification:** A table showing 'Left' and 'Right' ear tests for 'Hearing Aid' and 'Cochlear Implan' on 11/19/2006.

## Advanced Details

Record overview also provides access to a baby's Advanced Details. Clicking on the magnifying glass icon  will open a pop-up window showing Advanced Details.



Advanced Details is a helpful tracking tool because it provides access to each baby's Change Log and Match History. These two historical records are important tracking tools because most hearing screening programs merge their data from other sources such as equipment software or other HiTrack systems. These features let you see changes that other users make during these processes.

The following Baby fields are maintained in the Change Log:

Medical ID	gender
birth date	middle names
last name	alternate id
first name	deceased (Y/N)
test result edit	test result delete

The following Person/Contact fields are maintained in the Change Log:

last name	gender
first name	address
last name	phone
middle names	deceased (Y/N)

The Change Log entries displayed in the following format:

Edit:	The name of the field that was changed.
Prior Value:	This is the previous value for this field. Empty brackets mean that the value was empty (null) prior to being edited.
Prior change date:	This date in parentheses is the last time this was changed before this change.
New Value:	This value preceded by " to " is the new value of this field.
Date:	This is the date the field was changed

A baby's Match History is also provided on Advanced Details. A match history entry is recorded each time a record is combined or matched by HiTrack's [Merge process](#). These entries follow the following format:

Matched: Match confidence score #	This will be blank for Combined records, however for <a href="#">merged records it is a score calculated</a> using the number of matching fields.
User Match: Yes/No	This indicates whether or not HT automatically matched this record. A "Y" value will be shown if this match was done with the Combine feature or if this record was matched by hand during the merge process.
User	The user who combined or matched the records
Incoming	This shows the HT ID# from the discarded record for technical purposes only.
Date:	This is the date the match or combine occurred.



Classifications are a third noteworthy item included in Advanced Details. [Classifications](#) are calculated for each baby at each stage and represent how they are counted on reports. They are also use as criteria for generating letters.

## Tracking Options

The Tracking Options link is directly below the Record Overview link on the baby File Card. Clicking this link will display the following screen.

Tracking options for Mortonson, Emily

User: Chris Help

Actions | Stop Tracking

Current Action: **Need Outpatient Appt.**

Possible Actions

To advance this baby select an action below.  Suggested Actions

Action
Need to Locate Child
OP Appointment made
Recall for Inpatient Screening
Evaluate w/o Outpatient

History

From	To	Explanation	Entered
Need Inpatient Screenin	Need Outpatient Appt.	Missed Inpatient	10/11/201

Reset & Clear Tracking Options Edit

## Tracking Actions

This tab will allow you to advance a baby record from one tracking action to another. The top grid gives a list of suggested actions for a baby based on the protocol being used. (You will see these same actions when using the Quick Add "Finish" feature on the Baby List.)

Some possible Inpatient actions for the [Two-Stage Automated Pass/Manual Non-Pass](#) protocol are:

Finish IP without Passing	The baby was screened but the results were inconclusive or incomplete.
Refer for Outpatient	The baby was screened but the pass criteria was not met for one or both ears; refer for further screening.
Missed Inpatient Screening	The baby was not screened before the mother and child were discharged from the hospital.
Evaluate w/o Screening	No further screening was done because the baby was referred directly for Audiologic evaluation.
Refer for Intervention	The baby has completed evaluation and will be monitored as enrollment activities take place.

The list of suggested possible actions will automatically change as screening information is entered. The protocol being used at your facility will also determine what actions are available. You may override the choices at any time for special circumstances by deselecting this option.

The bottom grid provides a history of actions that have occurred for the selected baby. For record keeping purposes you may edit the explanation field by selecting a history row and then clicking the "Edit" button.

### **Stop Tracking**

This tab allows you to deactivate a baby record. Possible reasons why tracking might be stopped are:

Refused	Parental Refusal of screening or follow-up activities
Discontinued	Returned Mail, Moved out of State, Multiple failed attempts to contact, Unresponsive parents

Always be sure to record a note explaining the reason why tracking was stopped.

## Demographics

The Demographics link in the baby's File Card provides specific demographic data for a baby and their primary contact.

Demographics for Mortonson, Emily

Facility: Green River Regional Letters Reports Merge Tools Settings Log Out User: Chris Help

**Demographics for Mortonson, Emily**

Demographics | Contacts | Notes

**Baby**

Medical ID: 9568421 Birth Date: 10/ 8/2012 Birth Facility: Green River Regional

Name-Last: Mortonson Time: 08 52 Nursery: Well Baby

First: Emily Order: Single Location:

Middle: Weight: 658 g Gender: Female

Insurance: [blank] Gest. Age: Physician: [blank]

Race: [blank] ICU Days: Alt. ID:

Deceased:  Notes: [\[Click here to enter notes\]](#)

**Primary Contact**

Name-Last: Mortonson Address: 123 E 400 N Phone: 435-512-1267

First: Susan Ph 2:

Middle: Email:

Medical ID: 0094 City/State/Zip: Vernal UT 84321 County:

Birth Mother:  Birth Date: Language: ENGLISH

Title/Suffix: Race: [blank] Gender: Female

Education: [blank] Ethnicity: [blank] Deceased:

Save

Change Date: 10/11/2012 8:10:25 AM

**Mortonson, Emily**

Medical ID: 9568421 Birth Date: 10/8/2012 Action: Need Outpatient Appt. Status: In Progress

Record Overview

Tracking Options

Demographics

Screening

Risk Indicators

Diagnostics | Manage

Hearing Disposition

Amplification

EI Services

Letters

Transfer

**Recommended Actions:**

Schedule Outpatient Screening Manage

**Messages:**

Milestones: -/- PC: Susan Mortonson Phone: 435-512-1267

The demographics tab provides a preview of the latest note at the bottom of the baby demographics box. Clicking on this box will switch you to the Notes data entry screen.

A Change Date is also shown at the bottom of the demographics screen. This shows the last time this baby's demographics were changed.

Colored fields indicate mandatory data. Click on a field or press the Tab key (to move forward) or press Shift-Tab (to move backward) to select a field.

*To update a baby record:*

1. Enter the appropriate values into the fields shown.  
A drop-down field marked with a pad-lock icon indicates that this record is using Program Information that is [retired](#) or not available to your facility. Click on the pad-lock to unlock the drop-down and select a current Program Information entry.
2. Click the "Save" button to retain your changes.  
The Change Log is automatically updated when identifying fields are modified.

*To add additional contacts:*

1. Click on the "Contacts" tab above the demographics screen. (The Web Companion provides a separate "Contacts" link on the File Card for this feature.)

A list of all contact roles is shown. Persons may be shown more than once in this list because some contacts play more than one role. For example "Birth Mother" and "Primary Contact."

2. Click on "Add"

A blank Person Demographics screen is shown.

The screenshot shows a web form with three tabs: 'Demographics', 'Contacts', and 'Notes'. The 'Demographics' tab is active. The form is divided into two main sections: 'Contact Info' and 'Contact Role'. The 'Contact Info' section contains several input fields: 'Name-Last', 'First', and 'Middle' (all highlighted in light blue); 'Address', 'Phone', 'Ph 2', 'Email', 'Medical ID', 'City/State/Zip', 'County', 'Birth Date', 'Language', 'Title/Suffix', 'Race', 'Gender', 'Education', 'Ethnicity', and 'Deceased'. The 'Contact Role' section contains three columns of checkboxes for various roles: 'Primary Contact', 'Other Contact', 'Birth Mother', 'Father', 'Legal Guardian', 'Kin', 'Foster Parent', 'Parent', 'Primary Care Provider', 'Attending Physician', 'Attending Nurse', 'EI Service Provider', 'Other Health Care Provider', and 'Medical Home'. At the bottom of the form are 'Save' and 'Cancel' buttons.

3. Enter the contacts information in the appropriate fields and select appropriate contact role types from the checkboxes below the demographics fields. Colored fields indicate mandatory data.

Select "Primary Contact" (in addition to any other roles you select) if this is the person that should receive mail and contact for this baby.

4. Click the "Save" button at the bottom of the screen

*To enter additional notes:*

1. Click on the "Notes" tab above the demographics screen or click on the yellow and blue Note preview box.

This baby's Note history is shown. If no notes are present a "New Entry" field is selected and the cursor is set for you to type.

2. Type in your note. The date and time are automatically stored with each note entry. In HiTrack 4 each note should be entered as a separate new entry unless you are making revisions to a prior note.

Custom fields are also entered on this screen. Use the "Type: " drop-down to select from your predefined list of custom fields.

3. Click the "Save" button at the bottom of the screen to retain your changes to Notes.

## Screening

The Screening link on the File Card is where hearing screening information may be reviewed and entered for a baby.

*To update screening results:*

1. Click on the "Screening" link from the File Card for this baby.

The Screening screen is shown.

**Screening Tests for Johns, Joseph**

New Screening Result

Stage:  Type:

Date:  Time:

Result: Right:  Left:  Source:

Screener:

Equipment ID:

Screening History - Select an item to Edit or Delete the result

Test Date	Type	Ear	Result	Screener	Facility	Stage
6/12/2006 12:00 AM	OAE	RIGHT	Pass	GATES, CYNTH	Mountain View M	Inpatient
6/12/2006 12:00 AM	OAE	LEFT	Pass	GATES, CYNTH	Mountain View M	Inpatient

Most Conclusive Inpatient Results				Most Conclusive Outpatient Results			
Method	Result	Type	Date	Method	Result	Type	Date
IP Right	Pass	OAE	6/12/2008				
IP Left	Pass	OAE	6/12/2006				
IP Head	Pass	OAE	6/12/2006				

This screen is separated into three areas. The top area lets you add and modify screening results. The middle area shows the screening history for this baby. The Bottom area shows the current Most Conclusive Results.

2. Add a new result by selecting from the appropriate fields in the top "New Screening Result" area. Colored fields are mandatory.

To edit an existing result click on the result in the "Screening History" list. This will activate the top area and allow you to change fields for this existing entry.

3. Click on "Apply" to complete the entry for this result.

The "Save" button will begin flashing to remind you to save your changes before leaving this screen.

4. Click the "Save" button at the bottom of the history list.

Under rare circumstances you or an audiologist may want to override HiTrack's Most Conclusive Result rules. Use the "Override MCRs" button to hand-select results to be used for head results and ear results for each stage.

## **Risk Indicators**

Selecting the Risk Indicators link will display a list of risk indicators for a baby. You can use the drop downs to indicate whether a baby has any of the risks shown here. A history of risk indicators may also be accessed from this screen.

## **Diagnostics**

The Diagnostics link provides a quick add screen for entering Hearing Disposition, DxABR, OAE, Tympanometry and Behavioral results all in one step.



Diagnostics details are for documentation purposes only. There are no reports in HiTrack that use diagnostic details as criteria. However, a filter for “Has Diagnostics” is available as an Advanced Find filter on the Baby List.

To use this feature first fill in the top section that asks for Audiologist, Test Date, Facility and Test Stage. Next, fill in diagnostic details for each of the areas in the bottom. If you do not use an area use the click arrow on the left to hide that section.

## **Diagnostics – Manage**

The Manage link displays histories for DxABR, OAE, Tympanometry, and Behavioral and also allows you to add new results or edit existing ones.

A baby’s Individual Diagnostic Report can also be created from this screen. This report summarizes this baby's complete HiTrack record in a clean, printable format suitable for sharing with supporting organizations.

## **Baby Letters**

The Letters link on the File Card shows the pending letters for selected baby and allows you to request a specific letter. This link also provides access to letter history for the selected baby.

For more information on letters see [Step 2: Generating Letters](#).

## Recommended Actions

The Recommendation system in HiTrack works together with the File Cabinet to automatically create "To Do" lists for you and your staff. For example, babies who do not pass Inpatient screening are moved to the "Need Outpatient Appointment" folder. As you enter scheduled recommendations for Outpatient Screening babies will be automatically moved to "Need Outpatient Screening" folder. When parents return with their baby for Outpatient screening their records can be easily located in the "Need Outpatient Screening" folder.

This scheduling feature also applies to the Need Evaluation folders.

HiTrack will automatically list suggested recommended actions on the File Card for babies in "Needing Appointment" folders. Clicking on one of these shortcuts takes you to an "Add Recommendation" screen.

*To enter a suggested recommendation*

1. Click on a suggested recommendation link under the Recommended Actions on the File Card for this baby.

Add Recommendation Screen is shown.

Recommend: Needs Outpatient Screening

Target

Date:  10/18/2012  Age: 1 Months

Confirmed

Date:    Age: 1 Months

Status: Scheduled

Notes:

Responsible Party

Person  Organization

Last Name:  Find Person

First Name:

Recommended By

Last Name:  Find Person

First Name:

On: 10/11/2012

Save Cancel

**Target Date:** The optimal date of completion for this action. For common recommendations the target date is prefilled for you.

**Confirmed Date:** The date of an actual appointment for this action.

Instead of Target Date and Confirmed Date, Target Age and Confirmed Age can be used. This option is chosen when a user wants to specify the optimal age of completion for an action.

**Status:** Pick a status for this action. "Referral Made" can be used for actions that are carried out by a supporting organization.

**Notes:** Explain briefly the reason or the details of this action.

**Responsible Party:** This is the person or organization that is responsible for carrying out this action. The "Find Person" button opens up a pop-up window that allows you to select a person by using filter criteria.

**Recommended By:** This is the person who requested that this action occur.

**Recommended On:** This is the date the recommendation was made.

2. Select appropriate values for each field.
3. Click "Save" at the bottom of this screen. A message may be displayed if this recommendation affects where the baby is listed in the File Cabinet folders.

*To see all recommendations or make a recommendation not listed*

1. Click on the "Manage" link under the Recommended Actions on the File Card for this baby.

A recommendations list is shown.

Action	Target Date	Confirmed Date	Status	Responsible Person	Notes
Needs Outpatient Scree	30 Jul 2007		Scheduled	Granger, Miranda -	

Click on "Add", "Edit" or "Delete" to work with the recommendation list.

2. Enter the recommendation details in the appropriate fields. (See steps above.)
3. Click "Save" at the bottom of this screen. A message may be displayed if this recommendation affects where the baby is listed in the File Cabinet folders.

### **Other Important Recommendations**

The "Early Intervention Referral" recommendation is important because it affects the [CDC Survey report](#). The "Total Referred for EI" count is based on the number of babies with this recommendation.

## Hearing Disposition

There are two tabs on the Hearing Disposition screen. The first of these allows you to enter hearing disposition data and select the most conclusive outcomes. The second tab is for recording audiograms. (The Web Companion shows the Audiogram page as a separate link on the File Card instead of a tab within the Hearing Disposition screen.)

The Hearing Disposition Most Conclusive Outcomes are an important tracking tool for babies receiving evaluation because they determine how this baby is seen on reports.

### *To record hearing disposition*

You must have Audiologists entered in for your facility to begin entering hearing disposition.

1. Click on the "Hearing Disposition" link on the File Card for this baby.

The Hearing Disposition Screen is shown.

Ear	Degree of Loss	Type of Loss	Date Eval'd	Confirmed Date	Confirmed	Audiologist
-----	----------------	--------------	-------------	----------------	-----------	-------------

Most Conclusive Outcomes

Right:

Left:

SF:

This screen is separated into two main areas. The top area lists all of the ear disposition entries for this baby. The lower area allows you to specify which entries to use as the Most Conclusive Outcomes for this baby.

2. Click on "Add" to add disposition for this baby.

The "Add Hearing Disposition" screen is shown with entries for Left, Right and Sound Field.

The screenshot displays the 'Add Hearing Disposition' interface with three main sections:

- Right Ear (highlighted with an orange border):**
  - Degree of Loss: Mild-Moderate
  - Type of Loss: Fluctuating Conductive
  - Configuration: Cookie Bite
  - Audiologist: ISUPONUS, DAWN
  - Date Evaluated: 9/17/2006
  - Date Confirmed:  9/17/2006
- Left Ear:**
  - Degree of Loss: [blank]
  - Type of Loss: N/A (No Loss)
  - Configuration: [blank]
  - Audiologist: EDWARDS, JEANNETTE
  - Date Evaluated: 9/17/2006
  - Date Confirmed:  [blank]
- Sound Field:**
  - Degree of Loss: [blank]
  - Type of Loss: N/A (No Loss)
  - Configuration: [blank]
  - Audiologist: EDWARDS, JEANNETTE
  - Date Evaluated: 9/17/2006
  - Date Confirmed:  [blank]

At the bottom right of the form are 'Save' and 'Cancel' buttons.

3. For each ear or sound field enter in appropriate values. HiTrack has been designed based on audiologic research to not allow incorrect combinations of Type and Degree of Loss. To record Normal Hearing use Degree: "None (Normal Hearing)".

It is important to enter the date of loss Confirmation so that HiTrack can correctly report this baby on Reports. Confer with your audiologists to decide when a loss is considered confirmed.

3. Click the "Save" button at the bottom of the screen.
4. After saving HiTrack requires that you select Most Conclusive Outcomes. This is especially important for babies with multiple historical disposition entries.

Pick a disposition entry for each ear and click "Save Changes" to retain your Outcome choices.



If you report to the CDC please note that the CDC Survey reports both on the first and Most Conclusive Outcomes selected in Hearing Disposition. Infants that have confirmed initial entries, but unconfirmed most conclusive outcomes may adversely affect the numbers generated on HiTrack's CDC Report. Please contact the [Help Desk](#) for more information about how to check for these scenarios.

## Amplification

The Amplification link is where information may be entered regarding a baby's hearing devices.

## EI Services

There are three data entry screens in HiTrack that affect Early Intervention Tracking and [Reporting](#): EI Referral data entry, EI Enrollment and Status Update, and Stop Tracking.

To enter the EI Referral start by selecting "Manage" from the Recommended Actions section on the baby File Card. Next choose "Add" from the bottom of the page.

The screenshot shows a web form for entering an Early Intervention Referral. The form is divided into several sections:

- Recommend:** A dropdown menu set to "Early Intervention Referral".
- Target:** A section with a "Date" field (checkbox), an "Age" dropdown set to "1", and the word "Months".
- Confirmed:** A section with a "Date" field (checkbox) set to "1/ 2/2013", an "Age" dropdown set to "1", and the word "Months".
- Status:** A dropdown menu set to "Scheduled".
- Notes:** A large text input area.
- Responsible Party:** A section with radio buttons for "Person" and "Organization" (selected), and a "Name" dropdown menu set to "Creekside Hospital".
- Recommended By:** A section with "Last Name:" and "First Name:" input fields, a "Find Person" button, and an "On:" date dropdown menu set to "10/26/2012".
- Buttons:** "Save" and "Cancel" buttons at the bottom.

Callout boxes on the right provide instructions for several fields:

- Line 1: Select "Early Intervention Referral" as the type of recommendation.
- Line 2: Select an appropriate date in the future for when this should occur.
- Line 3: Select the Person or Organization to whom the referral is being made.
- Line 4: Select the person who is making the Referral.
- Line 5: Select the date the Referral was made. **This is used to calculate the Age at Referral for the CDC Survey.** 

To enter Enrollment details, including cases of Parent Opt Out, Out of Jurisdiction or Monitoring services, use the "EI Services" link on the baby File Card. From this screen use the "Add" button at the bottom to make a new entry.


The screenshot shows a form for entering enrollment details. The fields are as follows:

- Status:** A dropdown menu currently set to "Enrolled".
- Enrolled With:** A dropdown menu currently set to "Falk, Trudy - Rainbow Marathon".
- Date:** A dropdown menu currently set to "10/26/2012".
- Age:** A text field currently set to "4 months".
- Service Type:** A dropdown menu currently set to "Service Coordination".
- Activities:** A text area with a "500 chars max" limit.
- Private Enrollment (Non-Part C):** An unchecked checkbox.
- Buttons:** "Save" and "Cancel" buttons at the bottom.

Three callout boxes provide instructions:

- Top box: "Select the appropriate status. For babies that were not Enrolled choose an alternate status." (points to the Status dropdown)
- Middle box: "Select the program and date of Enrollment or non-enrollment status." (points to the Enrolled With and Date dropdowns)
- Bottom box: "Use this checkbox to indicate Non-Part C enrollment." (points to the Private Enrollment checkbox)

After saving a baby's enrollment entry you should indicate their eligibility by selecting an option from the drop down at the bottom of the EI Services screen. Yes, No and Unknown are the available choices.



To indicate that EI Follow-up should be discontinued on a baby, use the Tracking Options link on the File Card and then choose Stop Tracking. For this status to be counted as intervention-related the baby must be in the "Referred for Intervention" File Cabinet folder.

Medical ID: 0349854  
Birth Date: 6/11/2012  
Action: Referred for Intervention  
Status: In Progress

Record Overview

- Tracking Options
- Demographics
- Screening
- Risk Indicators
- Diagnostics | Manage
- Hearing Disposition
- Amplification
- EI Services

Letters  
Transfer

Recommended Actions:  
Manage

Messages:  
EVALUATION IN PROCESS.  
10/26/12  
Milestones: 1/-/-

Current Status: In Progress

Status History

From	To	During	Date
------	----	--------	------

Change Status to: Follow-Up Discontinue Date: 9/ 2/2012 Save

If the baby you are working on is not already in the "Referred for Intervention" folder, you will need to choose an Action from the "Actions" tab that will put them in this folder before stopping their tracking. For example if the record is currently in "Need Evaluation Appt" you would pick the "Refer for Intervention" Action Needed folder. For other folders you will need to first uncheck the "Suggested Actions" filter box to see more folder options. Contact the [Help Desk](#) for additional support information on this feature.

## Transfer

The Transfer feature works in one of two ways depending on your database layout. When using a Centralized database (i.e. Web Companion) the "Transferred" Baby list folder includes a read only view of babies that have transferred out of your Facility that are being listed and reported under a different facility. The Centralized Database setting can be changed in Advanced Settings.

If you access HiTrack using the Windows Client at your Facility and transmit data to a State or Regional EHDI coordinator you are typically configured with a Distributed database layout. In this case the records marked as being transferred are included in your periodic State/Regional transfer. These babies are shown in the "Transferred" Baby List folder and are marked as "Transferred Out" because they are not being tracked in your portion of the Distributed database.

The transfer link can be used to transfer a baby to another facility within the system. An option is also available to transfer the baby out of the system.

As [described earlier](#), the Centralized setting determines how transfer babies are handled in HiTrack, especially in regard to which facilities have access to the record after the transfer is entered.

Transfer history entries can be added, edited and deleted. HiTrack also allows you to determine the Responsible Facility transferred out of your Facility. If changes to Transfer history suggest that the baby should have a different Responsible facility, HiTrack will pre-select the Responsible facility drop-down and indicate that the user to save the change.

Responsible Facility from Baby records coming from HiTrack 3.5 (including merged ETF files) default to the baby Birth Facility unless there is a HiTrack 3.5 Transfer record to a different facility.

Note: Responsible Facility is not used on reports or any other tracking features.

## Tracking Tasks - Step 2: Generating Letters

Letters should be generated when you have finished reviewing the screening records. All baby records that need further screening or evaluation should be finished using Quick Add or Tracking Options. Contact information should also be provided for these records.

*To generate letters:*

1. Click on the "Letters" postage stamp icon from the Toolbar  
The Facility/Facility Group pending letter list tab will be shown.  
The baby, recipient, letter title and current screening outcomes are shown in columns on this list. Clicking on a column header will sort the list.
2. From this list you can review which letters will be generated.  
For Letters you do not wish to print select the "Never" print option.  
Use the shortcut buttons at the bottom of the screen to reset the print options for the whole list.
3. Click the "Process" button in the lower right when you are ready to continue.  
Your word processor will open this batch of Letters for your review.  
Use "File" and "Print" to send the letters to your printer.
4. After printing close your word processor and return to HiTrack.  
HiTrack will prompt you asking "Did the letters print successfully?"
5. If there was a paper jam or some other printing problem click "No". This leaves the letters in the list so you can reprocess them when you are ready.

A Letters Reload button and date entry gives you the opportunity to override the date logic that detects new possible for letters. In most cases this feature will not be needed; new letters will automatically be shown on the list in Step 1 above. Contact the HiTrack Help Desk if you have questions about Letter criteria.

The Letter History tab provides a list of all of the letters sent by HiTrack. Use the date filters at the top of the screen to help narrow down the list of letters. Use the Resend button to send an exact copy of the letters that were originally sent. To manually generate a new letter with updated content use the Letter Request feature from the Letters link on a baby's File Card.

The Letter History tab will also indicate the format of historical letters. A HiTrack database may have past letters in both the original Rich Text Format and Microsoft's DOCX format for the current Facility/Facility Group. Two letter documents will be generated if letters of different formats are selected to be resent.



In some cases the letter text will be missing from letters in Letter History. This scenario happens for patient records that have been merged in from older HiTrack databases or when a letter is removed from the Pending Letter list without being processed. New letter content is generated for that baby's letter if these letters are selected to be resent.

Pending Letters and Letter History for individual babies can also be accessed by using the Letters link on the baby File Card.

### Tracking Tasks - Step 3: Generating Reports

Program Coordinators should periodically use the HiTrack Reports to review the performance of their programs.

HiTrack's Reports are designed to provide statistics for Single and Two Stage protocols. The 'No Out Scrn ->Dx' ("Go directly to Diagnostics without Outpatient Screening") numbers account for non-Pass babies in Single Stage protocols on the Outpatient section of reports.

### General Filtering Options

Most reports in HiTrack can be filtered using the following criteria:

Date	Reports can be filtered by Birth date, Screening date or Enter Date.  "Enter Date" is the date when the baby record was hand entered or merged into the HiTrack database. This date does not change and is not the same as the last change date.  Open ended birth ranges can be used by un-checking the checkboxes to the left of the date entry boxes.
Location	Birth Facility, Inpatient Screening Site, Nursery are selectable as location filters. Facility and Nursery Groups filters can be used to combine different sets of infants for reporting.
Technology	Many of the reports allow filtering by the Technology used in screening. This allows program coordinators to show babies who were screened with only a specific technology or who were screened with both.  This helps answer the question: "How can I report screening effectiveness by technology?" 

These three filter options can be used in any combination to create a targeted report.



HiTrack will remember each report's latest criteria so remember to review the previously used filters when running reporting reports for different purposes.

Contact information can be included on reports that provide lists of babies.

Each report has its own relevant filter options in addition to the general filter options. [Appendix B](#) describes these options in detail for each report.

## Available Reports

The "Reports" icon from the Toolbar in HiTrack 4 provides access to the following reports:

Click on a report link listed below for detailed report usage and criteria.

<a href="#">Amplification Report</a>	Report amplification history for babies
<a href="#">CDC Survey Report</a>	A report providing CDC survey data
<a href="#">Data Cleaning Reports</a>	Generate reports to help in record consolidation
Duplicate Medical ID	Report babies with matching Medical ID
Probable Duplicate	Report babies with matching DOB and last name
Custom Duplicate	Report babies on a variety of matching field combinations
<a href="#">Early Intervention Report</a>	Report Early Intervention information for babies.
<a href="#">Hearing Status Report</a>	Report babies whose hearing disposition has been evaluated.
<a href="#">Needing Evaluation Report</a>	Report babies who need continued evaluation.
<a href="#">Needing Outpatient Report</a>	Report babies that need Outpatient Screening.
<a href="#">Outpatient Screening Report</a>	Report babies with Outpatient results.
<a href="#">Physicians Report</a>	Report babies grouped by Primary Care Provider
<a href="#">Risk Indicator Report</a>	Report babies with risk indicators.
<a href="#">Screener Profile Report</a>	Report screener response rates
<a href="#">Screening Results Report</a>	Report screening results for babies.
<a href="#">Flow Chart</a>	Generates a visual flow of screening and evaluation
<a href="#">Hospital Summary</a>	Reports screening and evaluation statistics by hospital
<a href="#">Quarterly Report</a>	Reports screening and evaluation statistics by quarter for hospitals
<a href="#">Technology Summary</a>	Reports screening outcomes by equipment technology
<a href="#">Milestone Summary</a>	Reports program timeliness for the CDC 1/3/6 timelines
<a href="#">User Defined Reports</a>	Create printable baby lists with user-specified columns and criteria

*To generate a program Summary report:*

1. Click on the "Reports" icon from the Toolbar  
The available reports will be shown in a list.
2. Select the "Summary Reports" entry from the list of reports.  
The report selection criteria screen is shown.  
Select the type of Summary report using the small list in the upper left.
3. Enter your filter criteria
4. Select a grouping option.  
Grouping the data as a single report shows all babies for the selection criteria on one report page. When grouped by Birth facility or IP Screening Site the total list of babies are divided into different counts based on where they were born or where they were screened. Note that changing the grouping options does not change the total number of babies included; it only affects how they are shown.
5. Click the "Generate" button.  
A progress bar is shown until the report is generated.  
After generating a new window is opened containing a print preview of the report.
6. To review the report page use the left and right wedge arrows on the toolbar.
7. When you are done reviewing and printing you may close the report with the X in the top right of the window.  
Report windows may be left open for side by side comparison.

*To generate a baby listing report:*

1. Click on the "Reports" icon from the Toolbar  
The available reports will be shown in a list.
2. Select a report from the list (excluding the "Summary Reports" entry) and click the "Select" button below.  
The report selection criteria screen is shown.  
In most cases the general criteria is listed in a column on the left and the report specific criteria is listed on the right.
3. Specify report criteria and click the "Generate" button.  
After generating a new window is opened containing a print preview of the report.
4. To review the report page use the left and right wedge arrows on the toolbar.
5. When you are done reviewing and printing you may close the report with the X in the top right of the window.

## Report Advanced Criteria

Most Reports in HiTrack can be filtered using the [Advanced Criteria](#) feature.

## Printing and Exporting Reports

Once a report preview window has been opened you print the report or export it in the following formats:

Crystal Reports (*.rpt)	A file that is Compatible with Crystal Report Designer tools
Adobe Acrobat (*.pdf)	A printable format that maintains page layout.
Microsoft Excel (*.xls)	An Excel format that attempts to maintain similar page formatting
Microsoft Excel Data Only (*.xls)	An Excel format that attempts to provide the data for spreadsheet processing. Empty rows and columns are included for report white space.
Microsoft Word (*.doc)	A Word format that attempts to maintain similar page formatting. Most text fields are converted to "frame" objects.
Rich Text Format (*.rtf)	A Word compatible RTF file that attempts to maintain similar page formatting. Report fields are shown in a single vertical column when opened in WordPad.

Please contact the [HiTrack Help Desk](#) if you would like help exporting a report in a different format.

## Milestones

The Milestone feature enables you to easily identify infants who are nearing or who have missed the CDC recommended dates for hearing screening (before 1 month), loss identification (before 3 months), and intervention (before 6 months).

Two features in HiTrack give you Milestone information, the Milestone Baby List and the Milestone Summary report.

When the Baby List Milestone feature is turned on, three new columns appear on the Baby List, Milestone, Status, and Target Date. In addition, a drop-down box appears directly to the right of the Milestones checkbox. The drop-down box allows you to select which milestone category to display. After selecting the milestone category, click the Refresh button to view the results in the baby list.

The screenshot shows the 'All Data for Green River Regional (15 of 15)' interface. It features a filter section with 'DOB From' (4/ 1/2010) and 'To' (10/ 9/2012) dropdowns, 'Advanced Find' and 'Milestones' checkboxes, and a 'Current Milestones' dropdown. Below this are input fields for 'Medical ID', 'Last Name', and 'First Name', a 'Birth Facility' dropdown, and 'Appointments Within' fields. A 'Result' dropdown is set to 'All', and 'Records' is set to 'All'. The main table has columns: Milestone, Standing, Target Date, Status, Baby Name, and Birth Facility. A filter dropdown is open over the 'Baby Name' column, showing 'Complete' and 'In Progress' options. Callouts identify 'Advanced Filter Criteria', 'Milestone Categories', 'Milestone Columns', and 'Auto Filters'.

Milestone	Standing	Target Date	Status	Baby Name	Birth Facility
ID by 3	Due soon	1/6/2013	In Progress	(Clear Filter)	Green River Regional
ID by 3	Met	1/4/2013	In Progress	Complete	Green River Regional
ID by 3	Due soon	1/3/2013	In Progress	In Progress	Green River Regional
Screen by 1	Due soon	11/2/2012	In Progress	Rose, Kevin	Green River Regional
ID by 3	Met	12/30/2012	In Progress	Stevens, Rebecca	Green River Regional
ID by 3	Met	12/30/2012	In Progress	Willard, Kevin	Green River Regional
ID by 3	Due soon	12/30/2012	In Progress	Smith, James	Green River Regional

The [Milestone Summary Report](#) is another powerful tracking tool. Please refer to the Report Criteria section in this document for information about this report.

## Tracking Tasks - Step 4: Transferring Data to a Central Database



This step creates an Encrypted Transfer File (ETF4) file for the purpose of reporting your HiTrack data to a state or regional Central coordinator. *This step is not required if you a) do not report your data to a Central coordinator or b) use the HiTrack Web Companion.*


This step requires you to install your state or region's Public Key. [Please refer to the Tools section for instructions for this one-time installation.](#)

The Standard transfer option creates a file containing new and changed records. Standard transfers should be created and sent to your Central coordinator on a regular basis.

*To create a Standard transfer:*

1. Select your individual facility from the Toolbar
2. Click on the "Tools" icon from the Toolbar. The Tools tabs will be shown.
3. Click on the "State/Regional Transfer" tab.

The Data Transfer Options sub-tab will be shown.



This tab shows the folder where the transfer file will be created and when the last transfer was made.

The number of records waiting to be transferred is shown (If this is your first data transfer all of your HiTrack data will be transferred.) This number is the total number of records including babies, medical persons and other Program Information lists.

If your Hospital uses a "Two-Stage Manual" or "Two-Stage Automated Pass/Manual Non-Pass" protocol, some infants may be listed under "Not sending". To include these infants in your Transfer, go to the Action Needed – Inpatient Screening folder and update the infants Tracking History. For assistance, go to Tracking Options page 21.

4. Click on the "Transfer" button. A success message will be shown after the ETF4 file is created.
5. Use your Email program to create a new message to your Central coordinator with this ETF4 file as an attachment. Alternatively, you could use the [HiTrack Automatic Email](#) found in User Settings to generate this email for you.

 You can open the folder where transfers are being made by double clicking in the box under "Transfer File to:" (This works with any folder field in HiTrack.)

### Types of Data Transfers

As mentioned above Standard transfers are the routine way of sending data to your Central coordinator. Periodic Standard data transfers keep your Central coordinators database up to date.

The Maintenance Transfer feature should generally be used at the direction of your Central coordinator or the HiTrack Help Desk. Maintenance transfers are needed to resend missing data, most often due to human error.



There are four types of HITACK data transfers; one routine and three are for maintenance:

<p><b>Standard Transfer</b> – A transfer file is created that includes new and changed baby records. (described above in Tracking Task 4)</p> <p>Standard transfers should be created and sent to your Central coordinator on a regular basis.</p>
<p><b>Transfer by Birth Date</b> – A transfer file is created that includes baby records with birth dates between two dates. All of the data fields in each baby's record are included in the transfer.</p> <p> <i>Use this transfer when you know which babies you want to resend to your Central coordinator.</i></p>
<p><b>Transfer by Change Date</b> – A transfer file is created that includes baby records that have changed between two dates. Only parts of the record that changed between these two dates are included in the transfer. (Demographics are always included.)</p> <p> <i>Use this transfer to resend a lost period of changes independent of how old the babies are.</i></p>
<p><b>Transfer entire current facility</b>– A transfer file containing all of your Facility records is created. Every part of each baby's record is included in the transfer.</p>
<p><b>Baby Basket Transfer</b>– A transfer file containing all of your babies' records in your Baby Basket is created.</p>

The following example explains the differences between the Birth Date and Change Date transfers:

Three babies were born at your hospital during the month. Baby "A" was born on the 6<sup>th</sup>, Baby "B" was born on the 15<sup>th</sup> and Baby "C" was born on the 17<sup>th</sup>.

**Birth Date Calendar**

			1	2	3	4
5	6 Baby A Born & Entered	7	8	9	10	11
12	13	14	15 Baby B Born & Entered	16	17 Baby C Born & Entered	18
19	20	21	22	23	24	25
26	27	28	29			

Three baby records will be created in HiTrack as their information is manually entered on their birth date. A Birth Date transfer with a date range from the 6<sup>th</sup> to the 15<sup>th</sup> would include Baby "A" and Baby "B" but not Baby "C". As explained above all of the data fields for Baby "A" and Baby "B" will be included in the transfer file.

After entering these babies into HiTrack updates were made to their records. These changes began on the 21<sup>st</sup> and were finished on the 28<sup>th</sup>.

**Record Change Date Calendar**

			1	2	3	4
5	6 Baby A Born	7	8	9	10	11
12	13	14	15 Baby B Born	16	17 Baby C Born	18
19	20	21 Baby A Edited	22	23 Baby B Edited	24	25
26	27	28 Baby C Edited	29			

A Change Date transfer with a date range from the 20<sup>th</sup> to the 24<sup>th</sup> would include data fields that changed for Baby "A" and Baby "B" in the transfer. Baby "C" is not included because their last change occurred on the 28<sup>th</sup>.



Note that a Change Date transfer using the prior date range of the 6<sup>th</sup> to the 15<sup>th</sup> would not include any records because no records were changed during that period of time.

The biggest difference between the Birth Date and Change Date transfers is that Change Date transfers include records based on the when the edits and updates were made in HiTrack and Birth Date transfers include records based on the date of birth of the baby.

### Remake/View Transfer History

The best way to re-create a lost transfer file is to use the "Remake/View Transfer history" feature. After clicking on this button a historical list of data transfers will be shown. To remake a transfer click on the row that you would like to remake. The Data Transfer options screen will be returned and the date range criteria for a maintenance transfer will be preselected for you.

### Transfer File Name

ETF4 transfer files use the following naming format:

The first 10 Letters of your Facility name (underlines for spaces) & ...
4 Characters of the HiTrack ID for that Facility (This is part of an internal unique identifier and is not the same as HiTrack 3.5 "State Transfer ID") & ...
Date of Transfer: Day, Month (Three letter abbreviation), Year (Two Digit) & ...
A Three or Two letter "Type of Transfer" description: --"STD" is Standard Transfer, (New and changed records) --"DB" is Date of Birth Range Maintenance Transfer --"CD" is Change Date Range Maintenance Transfer --"EN" is Entire Facility Maintenance Transfer --"EP" is Program Information export

Be aware that subsequent transfer files will be named identically if you create more than one transfer of a specific type on a single day. If you use email to send a maintenance transfer to your Central coordinator be sure to explain the type of transfer and the date range used. Otherwise, they may skip your file because it is named the same as a previous maintenance transfer.



## **Maintaining the HiTrack Database**

HiTrack provides several tools to help you maintain your database. The Tool box icon opens the Tools screen with several tabs along the top. A brief description of each tab is given below.

<b>Program Information</b> – HiTrack program information can be exported to a specified file to be imported by screening equipment software.
<b>Change Logs and Audits</b> -- Allows users to view changes made to baby data after a selected start date.
<b>State/Regional Transfer</b> – Provides an encrypted data transfer file for use in another HiTrack system at a State or Regional office.
<b>Back-ups/Conversion</b> – Backup or restore a HiTrack database. Convert existing data from HiTrack 3.5 Desktop or Enterprise editions.
<b>Key Maintenance</b> – Install encryption keys for State/Regional transfer data and other encryption key Tools.
<b>Advanced Tools</b> – Other tools for maintaining your database.

### **Program Information**

This tool is used to export your Facility's lists of birth locations, physicians, audiologists and screeners for use in [screening equipment software](#) when HiTrack is not installed on the screening equipment computer. (If HiTrack's Windows user interface is installed on the same computer as your screening software you will not need this tool.)

To export your lists specify a folder or drive and then click the "Export Program Info" button. This creates a series of .TXT and .DBF files used by screening equipment software. An XML formatted version of your Program Information is also exported.

### **Change Log and Audits**

This feature allows you to view the changes made to baby and person records at your Facility. Use the checkbox and date selector to filter what changes are shown. The Change Log can be browsed on-screen or exported to a file for future analysis.

Read [Tracking Tasks – Step 1](#) for more about the fields maintained in the Change Log.

### **State/Regional Transfer**

This tool is outlined in the [Tracking Tasks - Step 4](#).

## Back-ups and Conversion

Backing up data is one of the most important elements of good data management. Although the HiTrack Back-up feature is not a substitute for a permanent back-up system, this utility does provide a quick and convenient way for you to copy infant records and critical program information to a diskette or zip drive. If you do not have another back-up system, take advantage of this feature to back up your data on at least a weekly basis.

Backups can only be made using HiTrack if you are using the computer that is serving the database. Backup features are disabled if your computer is only installed with the application component of HiTrack.

*To make a backup of your HiTrack database:*

1. Click on the "Launch Database Manager" button from the "Backups/Conversion" tab.  
The Database Manager will be launched. Three buttons will be shown: "Backup Database", "Restore Database" and "Update database".
2. Click on the "Backup Database" button.  
A Browse for Folder window will be shown. Choose a folder where the backup should go and then click "OK".
3. You will then be asked for an optional password. This allows you to protect your database file. If you do not want to put a password on the backup you may simply click Proceed. Forgotten passwords on back up files cannot be retrieved.

After processing a "Database Backup Completed Successfully" message will appear.

The name of the backup file created will start with "HT4Backup\_" and then be followed by the date and time the backup was created. It will have a .zip extension because it was compressed to save file space.



HiTrack also offers an automated backup feature. This will back up your HiTrack database at a specified time or event. To use this feature, open up the Database Manger and click on the "Add Auto Backup" button. Now select a folder where your HiTrack backups will be placed and then complete the rest of the backup wizard.

Upon completion of the wizard, the "Add Auto Backup" button reads "Edit Auto Backup". If you click on this button, you can change the settings for your automated backup of HiTrack. If at any time you want to stop doing automated backups, click on the "Delete Auto Backup" button and your HiTrack database will no longer be automatically backed up.

*To restore an existing HiTrack 4 backup:*

**THESE STEPS WILL OVERWRITE YOUR CURRENT HiTrack DATABASE**

1. Click on the "Launch Database Manager" button from the "Backups/Conversion" tab.  
The Database Manager will be launched. Three buttons will be shown: "Backup Database", "Restore Database" and "Update database".
2. Click on the "Restore Database" button.  
A small Restore Database window will be shown.

Use the small folder icon to the right to locate a HiTrack 4 backup file.

3. Click on the "Proceed" button to restore the data and overwrite your current HiTrack database.

The Database Manager can also be access by clicking on the "plug-in" icon from HiTrack's opening page.

Please contact the [HiTrack Help Desk](#) for Conversion support.

## Key Maintenance

This tool helps you maintain your public and private keys used for data encryption in HiTrack.

If you report your HiTrack data to a state or regional coordinator you will need to install their Public Key. In most cases this is provided to you as an email attachment called "PublicKey.pbk". Use your email program to save this attachment to your computers hard drive and then use this screen in HiTrack to install it. Use the "Browse" button on the left side of the screen to locate the PublicKey.pbk file.



(Refer to [Tracking Tasks - Step 4](#) for steps on creating a Standard data transfer.)

If you are a state or regional Central HiTrack user you will need to create a Decryption Key Set and then export your Public Key for the users who transfer data to you. Use the buttons on the right side of this screen to perform these tasks.

After exporting the Public Key it can be safely emailed to all of your users. Public Keys can be shared over non-secure networks because they can only be used to encrypt data files. *Private keys are required to decrypt data files and should not be shared over non-secure networks.*

## Advanced Tools

Use Advanced Tools under the direction of the HiTrack Help Desk.

### Classification Browser Tool

The Classification Browser tool in HiTrack allows program coordinators to browse baby records by creating combinations of [Classifications](#). This tool is available by clicking on the Tools link in the tool bar and then clicking on the Classification Browser button at the bottom of the Advanced Tools tab.

The pick lists on the left allow for filtering by stage specific classifications. The "NOT" checkboxes provide the ability to exclude the selected classification instead of including it. Any combinations of classifications can be used across the four stages. A birth date filter is also provided to narrow down the results.

The Refresh button queries the database and reloads the list according to the selected classification filters. After the list loads the number of baby records is reported and highlighted in a box to the right.



The Classification browser helps answer the following questions: "How do I see a list of infants show on the Flow Chart or other Summary report?" and "How can I improve my program by accounting for infants with atypical diagnosis?"



The Classification Browser grid includes a row for each baby's classification in the list. This shows how the baby is being reported across all four stages. This list can be sorted and exported for use outside of HiTrack. The baby's internal HiTrack ID number is also provided in this list.

## Data Integration



HiTrack is designed to integrate with other database systems in several ways. Incoming data is imported through HiTrack's advanced Merge system so that it can be checked for completeness and existing records can be are matched based on rule-based algorithms. It also means that questionable matches are left for manual review.

The following methods are available for populating and querying HiTrack:

<b>File Exchange (XML, TXT, CSV)</b> – HiTrack can import data from XML and Comma Separated Value formatted files. HiTrack data can also be exported in these formats.
--

<b>Stored Procedures (SQL)</b> – This method allows batches of data to come into HiTrack using HiTrack's Merge Database Stored procedures. HiTrack also provides stored procedures for other systems to query.
--

<b>Web Services (SOAP)</b> – HiTrack provides several industry standard Web Services that can be used to add or query data using web protocols.
---

<b>Electronic Messaging (HL7)</b> – HiTrack provides HL7 V2 compatibility using integration software called Mirth Connect. Both incoming and outgoing data channels are available to use with HiTrack data.
---

Please contact the [HiTrack Help Desk](#) to discuss data integration installation and configuration.

## ***Helpful Hints for Data Entry in HiTrack***

This section contains data entry shortcuts to make entering data quick and easy.

### **Navigating screens and fields**

Keyboard shortcuts for buttons are shown on the Windows user interface with an underline under the shortcut letter. For example if a button reads "Cancel" you can press Alt+C on the keyboard to cancel the screen.

To move through the data entry fields, click the field, press Tab to move forward or press Shift+Tab to move backward.

HiTrack Windows now supports additional standard Keyboard shortcuts, including:

- ESC ..... Leave the screen you are on
- ENTER .. Save the changes you have made on this screen
- F1 ..... Shows support options
- F2 ..... Opens the Action Needed Cabinet
- F3 ..... Activates Advanced Find feature
- F4 ..... Shows/Hides the File Card
- F5 ..... Reloads the Baby List
- F6 ..... Shows/Hides the Baby Basket
- F7 ..... Highlights the Facility pick list on the Tool bar
- F8 ..... Opens the Letters Screen
- F9 ..... Opens the Reports Screen
- F10 ..... Opens the Merge Screen
- F11 ..... Opens the Tools Screen
- F12 ..... Opens the Settings Screen

HiTrack Web now supports access keys. Depending on the browser, they can be used by holding down Alt and pressing the access key, or by holding down Alt+Shift and pressing the access key. The supported access keys are:

- 1 ..... Shows support options
- 2..... Opens the Action Needed Cabinet
- 3 ..... Activates Advanced Find feature
- 4 ..... Shows/Hides the File Card
- 5 ..... Reloads the Baby List
- 6 ..... Shows/Hides the Baby Basket
- 7 ..... Shows a list of the appointments happening today
- 8 ..... Opens the Letters Screen
- 9 ..... Opens the Reports Screen
- 0 ..... Opens the Merge Screen
- - ..... Opens the Tools Screen

- = ..... Opens the Settings Screen
- B ..... Navigates back to the Baby List

### Opening Records

In many case you can open a record for review and editing by double-clicking the mouse on the row of information you are interested in.

### Entering Dates

Each date entry control allows you to enter the month, day and years using the keyboard or a pop-up calendar. Clicking on the icon to the side of the date field will open the pop-up calendar.

In some cases dates are optional on several screens, including report filters. Optional dates will have a check box to the left of the date. Uncheck the box to exclude a date or leave it blank.

### Using drop-down menus

Drop-down menus offer preset options that save time and standardize data entry. You can also use the keyboard to help you select items from drop-down menus. When a drop-down is selected you can jump to existing presets by typing the first letter of the entry.

### Entering File Locations

Each Windows user interface screen that requires you to fill in a computer file location includes a folder icon. Clicking this folder icon will help you select a location from your computer or local network. It also allows you to create a new folder and select it in one step. The folder shown in any file location box can also be automatically opened by double-clicking in white area of box.

### Data Grids

The contents of most data grids in the Windows user interface may be printed or exported by right-clicking in the grey area of a grid. After right-clicking you will have these options:

<b>Copy Selected Row</b> – Copies the selected row into the Windows clipboard. You can then "Paste" this entry as text into another program.
<b>Print</b> – Opens a Print Preview screen for the grid you have selected.
<b>Save as XML</b> – Allows you to save this grid as an "XML" file with accurate database structure representation.
<b>Save as CSV</b> – Allows you to save this grid as a "Comma Separated Values" file that is compatible with Microsoft Excel.

## HiTrack's Web Companion

This section describes HiTrack's Web Companion. The Web Companion is a web-based user interface designed to operate within the modern versions of Internet Explorer or Firefox web browsers.

As [described above](#) your HiTrack database may be configured with a Distributed or Centralized layout. In most cases the Web Companion provides access to a Centralized state or regional database. Please contact the [HiTrack Help Desk](#) with Web Companion usage questions or for instructions on how to install the Web Companion IIS Server application (IT required).

To use the Web Companion you will need to open your web browser to a predetermined web site provided by your Central coordinator (i.e., "https://hitrack.state.gov/"). This page will open the HiTrack login page. You should see a yellow security padlock in your web browser when this page loads. (If you do not see a yellow security padlock please contact your Central coordinator before using HiTrack.)

### Web Browser Settings

To use the Web Companion the following settings must be configured in your web browser:

<b>Internet Explorer</b> 7 or newer (recommended), Firefox 2 or newer. Other web browsers (i.e. Opera, Netscape) are NOT supported.
<b>Trusted Sites</b> – The HITRACK web application site must be added to the Internet Explorer "Trusted Sites" list and "Automatic prompting for file downloads" must be enabled for Trusted Sites.  To do this in Internet Explorer go to the HiTrack Web site and select Tools and then Internet Options. From here select the Trusted Site icon from the Security tab. After selecting the icon click the Sites button to open the Trusted Sites list that allows you to add the HiTrack site.  To allow prompting for downloads click on the Custom Level... button for the Trusted Sites zone. Enable automatic prompting by filling in the radio button under the Downloads section.
<b>JavaScript</b> must be enabled for the HITRACK web application site.
<b>Popup windows</b> must be allowed for the HITRACK web application site.  Most browsers provide feedback if a pop-up is blocked. The option to "Always allow for this site" should be used.
<b>A PDF viewer</b> (i.e. Adobe Reader) if using the Reports feature.
<b>Microsoft Office Word</b> 2003 (or later) is required for Letters
<b>Microsoft .Net Framework 2.0</b> is required on client computers where users choose to use the Merge Companion upload tool.  The Merge Companion provides automated upload features and renames screening equipment merge files as they are processed.

## Differences in Web Companion

A few main differences between HiTrack's Windows user interface and the Web Companion are:

<p><b>Automatic Session Timeout</b> –Web Companion is a web application that can be used by hundreds of users. To conserve computing resources the system provides a built-in session timeout feature. This feature (provided by the web server system) will end your session if it is unattended for a certain amount of elapsed time. Typically this time period is set to 20 minutes of inactivity.</p>
<p><b>Software Installation and Maintenance</b> – The Web Companion server software can only be installed and maintained with IT support.</p>
<p><b>Letter Templates</b> – Letter Templates for your Facility or Facility Group cannot be edited using the Web Companion. If you need to customize the templates for your facility you can submit your suggested letter text to your Central coordinator or install the Windows user interface on one computer at your site to access the template editor. Please call the <a href="#">HiTrack Help Desk</a> for letter template support.</p>

Other functional differences between the Windows user interface and the Web Companion user interface include how data is imported and exported.


<p><b>Collecting Screening Information from Screening Equipment Software</b> –</p> <p>Web browser software applications such as the HiTrack Web Companion are not allowed to access files on your computer, including merge equipment files. To account for this constraint the Web Companion provides both a <b>manual</b> and an <b>automatic</b> upload method for merging equipment files. These two options are described in the next section.</p>
<p><b>Program Information</b> – Instead of exporting separate program information files to a given folder the Web Companion provides a download link for a single ZIP file that contains all of the Program Information files for your facility. After downloading, the contents of this Compressed Zipped file can be extracted <a href="#">for your screening equipment</a>.</p>
<p><b>Processing Letters</b> – After clicking "<a href="#">Process</a>" on the pending Letter list an RTF word processing file containing your letters will be downloaded from the Web Companion site. In some cases your web browser will automatically show this file, in other cases you will need to download the file and open it manually for printing.</p>

## Manually Uploading Equipment Files to the Web Companion

One advantage to the manual upload method is that it does not require any additional tools or downloads. The main **disadvantage** is that you must manually rename or move the equipment file on your computer after the merge completes.

*To manually upload and merge an equipment file:*

1. Click on the Merge icon from the toolbar. The Merge screen is shown with the Manual Upload tab selected by default.



The screenshot shows a web interface with two tabs: "Manual Upload" (selected) and "Automated Upload". Below the tabs, there are three steps:


- Step #1: Select a source.** A dropdown menu is open, showing "HT4 XML Data File 2011".
- Step #2: Use 'Browse' to find a merge file.** A "Choose File" button is highlighted, with "No file chosen" text to its right. Below it is an "Upload File" button.
- Step #3: Click 'Merge Data' to Merge this file.** A "Merge Data" button is highlighted.

Select your equipment [Merge Source](#) from the drop down in Step #1.

- 2a. Click the "Choose File" button. A "Choose File" pop-up window will be shown. Locate the equipment file on your computer and then click Open.

The equipment file path will be shown in the small field to the left of the Browse button.

- 2b. Click the "Upload File" button to send a copy of this file to the server computer where the Web Companion site runs. This screen will be updated showing the date that file was uploaded.



The screenshot shows a confirmation message: "Uploaded file waiting to be merged: infoscrn.dbf 6/27/2008". Below the message is a "Restart Upload" button.

3. Click the "Merge Data" button to begin the merge process after uploading.
4. After the merge completes you will be reminded to rename or move your equipment file.



If you use the Manual Upload feature it is important that you rename or move your equipment file after each merge. Failure to do so may cause your screening equipment software to append the next export to the information you previously merged. Re-merging the same baby records requires more staff effort and duplicates screening test outcomes unnecessarily.

## Automatic Uploading Equipment Files using the Merge Companion


The HiTrack Merge Companion is a small tool (100KB) that can be downloaded onto your computer to automate screening equipment file merges. The Merge Companion uses a technology called "Web Services" to upload and merge equipment files. To use this tool you must have Microsoft's .Net Framework 2.0 installed on the computer you use to access HiTrack. (Please contact the [HiTrack Help Desk](#) for more information about this requirement.)

Merge Companion must be downloaded and configured the first time you use it.

*To download and configure the Merge Companion (one time only):*

1. Click on the Merge icon from the toolbar. The Merge screen is shown with the Manual Upload tab selected by default.

Select the Automated Upload tab.




From this page you can download the Merge companion by clicking on the [Merge Companion](#) link with the Right mouse button and selecting "Save Target/Link as...".

Download the exe file to your Desktop (or any other easy to remember folder location).

2. After downloading the Merge Companion exe locate the file from the Desktop and run it by double clicking on the icon.

The Merge Companion screen will be shown and the following message will appear.



*(continued)*

3. Click on the "Configuration Options" link. The Configuration options will be shown.

The screenshot shows a window titled "HITRACK 4 Merge Companion" with a "Configuration Options" tab. The form contains the following fields and controls:

- User:** Text box containing "Dora".
- Facility:** Text box containing "Green River Regional".
- URL\*:** Text box containing "https://state.hitrack.gov/Hitrack\_4/hitrack\_4\_ws/WSMergeCompanion.asmx". Below this is a note: "\*Please provide the URL to the HITRACK 4 Merge Companion Web Service at your site." and a "Save" button.
- Merge New Data:** A section containing:
  - Select a Source:** A dropdown menu showing "Natus ALGO 2e/3 (ALGO.TXT) (algo.txt)".
  - Source Location:** Text box containing "algo.txt".
  - Password:** Text box with a "Start" button.
- Status:** A text area showing "Checking for file..." and "Not found." with a "Close" button at the bottom right.

Fill in the HiTrack User name, Facility and Web Service URL fields.


The URL is typically the same as your Web Companion site address followed by: "/HiTrack\_4/HiTrack\_4\_WS/WSMergeCompanion.asmx".

4. Click on the "Save" button. The Configuration Options will close and the Merge Companion will load your Facility's [Merge Sources](#).

The Merge Companion will save your configuration for next time and does not need to be downloaded or configured again.

To merge an equipment file using the Merge Companion:

1. Double click on the Merge Companion icon from the Desktop (or wherever you saved it after downloading). The Merge Companion screen will be shown.



2. Select your equipment [Merge Source](#) from the drop down in and then type in your HiTrack password.

The "Start" button is activated if a file from the selected equipment is found. If it is not found, the "Start" button's text will become "Search" indicating that a search is required to find your specified file. Clicking on "Search" will automatically search the root directories of your system for the file. If it is not found, it will allow the user to search manually for the file's location. This search feature does not allow users to search for XML, ETF, or ETF4 files.

3. Click on "Start". The equipment file will be uploaded and the merge process will begin. After merging, the equipment file on your computer will be automatically renamed to indicate that the file has been processed. For more information refer to the [File Location details in Appendix A](#).
4. Any pending merge batches should be resolved by logging into the Web Companion and clicking on the Merge icon from the toolbar.

# Appendix A: Linking HiTrack with Equipment Manufacturer Software

## ***Manufacturer Software Instructions***

This section discusses linking HiTrack with a few specific screening systems. Screening systems software changes frequently so please contact the screening manufacturer if there are any questions beyond what is covered here in this manual.

The instructions below refer to the HiTrack Windows user interface. For information about how to Merge into HiTrack's Web Companion see the [documentation earlier in this document](#).

Four commonly used screening software packages are detailed in this section, Bio-logic ABAer and AOAE, ALGO 2e Color, ALGO 3 and Otodynamics EZScreen2.

## **About File Locations**

For any screening software link you need to establish a file location that is entered into both HiTrack and the screening software. This location will be used by the screening software to write the file and used by HiTrack when merging the file.

Any file path on the computer may be used when HiTrack and the screening software are installed on the same computer as long as your Windows User account has read and write access to it. In most cases a new folder such as "C:\merge files\" can be created as used as an export location.

When HiTrack and the screening software are not on the same computer, removable media may be used. In the following directions, note that any removable media such as USB Flash Memory drives may be used in cases where terms "A:" and "diskette" are used. File paths on network drives can also be used for equipment data files when HiTrack is on a different computer. Please call the [HiTrack Help Desk](#) for additional file location support.

## **File Processing**

After HiTrack is able to successfully read in records from a merge, the file is renamed with a date and time "stamp". This prevents future screening activity from writing over the completed file and also prevents unnecessary re-merging of completed files.

The date and time stamp is added to the end of the file name. For example if you merge an "ALGO.TXT" file on March 21<sup>st</sup> 2006 at 11:02 it will be renamed to "ALGO\_TXT\_2006-03-21-11-02-52.old" after merging.

## **Program Information Sharing**

[HiTrack's Program Information](#) can be exported automatically when the HiTrack Windows user interface is installed on the same computer as the screening equipment software. To enable this setting turn on the "Auto Export PI" (Program Information) feature listed under [Other Settings](#) for your Facility. When this setting is enabled changes made to Medical People or Pick lists are automatically exported to HiTrack's "ProgInfo" folder (the default path is "C:\Program Files\NCHAM\HiTrack 4 for Windows\ProgInfo"). Some screening equipment software can then be configured to use this folder to import the lists.

When HiTrack is not installed on the same computer the [Program Information tab in Tools](#) can be used to export the lists to a removable storage device.

### **Free Text Program Information Support**

Instead of importing [HiTrack's Program Information](#) to use in drop-downs some equipment manufacturers allow "Free Text" entries collected from their software. Although this method of exporting data is not recommended by the HiTrack Help Desk it is a way to obtain patient related entries from rudimentary software.

"Placeholder Created" warnings are shown during the Merge process when HiTrack encounters a Free Text entry it does not recognize. The Fix Unknown Code Placeholders tool in HiTrack Advanced Tools should be used regularly to keep your lists of Program Information clear of duplicate entries.

### **Automating the HiTrack Merge**

With the help of IT or Data support staff the HiTrack Merge can be automated to run on a recurring schedule without user interaction. As with [other merges](#), any records that have questionable matches or need manual attention are left in the Pending Data merge screens.

A Windows "bat" file should be created with the following data fields in sequence (separated by spaces):

1. The location of the HiTrack 4 executable.
2. The action to perform (MERGE)
3. The Merge Source ID to merge as (found on the Facility Merge Sources screen)
4. The HiTrack username to merge with
5. The password for the HiTrack 4 username
6. The Facility ID (GUID) to Merge for (can be found in Settings → Select the hospital → Edit → Advanced)
7. The Facility or Facility Group ID to use to obtain Letter settings.

Sample "bat" file contents:

```
"C:\Program Files\NCHAM\HiTrack 4 For Windows\HiTrack_4_Win_PL.exe" MERGE 10 usern passw  
56045B2B-DB17-42DE-BBDC-735F8EA3B5C9 56045B2B-DB17-42DE-BBDC-735F8EA3B5C9
```

The IT or Data support staff should contact the [HiTrack Help Desk](#) to complete the setup for automated merges.

Please review the [Data Integration section](#) for more information about automating data merges.

## **Bio-logic System Corp ABaer and AOAE**

### **(Version 3.17.0.i.020)**

Demographic data for babies to be screened can be entered into the HATS™ (Hearing Assessment Tracking System) that comes with the Bio-logic equipment. To set the exporting file path click Setup on the main screen of the Bio-logic® Hearing Screening module (shown below), select Export Path and then click HiTrack. If the Bio-logic software and HiTrack are on the same computer, enter C:\Program Files\Bio-logic Systems Corp\ABaer. If they are on different computers, enter A:\.

Before starting verify that Biologic AudX/ABaer was added to your Facility [Merge Sources](#) and that the file location (i.e., the location from which the data will be coming) is set to c:\Program Files\Bio-logic Systems Corp\ABaer if Biologic is installed on the same computer as HiTrack or to A:\ if Biologic is installed on a different computer.

### **Link Setup**

*Use the following steps to set your [file location](#) and program information in the Biologic software:*

1. Click the Biologic Hearing Screening icon from the Desktop
2. From the menu bar select Setup > Export Path > HiTrack  
(If Export Path does not show in your menu use the Switch Profile to change to the Administrator user.)  
A pop-up with an “Export to” path is shown.
3. Use the “...” button or type in the file location path and then click the "Ok" button.
- 4a. To import Program Information (in later versions of HATS) when HiTrack is installed on the same computer: Use Setup > HiTrack List Import. Browse for "C:\Program Files\NCHAM\HiTrack 4 For Windows\ProgInfo". Also fill in the "Hitrack 4" checkbox.
- 4b. To import Program Information (in later versions of HATS) when HiTrack is installed on a different computer you will need to first export the Program Information from HiTrack onto a floppy disk or removable flash drive. Next, use Setup > HiTrack List Import, browse for "A:\" (or the drive letter of the flash drive) and fill in the "Hitrack 4" checkbox.

## Data Entry

Data can be entered into the HATS system in one of two places. The first option is when the screening tests are actually conducted on the front screen. A limited amount of data can be entered in the Patient Information form shown below.

The second option is to click the HATS Access button in the bottom left hand corner of the above screen and use the somewhat more comprehensive data collection form shown below.

Any data entered on the above form is automatically copied into the Hearing Screening form when screening is done. Regardless of which form is used, when a hearing screening test is completed, the demographics and test results for each ear are stored in a database from which it can be exported to a file for merging into HiTrack. Some additional information fields can be added to HATS using instructions found in the HATS manual.

The following table lists the HATS data that can be transferred from Bio-logic equipment and the HiTrack field to which it is imported:

Data enter into HATS as:	Is imported into HiTrack as
ID number*	Child's Medical ID
Last Name*	Child's Name-Last

First Name	Child's First
Birth Date*	Child's Birth-Date
Sex	Child's Gender
TOB	Time
Birth Hospital	Birth Facility
Doctor	Physician
Alternate ID#	Alt Rec #
Multi-Code	Order
Nursery Type	Nursery
Inpatient/Outpatient	Level
Weight (in grams)	Child's Birth-Weight
Notes	Notes
Risk Factors -AP	Risk Indicators Apgar scores of 0-4 at 1 minutes or 0-6 at 5 minutes
Risk Factors -BW	Risk Indicators Birth Weight
Risk Factors -CF	Risk Indicators Cranio Facial Anomalies
Risk Factors -CI	Risk Indicators Congenital (in utero) infection
Risk Factors -FX	Risk Indicators Family History
Risk Factors -HL	Risk Indicators Stigmata
Risk Factors -HY	Risk Indicators Hyperbilirubinemia
Risk Factors -ME	Risk Indicators Bacterial Meningitis
Risk Factors -MV	Risk Indicators Mechanical Ventilation
Risk Factors -OD	Risk Indicators Ototoxic Medication
Mother Last Name	Mother's Name-Last
Mother First Name	Mother's First
Mother DOB	Mother's BirthDate
Mother ID #	Mother's Contact ID
Mother Language	Mother's Language
Street Address	Mother's Address line 1
Apartment Number	Mother's Address line 2
City	City

State	State
Postal Code	ZIP
Telephone	Phone

\*Indicates fields which must be added to the basic HATS screen.

The following screening information is imported into HiTrack: Screening type (OAE or ABR), Left & Right Results, Test Date, Test ID, Technician-Screener, Screening Site, Level Inpatient/Outpatient

Currently, the Bio-logic system does not transfer some data entered into HATS (e.g., Patient Outcome at Discharge, Status, History, Testers' Notes) to HiTrack. Also, only one Risk Factor may be selected, if more than one is selected none transfer over to HiTrack. These may change in the future.

*To import data that has been collected in HATS into HiTrack:*

1. Enter the child's demographic data and collect screening data using the Bio-logic Hearing Screening module or the HATS module
2. From the Bio-logic Hearing Screening module, select Patients >Export>HI\*TRACK.
3. Minimize or Exit from the Bio-logic Hearing Screening module.
4. On the HiTrack main screen, click the Merge button.
5. From the Merge Type drop-down menu, select Biologic AuDX/ ABAer.
6. Click the Merge button.

All of the data collected using the ABAer and AOAe equipment for the date range specified will be imported to HiTrack. When the merge is completed, a message appears telling you how many records were merged successfully.

### **Export tips**

You probably want to export (download) infant records from HATS and merge them into HiTrack on at least a weekly basis in order to generate timely letters and up-to-date reports. During the export process, a copy of the data is sent to HiTrack. The original data is retained in the HATS database. Exporting data from HATS and merging it into HiTrack is always a two-step process:

1. Records from HATS are automatically stored at the specified location as data for each child are collected. These records remain at that location and new records are added to the database as infants are screened. When an Export is done, the information is copied to the file location specified. Caution: Any data exported but not merged into HiTrack will be **OVERWRITTEN**.
2. The data from the exported file is merged into the HiTrack database by clicking the Merge icon from the Toolbar and clicking the "Merge" button. Refer to the [Merging process documentation](#) in this manual for more information. After a completed merge the export file will be renamed with a date and time stamp.

## **Natus ALGO 2e Color or 3 with DataBook**

Because the ALGO equipment cannot be used for purposes other than screening, HiTrack must be installed on a different computer. DataBook is not designed to import your Program Information from HiTrack or to make that available to screeners during data entry. Therefore, screeners must take great care to follow the data entry conventions described below.

Before starting verify that Natus ALGO 2e Color/3 was added to your Facility [Merge Sources](#) and that the file location (i.e., the location from which the data will be coming) is set to a folder on the local hard drive if HiTrack is installed on the same computer as or to A:\ if HiTrack is installed on a different computer.

### **Link Setup**

*To set up your DataBook and HiTrack system:*

1. Open HiTrack.
2. Click the Settings button and choose the Medical Persons category on the left.
3. Fill in the "Only List People marked as Show" option and click the "Refresh" button.
4. Using the right side button on the mouse right-click and select Print Grid.

(Repeat these steps for other Physicians and Audiologists.)

This step print lists of the HiTrack screener and physician codes that must be used during data entry. See "Data entry and export tips" below for detailed information.

5. Open the ALGO 2e Color or 3 Databook
6. On the User Preferences screen, the following fields should be checked and the date format set as follows:

Medical Record No. required

Date of Birth required

Last Name required

Gender required

Date format: mm/dd/ccyy (Setting the date format may not be sufficient to ensure proper data entry—it is recommended that you post the data entry format for screeners directly on the screening equipment to avoid merge errors.)

## Data entry and export tips

Programs using DataBook software for infant hearing screening must follow the data entry guidelines provided below to ensure smooth transfer of information into HiTrack. Screening data that has not been entered in this specific format is not recognized by HiTrack.

*Use the following conventions during data entry into DataBook:*

1. Medical Record No: Enter up to 10 digits or characters.
2. Name, Last: Enter infant's last name.  
Name, First: Enter infant's first name or BB or BG (upper or lower case) to designate baby boy or baby girl.
3. Gender: Check Male or Female
4. Pediatrician: Enter Physician code as designated on the HiTrack code list.
5. Date of Birth: Enter as mm/dd/ccyy
6. Birth Location: Enter I, II, or III as Nursery level in the initial spaces. HiTrack ignores additional information.
7. Other: Enter O for Outpatient screening; enter I for Inpatient screening; HiTrack will ignore additional information.
8. Comment: Field not read by HiTrack.
9. User ID: Enter Screener code as designated on the HiTrack code list.
10. Risk Indicators: Leave as U if unknown; enter Y if present; enter N if certain it is not present. (Note: Severe asphyxia appears in HiTrack's Other Risk Indicator field.)

You probably want to export (download) infant records from DataBook and merge them into HiTrack on at least a weekly basis in order to generate timely letters and up-to-date reports. During the export process, a copy of the data is sent to HiTrack. The original data is retained in the DataBook system. Exporting data from DataBook and importing (merging) it into HiTrack is always a two step process:

1. From the Databook module, infant records are exported to a diskette by date range (refer to DataBook documentation for more explicit instructions). Always name the download file ALGO (DataBook automatically adds an extension so the file on the diskette is named ALGO.TXT). Do not download subsequent ALGO data onto the same diskette until you have completed the merge process, if you do, the first data file is OVERWRITTEN.
2. The data from the exported file is merged into the HiTrack database by clicking the Merge icon from the Toolbar and clicking the "Merge" button. Refer to the [Merging process documentation](#) in this manual for more information.

After a completed merge the export file will be renamed with a date and time stamp.

## **Natus ALGO 2 or 2e with DataBook**

Because the ALGO equipment cannot be used for purposes other than screening, HiTrack must be installed on a different computer. DataBook is not designed to import your Program Information from HiTrack or to make information available to screeners during data entry. Therefore, screeners must take great care to follow the data entry conventions described in this section.

Before starting verify that Natus ALGO 2/2e was added to your Facility [Merge Sources](#) and that the file location (i.e., the location from which the data will be coming) is set to a folder on the local hard drive if HiTrack is installed on the same computer as or to A:\ if HiTrack is installed on a different computer.

### **Link Setup**

*To set up your DataBook and HiTrack system:*

1. Open HiTrack.
2. Click the Settings button and choose the Medical Persons category on the left.
3. Fill in the "Only List People marked as Show" option and click the "Refresh" button.
4. Using the right side button on the mouse right-click and select Print Grid.

(Repeat these steps for other Physicians and Audiologists.)

HiTrack prints a lists of the screener and physician codes that must be used during data entry

### **Data entry and export tips**

Programs using DataBook software for infant hearing screening must follow the data entry guidelines provided below to ensure smooth transfer of information into HiTrack. Screening data that has not been entered in this specific format is not recognized by HiTrack. Because DataBook software does not provide separate data entry fields for vital information such as infant birth date, or inpatient screening vs. outpatient rescreening, screeners must modify the way that data is entered into available fields. An example of correct data entry on the DataBook screen is shown and explained below.

Also note that DataBook does not have a function to automatically recall infant information entered during a previous screening session so screeners must enter all information, following this specific format, every time a screening session is conducted even when successive screening attempts are made on the same infant.

Table 3: Example of Data Entry on the ALGO 2 DataBook Screen

Screen Prompt	What screener enters
Type first name, then press Enter:	07/01/97 bb1 (Birth date, sex, 1 = inpatient; 2 = outpatient) There must be 6 digits entered in the date field. If the month or day is a single digit, add a zero as shown above. There must be a space between the year and the sex of the infant.
Type last name, then press Enter:	Smith
Type patient ID number, then press Enter:	1254324
Type operator initials, then press Enter:	MDL (Designated screener codes/initials)

Screeners should use the code list created in the setup steps as they enter information into the Algo.

- 1 In the first name field, do not enter a first name; instead, enter the infant's birth date, as shown above, with the month/day/year separated by either slashes or hyphens, followed by a space. Then enter the letters bg (baby girl) or bb (baby boy), in either upper (BB) or lower case (bb) letters.
- 2 In the last name field, enter the infant's last name in either upper or lower case letters.
- 3 In the patient ID number field, enter the infant's medical record or patient ID number. (Most hospitals assign unique medical identification numbers to each twin, triplet, etc. If your hospital does not, you should add a multiple birth code to the infant's patient ID number, for example: 1254324A, 1254324B, etc.)
4. In the operator initials field, enter the three character or three digit code for the screener as shown on the HiTrack code list. These codes will be matched in HiTrack when the DataBook records are merged into HiTrack

You probably want to export (download) infant records from DataBook and merge them into HiTrack on at least a weekly basis in order to generate timely letters and up-to-date reports. During the export process, a copy of the data is sent to HiTrack. The original data is retained in the DataBook system. Exporting data from DataBook and importing (merging) it into HiTrack is always a two step process:

1. From the Databook software, infant records are exported to a diskette by date range. (Refer to DataBook documentation for more explicit instructions.) Always name the download file ALGO.DAT because HiTrack looks for that file name when attempting to merge new data into your tracking database. Do not download subsequent ALGO data onto the same diskette until you have completed the merge process, if you do, the first data file is OVERWRITTEN.
2. The data from the exported file is merged into the HiTrack database by clicking the Merge icon from the Toolbar and clicking the "Merge" button. Refer to the [Merging process](#)

[documentation](#) in this manual for more information.

After a completed merge the export file will be renamed with a date and time stamp.

## Otodynamics EZ-Screen2

### (Version 6.25.0.0)

The EZ-Screen2 software can be used with any of the Otodynamics equipment that has a keyboard for entering demographic data (e.g., ILO 88, Echoport plus, ILO 288, etc.).


Before starting verify that Otodynamics EZ-Screen was added to your Facility [Merge Sources](#) and that the file location (i.e., the location from which the data will be coming) is set to a folder on the local hard drive if is HiTrack installed on the same computer as or to A:\ if HiTrack is installed on a different computer.

### Link Setup

To link EZ-Screen2 with HiTrack, specify the following settings in EZ-Screen2:

1. On the EZ-Screen2 main screen, select the File menu, then Options. Select the Exports/Imports tab.

The following screen appears.



2. In the Screening data transfer area, field enter the merge file location that was entered in HiTrack for your EZScreen Merge Source. (The default is C:\Program Files\Otodynamics\ILO4win\Transfer.)
3. From the Export data to: drop-down menu, select NCHAM, HI\*TRACK.
4. Click the OK button to retain these changes.

If HiTrack and EZScreen2 are on the same computer, the user can configure the shared Program Information before merging into HiTrack.

*To configure Program information:*

- 1 Open EZScreen2 and make sure to be log in as an ILO Administrator.
- 2 Select Data Management > User Password > Configure.

Listed here will be all the users for EZScreen2 and screeners and hospitals for HiTrack.

- 3 Use the arrows to assign HiTrack codes to their EZScreen2 equivalents.


If the user changes the hospital affiliation, then all babies at this point and on will have the new hospital information. This does not change the hospital information for existing babies in HiTrack.

4. Click the OK button to retain these changes.

### Data Entry

Refer to the EZ-Screen documentation to select mandatory fields for data entry. Most, but not all, of the demographic data collected by EZ-Screen2 can be imported or merged with HiTrack.

The EZ-Screen2 data entry screen is shown below:



When demographic data is entered into EZ-Screen2 and a test is completed, a copy of the data and test results for each ear are stored in the location specified in the EZ-Screen2 configuration. The default location is C:\program files\Otodynamics\ILO4win\Transfer and the file name is EZScreen.txt.

The following table lists what EZ-Screen2 data can be transferred and the HiTrack field to which it is imported:

Data entered in EZ-Screen2 as:	Imported in HiTrack as:
Family Name	Child's Name-Last
First name	Child's First

ID number	Child's Medical ID
Mother's Name	Mother's Name-Last
Weight (kg)	Child's Birth-weight
Date of birth	Child's Birth-Date
Time of birth	Time
Male/Female	Gender
Status	Nursery
Multiple births	Order
Notes	Notes
Secondary ID	Alt. Rec. #
Mother's Forename	Mother's First
Mother's Maiden Name	Mother's Maiden
Mother's ID Number	Mother's Contact ID
Mother's Address (house/street)	Mother's Address line 1
Mother's Address (town)	Mother's Address line 2
Mother's Address (city)	Mother's City
Mother's Address (state/province)	Mother's State
Mother's Address (postcode/zip)	Mother's Zip
Mother's Phone number	Mother's Phone
Mother's Date of Birth	Mother's Birthdate
Guardian's Family Name	Mother's Alt Name-Last
Guardian's Forename	Mother's Alt First
Guardian's Address (house/street)	Mother's Alt Address line 1
Guardian's Address (town)	Mother's Alt Address line 2
Guardian's Address (city)	Mother's Alt City
Guardian's Address (state/province)	Mother's Alt State
Guardian's Address (postcode/zip)	Mother's Alt Zip
Guardian's Phone number	Mother's Alt Phone
Risk Factors Family History	Risk Indicators Family History
Risk Factors Low Birth Weight	Risk Indicators Birth Weight
Risk Factors Bacterial Meningitis	Risk Indicators Bacterial Meningitis

Risk Factors Ototoxic Medication	Risk Indicators Ototoxic Medication
Risk Factors Syndrome	Risk Indicators Stigmata
Risk Factors Hyperbilirubinemia	Risk Indicators Hyperbilirubinemia
Risk Factors Congenital Infection	Risk Indicators Congenital (in utero) infection
Risk Factors Cranio Facial Anomalies	Risk Indicators Cranio Facial Anomalies
Risk Factors Mechanical Ventilation	Risk Indicators Mechanical Ventilation
Risk Factors Abnormal Apgar Scores	Risk Indicators Apgar scores of 0-4 at 1 minute or 0-6 at 5 minutes
Risk Factors Other Risk Factors	Risk Indicators -Other tab (brings over a 'Y' if selected in EZScreen)

The following screening information is imported into HiTrack: Screening type (OAE or ABR) , Left & Right Results, Test date, Test ID ,Screener ID (if pick lists are set up).

Even though data for other fields can be entered into EZ-Screen2 (e.g., Middlename, Alternate Contact ID # etc.), they are not presently transferred to HiTrack.

All of the data collected during screening since the last time you merged will be imported to HiTrack. When the merge is completed, a message appears telling you how many records were merged successfully.

### Export tips

You probably want to export (download) infant records from EZScreen2 and merge them into HiTrack on at least a weekly basis in order to generate timely letters and up-to-date reports. During the export process, a copy of the data is sent to HiTrack -the original data is retained in EZ-Screen2. Exporting data from EZ-Screen2 and importing (merging) it into HiTrack is always a two-step process:

1. Records from EZ-Screen2 are automatically stored at the specified file location as data for each child are collected. This data remains at that location and new data is appended to it until a merge is done.
2. The data from the data file is merged into the HiTrack database by clicking the Merge icon from the Toolbar and clicking the "Merge" button. Refer to the [Merging process documentation](#) in this manual for more information.

After a completed merge the export file will be renamed with a date and time stamp.

## Appendix B: Report Criteria Documentation

### Classifications

Classifications are used to categorize babies on reports including the Flow Chart and other summary reports. Each baby may be classified for each [EHDI Stage](#): Inpatient, Outpatient, Evaluation and Intervention. A baby may have at most four classifications, one for each stage.

Classifications are based on Most Conclusive Results and other parts of the baby record. Most Conclusive Results are computed for each ear at Inpatient and Outpatient stages. Most Conclusive Outcomes are also defined for Hearing Disposition but must be manually selected for each infant after any type and degree of hearing loss is entered or changed.

The criteria lists below are presented in the order of consideration. Once a qualifying classification is found for each stage no other criteria is considered.

Updated 5/19/2008.

### Inpatient Classifications

Last updated: 10/9/2007

	<i>Classification</i>	<i>Criteria</i>
1	<i>Deceased before/during Inpatient Screening</i>	Infant is marked as deceased and no screening results exist
2	<i>Inpatient Pass</i>	Both Inpatient Most Conclusive Results are Pass
3	<i>Inpatient Referred</i>	Both or one Inpatient Most Conclusive Result is Refer
4	<i>Inconclusive Inpatient Results</i>	Other Inpatient screening results exist that are not Refer
5	<i>Missed Inpatient Screening</i>	Infant's most recent Tracking Action for Inpatient Screening is "Missed" and the baby has no Transfer entries.
6	<i>Refused Inpatient Screening</i>	Infant's most recent Tracking Status for Inpatient Screening is "Refused".
7	<i>Transferred Out before/during Inpatient Screening</i>	Infant's most recent Tracking Status for Inpatient Screening is "Transferred Out" or a Transfer entry exists.
8	<i>No Inpatient Information</i>	(default Classification)

## Outpatient Classifications

Last updated: 10/9/2007

	<i>Classification</i>	<i>Criteria</i>
1	<i>Deceased before/during Outpatient Screening</i>	Infant is marked as deceased, Inpatient screening results exist and no Outpatient screening results or Hearing Disposition is entered.
2	<i>Outpatient Pass</i>	Both Outpatient Most Conclusive Results are Pass.
3	<i>Outpatient Refer</i>	Both or one Outpatient Most Conclusive Result is Refer.
4	<i>Outpatient Inconclusive Results</i>	Other Outpatient screening results exist that are not Refer.
5	<i>Refused Outpatient Screening</i>	Infant's most recent Tracking Status for Outpatient Screening is "Refused".
6	<i>Scheduled for Outpatient Screening</i>	Infant's folder is "Need Outpatient Screening" and no Outpatient results exist.
7	<i>Follow-Up Discontinued before/during Outpatient Screening</i>	Infant marked Discontinued while in "Need Outpatient Appt.", "Need Outpatient Screening" or "Need to Locate" folders. (Excludes the infant if they need to be located for Evaluation.)
8	<i>Need to Locate for Outpatient Screening</i>	Infant does not have Outpatient screening results and their folder is "Need Outpatient Appt." or "Need to Locate".
9	<i>Go directly to Evaluation without Outpatient Screening</i>	Infant does not have Outpatient screening results and any of the following: a. folder is "Need Evaluation Appt.", "Need Evaluation", or "Completed Evaluation". b. "Evaluate w/o Outpatient" Action exists in Tracking History. c. Hearing Disposition has been entered.
10	<i>No Outpatient Information</i>	Infant does not have Outpatient screening results

## Evaluation Classifications

Last updated: 4/28/2008

	<i>Classification</i>	<i>Criteria</i>
1	<i>Deceased before/during Eval</i>	Infant is marked as deceased, Outpatient screening results or Hearing Disposition is entered and no EI Services are entered. (Single Stage Protocols check for Inpatient Screening.)
2	<i>Confirmed Loss after Evaluation</i>	Infant has Most Conclusive Hearing Disposition where any ear is entered with a confirmed type of Loss or Sound Field is entered and confirmed.
3	<i>Normal Hearing after Evaluation</i>	Infant has Most Conclusive Hearing Disposition where a. both ears are confirmed Normal ("N/A (No Loss)" or "None" (for HT3.5 converted data)) d. one ear is entered and confirmed normal c. sound field is entered and confirmed normal
4	<i>Refused Evaluation</i>	Infant's most recent Tracking Status for Evaluation is "Refused".
5	<i>Follow-Up Discontinued before/during Evaluation</i>	Infant's most recent Tracking Status for Evaluation is "Discontinued".
6	<i>Out of Jurisdiction – No Evaluation</i>	Infant with a "Transfer Out of System" History record entered while the record is in "Need Evaluation Appt" or "Need Evaluation" and no Hearing Disposition entries exist. (Transfer date must be after DOB)
7	<i>Evaluation in Process</i>	Infant has any of the following a. Most Conclusive Hearing Disposition with unconfirmed type and degree of Loss entered b. An "Evaluate without Screening" Tracking Option has been used.. c. Outpatient Most Conclusive Result of Refer, Invalid or Could not Test.
8	<i>Need to Locate for Evaluation</i>	Infant's folder is "Need Evaluation Appt." or "Need to Locate" and a Tracking Action of "Need to Locate" during Evaluation exists.

## Intervention Classifications

Last updated: 10/26/2012

	<i>Classification</i>	<i>Criteria</i>
1	<i>Deceased before/during Intervention</i>	Infant is marked as deceased and an EI Service enrollment has been entered.
2	<i>Refused Intervention</i>	All EI Service entries are marked as "Parent Opt Out" or the infants most recent Tracking Status for "Referred for Intervention" is "Refused".
3	<i>Follow-Up Discontinued before/during Intervention</i>	Infant's most recent Tracking Status for "Referred for Intervention" is "Discontinued".
4	<i>No Intervention Info</i>	Infant has no EI Service enrollments entered that are "Pending", "Enrolled", or "Monitoring"
5	<i>Enrolled in Intervention</i>	Infant has at least one EI Service enrollment entered that is "Pending", "Enrolled", or "Monitoring"

## ***Amplification Report***

This report lists babies who have Amplification entries and provides calculations for the JCIH 2007 Benchmark for Amplification. Deceased babies are not included in this report.

### **Report Specific Filtering Options**

This report uses the general filtering options.

### **JCIH 2007 Benchmark**

This report includes the JCIH Benchmark for amplification:

“For families who elect amplification, the percentage of infants with confirmed bilateral hearing loss who receive amplification devices within 1 month of confirmation of hearing loss; the recommended benchmark is 95%.”

Although the benchmark applies only to infants with bilateral hearing loss, information has also been provided on infants who have unilateral hearing losses or where laterality is unknown.

In order to determine whether a baby meets the benchmark or not, two pieces of information are required: 1) the fitted date for the amplification device, and 2) a hearing disposition entry. Babies who do not have these two pieces of information will not be counted on the JCIH 2007 Benchmark report, but they will be counted on the Amplification report as long as they have an amplification entry.

A count of babies without a fitted date is included on the JCIH 2007 Benchmark report to show how many babies are not included on the report because the fitted date is missing. The Total Babies count includes babies without fitted dates, babies without hearing disposition entries, and babies with bilateral, unilateral, and laterality unknown hearing losses.

## CDC Survey Report

The HiTrack 4.6.3 CDC Summary Report follows the guidelines given in the "CDC-EHDI Hearing Screening and Follow-up Survey (HSFS): 2016 Explanations." This report also contains quality assurance feedback statistics not included by the CDC explanations. These additional numbers help program coordinators check the quality of their data.

### Report Specific Filtering Options

Birth date range	Specifies the Birth Date Range for the report.
Only include Baby Basket babies	Allows the CDC Survey report to be run for a small subset of records.
Home Birth Facility	Used to indicate which HiTrack Facility was used in records for babies born at home. Create a Facility Group of Home Births if multiple HiTrack Facilities have been used for births at home.
State Abbreviation for Non-resident calculation	The entry in this option is used to calculate non-resident totals.

### Record Browse List Feature

When generated in Browse List format HiTrack's CDC Summary Report allows users to see lists of babies that make up the counts for most of the survey categories. The records from these filtered lists can then be added to the Baby Basket for later review.

### Lost to Follow-Up % and Number (Diagnosis)

This statistic is now included at the end of section 1 of the CDC Report.

### Age Calculations and Types of Loss Considered



The CDC defines the date ranges as follows:

0 to 1 month = 0 to 30 days
0 to 6 months = 0 to 180 days
6 to 12 months = 181 to 360 days
12+ months = 361+ days

As prescribed by the CDC, ears Fluctuating Conductive losses are counted as Normal hearing ears.

### Non-Resident CDC Counts

Starting with the 2012 instructions the CDC has added a new count to each section of Part 1 that allows users to account for No Info/Unknown babies by categorizing them as Non-Resident. To determine if a baby is Non-Resident HT uses the following rules: If a state is specified it must

be a valid abbreviation or fully spelled out state name. If it is not a valid spelled out name or two letter abbreviation it is considered Non-Resident. Blank entries are counted as Residents. Only records without screening or diagnosis (depending on which count) qualify for the Non-Resident count.

A count of invalid Address entries is included in the Warnings List (see below).

### **Adoption Counts**

The CDC now allows “Adopted” to be used as an explanation for why Screening, Diagnostic or Intervention was not completed. HiTrack reports this number by checking the “Adoptive Parent” contact role specified in Demographics – Contacts. The Primary Contact must have this checked to qualify for these counts.

### **Hearing Loss Type and Severity Summary**

This section of the survey report summarizes types and degrees of loss for the report selection criteria. As required by the CDC, the hearing loss summary grid reports statistics based on the most complete diagnosis. To meet this requirement HiTrack uses selections based on the Most Conclusive Outcomes selected on the Hearing Disposition screen.

Each Type and Degree of loss category reports the number of ears recorded in HiTrack with that disposition. When ranged Degrees of Loss are used in records they are included in the counts with the most severe end of that range. (For example, Mild-Moderate is counted as Moderate and Severe-Profound is counted as Profound.) Babies with a type of loss that did not follow a normal screening protocol are not counted in this section, including babies for whom tracking was discontinued and babies that passed their screening results.

This summary also includes other statistics to help program coordinators understand their data, especially numbers of babies being excluded by the CDC survey reporting instructions.

There may be a difference in the number of babies with hearing loss on the summary grid when compared to counts on other parts of the survey. Parts 1 and 3 report statistics based on the initial diagnosis. The report includes two numbers that help account for this difference. "Total Cases Resolved" counts babies who initially had a type of loss that resolved to Normal hearing. These babies do not show on the type of loss grids. The second number, "Initial Normal Cases to Identified later as a Loss" counts babies who initially had confirmed Normal hearing, but developed a confirmed loss. This second number is provided for the Comments section of the survey.

### **Not Pass Babies with Type of Hearing Loss**

This number provides the number of babies represented by the ear counts in the Type and Degree grid. Types of loss for Late Onset/Atypical babies are excluded from this grid.

### **Total Number of Babies with Hearing Disposition Data**

This number provides the number of babies with any entries in Hearing Disposition. This can be used as a quality assurance feedback especially if it differs greatly from the number of babies reported on the grid. Differences in the number could be due to Sound Field entries or losses entered as "Undetermined" type of loss.

## Total Ears with Undetermined Type of Hearing Loss

This number provides the number of ears entered with Undetermined as a quality assurance feedback.

## Quality Assurance Tools

HiTrack provides several quality assurance tools to assist you in preparing your CDC Survey.

User Defined Report Columns	CDC Screening Outcome, Type of Loss (Initial), Degree of Loss (Initial)
User Defined Report Auto Filters	The output grid of all User Defined Reports can be dynamically filtered further depending on the values obtained from the report.
Report Advanced Filter Criteria	"CDC Initial Dx is Different" Yes/No "In Baby Basket" Yes/No
Baby Advanced Details	CDC Screening Outcome
Baby Basket	The baby basket can be used to keep an intermediate list of CDC qualifying babies to use for other Reports.

Please contact Help Desk to brainstorm additional quality assurance checks.

## Other Notes

Note that the table arrangement on the HiTrack report provides types of loss as columns although the 2011 CDC report shows types of loss as rows. This alternate arrangement does not affect the data provided in the table.

The "Profound" degree under Conductive loss is shaded grey to indicate that this is not a valid category. The number of ear losses for this category is provided as a quality assurance feedback.

When laterality is unknown the ear with a specified type of loss will be represented in the Laterality Unknown count.

Unknown Ear counts are not provided on this section because HiTrack requires left/right ear identification on audiometric data entry.

## Hearing Status for CDC Reporting

The CDC requests Initial Hearing Status for Parts 1 and 3 of their survey and the most current status for Part 2. The current status uses the Most Conclusive Outcomes selected by the user. The initial hearing status defaults to the first chronological entries. Beginning with HiTrack 4.6.2 users in the System Administrator group may override the default initial disposition entries for each baby. This must be done on a case by case basis.

In bilateral cases where the type of loss is different in each ear the following rules are used:

Ear 1 Type	Ear 2 Type	"Head" Type of Loss
Sensorineural	Conductive, Mixed, Neural, Fluctuating*	Sensorineural
Neural	Anything except Sensorineural.	Neural
Conductive	Mixed	Conductive
Conductive	Fluctuating	Conductive
Fluctuating	Mixed	Mixed

### Missed ears

CDC and JCIH have recently recommended that both ears should always be rescreened at outpatient screening appointments. From the JCIH "Year 2007 Position Statement" Statement:

*"For rescreening, a complete screening on both ears is recommended, even if only 1 ear failed the initial screening."*

The following table summarizes how the HiTrack 4 CDC Survey reports babies with only one ear screened.

Inpatient	Outpatient	CDC Report Outcome
Single ear Pass, Refer, Invalid or Could Not Test.	[blank]	Missed
Inconclusive (both ears screened or single Pass ear)	Single ear Pass, Refer, Invalid or Could Not Test	Not Passed (reflects Inpatient result)
Refer (both ears screened)	Single ear Pass, Refer, Invalid or Could Not Test	Refer (reflects Inpatient result)

### Warning List Feature

After generating the CDC Report, HiTrack provides a list of babies and issues that may be of interest for data entry completion purposes. The following categories are reported:

Names of babies who have an Initial Loss reported on Part 1, but not on Part 2
Names of babies who have Unknown Diagnosis
Names of babies who were Not Referred to EI
Names of babies with have Unknown EI Services
# invalid address State Fields
# of NICU Days entries

Babies from the Warnings List can be added to the baby basket with using the baby basket button on the warning list screen. Once in the baby basket, this list can be used as a to-do list for data entry completion or deeper reporting purposes.

### **User Defined Report Columns**

Several are available in User Defined Reports to help locate and verify data to be submitted:

Evaluation Head Result – Hearing Disposition
EI Services Outcome – EI Services
Primary Contact State – Contacts
Primary Contact Zip – Contacts

## **CDC Survey Report Counts**

The following sections of this document explain how babies are counted in each survey question.

Based on CDC reporting instructions, Inpatient and Outpatient hearing screening result outcomes have been combined into overall outcomes of Pass, Not Pass/Inconclusive, Refer, Refused, Missed, and Deceased. For more information about these combinations see the [Appendix](#).

### **Survey Counts**

#### **Total Occurrent Births**

This reports all baby records that are known to and born at the current Facility or Facility Group including Home Births.

#### **Total Screened**

*Calculated: Total Passed + Total Not Pass*

#### **Total Passed**

This reports babies with Pass as their [combined Screening outcome](#).

#### **Total Not Pass**

This reports the number of babies with Refer or Not Pass/Inconclusive as their [combined Screening outcome](#). Babies who did not Pass Inpatient and were screened in only one ear and passed are also included in this count. This count is used later as a basis for babies that should be receiving diagnosis and is typically used as a denominator for the Lost to Diagnosis percent.

#### **Detailed Screening Results**

Counts are for the requested Inpatient and Outpatient screening result combinations.

#### **Total Not Screened**

*Calculated: Died + Non-Resident + Medical Reasons + Transferred + Adopted + Homebirth + Unresponsive + Unable to Contact + Unknown Screening*

#### **Not Screened – Infant Died/Parents Declined Services**

This reports babies with Deceased or Refused as their [combined Screening outcome](#).

#### **Not Screened Non-Resident**

This reports babies with No Information or Missed as their [combined Screening outcome](#) that also have a Primary Contact address entered with a different State Abbreviation than what was requested in the criteria. (Record with blank addresses are not considered Non-Resident.)

#### **Unable to be Screened Due to Medical Reasons**

This reports babies with No Information or Missed as their [combined Screening outcome](#) that have been moved to the “No Scrn. Medical Reasons” tracking folder.

### **Infant Transferred**

This reports babies with Transferred Out as their [combined Screening outcome](#).

### **Infant Adopted**

This reports babies with No Information or Missed as their [combined Screening outcome](#) that also have a contact entered as “Adoptive Parent”.

### **Homebirth**

This reports babies with No Information or Missed as their [combined Screening outcome](#) that also have Home Birth (as indicated by the CDC Report Criteria) as their Birthing Facility.

### **Parents / Family Contacted (for Screening) but Unresponsive**

This reports babies with No Information or Missed as their [combined Screening outcome](#) that have been moved to the “Contacted but Unresponsive for Screening” tracking folder.

### **Unable to Contact (for Screening)**

This reports babies with No Information or Missed as their [combined Screening outcome](#) that have been moved to the “Unable to Contact for Screening” tracking folder.

### **Unknown Hearing Screening**

This reports babies with No Information or Missed as their [combined Screening outcome](#).

### **Not Screened - Other**

This number is not reported from HiTrack.

### **Total with Documented Diagnosis**

*Calculated: Total with No Hearing Loss + Total with Permanent Hearing Loss*

### **Total with No Hearing Loss (Initial diagnosis)**

This reports the number of Not Pass babies with confirmed normal hearing. This includes babies with fluctuating loss in on ear and confirmed normal hearing in the other ear.

The computed age subcategories are based on the date of the Initial Outcomes selected on the Hearing Disposition screen.

### **Total with Permanent Hearing Loss**

This reports the number of Not Pass babies with confirmed hearing loss as a subset of Total Not Pass. The computed age subcategories are based on the date of the Initial Outcomes selected on the Hearing Disposition screen. Non-permanent Fluctuating losses (both bilateral and unilateral) are excluded from this count but included in the No Hearing Loss counts on the CDC Survey Report.

### **Total with No Documented Diagnosis**

*Calculated: Audiological Diagnosis in Process (Awaiting Dx) + All No Diagnosis counts*

### **Audiological Diagnosis in Process (Awaiting Diagnosis)**

This reports the number of babies with all of the following criteria:

Classified as Evaluation In Process (Outpatient Refer or Evaluate without Screening)
One ear-specific Hearing Disposition entry
One diagnostic Recommended Action six months since the last Hearing Disposition entry
One currently scheduled diagnostic Recommended Action at the time the CDC Survey Report was generated.

### **No Diagnosis – PCP/ENT Did Not Refer Infant for Diagnostic Testing**

This reports the number of Not Pass babies without Hearing Disposition entered that have been moved to the “PCP/ENT did not Refer for Evaluation” tracking folder.

### **No Diagnosis – Infant Died/Parent Declined Services**

This reports the number of babies who are classified as "Deceased during/before Evaluation" or "Refused evaluation"

### **No Diagnosis – Unable to Receive Diagnostic Testing Due to Medical Reasons**

This reports the number of Not Pass babies without Hearing Disposition entered that have been moved to the “No Eval. – Medical Reasons” tracking folder.

### **No Diagnosis – Non-resident or Moved Out of Jurisdiction**

This reports the number of babies who are entered as Transferred Out of System while in an Needing Evaluation folder.

### **No Diagnosis – Infant Adopted**

This reports the number of Not Pass babies without Hearing Disposition entered that also have a contact entered as “Adoptive Parent”.

### **No Diagnosis – Unresponsive/Unable to Contact**

This reports the number of Not Pass babies without Hearing Disposition entered that have been moved to the “Contacted But Unresponsive” or “Unable to Contact” Evaluation tracking folder.

### **No Diagnosis – Unknown**

This reports the number of Not Pass babies without Hearing Disposition entered that have no other documentation.

### **Early Intervention Counts**

The Early Intervention counts are based on several data entry fields. A baby must qualify for each of the items listed as criteria for the EI counts.

### **Total Referrals to Part C EI**

Refer or not pass the hearing screening results
---

Must have a confirmed hearing loss.
Must have a Recommended Action entry for Early Intervention Referral.

**Total Not Referred to Part C EI and Unknown**

Refer or not pass the hearing screening results
Must have a confirmed hearing loss.
Must not have a Recommended Action entry for Early Intervention Referral.

**Total Enrolled in Part C EI**

Refer or not pass the hearing screening results
Must have a confirmed hearing loss.
On the EI Services screen, there must be a record entered where the status is not "Parent Opt-out", "Not eligible/closed" or "Out of Jurisdiction." The checkbox must not be checked for "Private Enrollment (Non-Part C)." If there are multiple enrollments for a baby and some of them are Non-Part C enrollments while others are Part C, then the baby will still be counted in this section.
Each baby enrolled in EI Services must also have a separate IFSP EI Service recorded to meet the CDC's "signed IFPS" requirement

**Enrolled in Part C Before 6 Months of Age, Enrolled in Part C After 6 Months but Before 12 Months of Age, Enrolled in Part C After 12 Months of Age, Enrolled in Part C Age Unknown**

Refer or not pass the hearing screening results
Must have a confirmed hearing loss.
On the EI Services screen, there must be a record entered where the status is not "Parent Opt-out", "Not eligible/closed" or "Out of Jurisdiction." The checkbox must not be checked for "Private Enrollment (Non-Part C)." If there are multiple enrollments for a baby and some of them are Non-Part C enrollments while others are Part C, then the baby will still be counted in these sections.
A baby's age is determined by the difference between the birth date and the enrollment date. If there are multiple enrollments for a single baby, the earliest enrollment date is used in determining the age. If the birth date or the enrollment date is blank, then the baby will be counted in the "Enrolled in Part C Age Unknown" category. The date ranges for the CDC report categories are calculated by the number of days, but displayed on the report as months.

**Total Receiving Only Monitoring**

This reports the number of enrolled babies that only have Monitoring EI Service entries.

### **Total Services from Non-Part C EI Services Only**

Refer or not pass the hearing screening results
Must have a confirmed hearing loss.
On the EI Services screen, there must be a record entered where the status is not "Parent Opt-out", "Not eligible/closed" or "Out of Jurisdiction." The checkbox must be checked for "Private Enrollment (Non-Part C)." If a baby has multiple enrollments, they must all be non-part C enrollments in order to be counted in this section.

### **Non-Part C EI Services Before 6 Months of Age, Non-Part C EI Services After 6 Months but Before 12 Months of Age, Non-Part C EI Services After 12 Months of Age, Non-Part C EI Services Age Unknown**

Refer or not pass the hearing screening results
Must have a confirmed hearing loss.
On the EI Services screen, there must be a record entered where the status is not "Parent Opt-out", "Not eligible/closed" or "Out of Jurisdiction." The checkbox must be checked for "Private Enrollment (Non-Part C)." If a baby has multiple enrollments, they must all be non-part C enrollments in order to be counted in this section.
A baby's age is determined by the difference between the birth date and the enrollment date. If there are multiple enrollments for a single baby, the earliest enrollment date is used in determining the age. If the birth date or the enrollment date is blank, then the baby will be counted in the "Enrolled in Part C Age Unknown" category. The date ranges for the CDC report categories are calculated by the number of days, but displayed on the report as months.

### **No EI Services Infant Died/No EI Services Parents Declined Services**

Refer or not pass the hearing screening results
Must have a confirmed hearing loss.
If a baby has a tracking folder of Receiving Intervention prior to stopping tracking (from the Tracking Status screen) because of parent refusal, the baby will be counted in this category.
If a baby has one or more enrollment entries from the EI Services screen where the Status is "Parent Opt-out," then that baby will be counted in this category. However if that baby also has enrollment entries where the status is not "Parent Opt-out," then that baby will not be counted in this section of the report but will be counted in one of the "Enrolled" sections.
If a baby has a tracking folder of Receiving Intervention and an enrollment entry from the EI Services screen prior to marking the baby as Deceased, that baby will be counted as Deceased in this category.

### No EI Services Non-resident or Moved Out of Jurisdiction

This reports the number of hearing loss infants whose latest EI entry is marked as Out of Jurisdiction or whose contact address is in a state other than the state used for reporting.

### No EI Services Unresponsive/Unable to Contact/Medical/Reasons/Unknown

Refer or not pass the hearing screening results
Must have a confirmed hearing loss
If a baby has a tracking status of Receiving Intervention prior to stopping tracking by changing the status to "Follow-Up Discontinued," the baby will be counted in this category.
If a baby has a tracking status of Receiving Intervention and no enrollment entries on the EI Services screen, the baby will be counted in this category.
If a baby has one or more enrollment entries from the EI Services screen where status is "Not eligible/Closed," then that baby will be counted in this category. If that baby also has enrollment entries from the EI Services screen where the status is "Parent Opt-out," then that baby will still be counted in this category.
Records who have had their Tracking folder changed to No Services - Parent/Family Unresponsive or No Services - Unable to Contact or No EI Medical Reasons will count in these respective folders. Otherwise they will be counted as Unknown.

### Additional Cases Not Reported

#### Additional Cases of Hearing Loss not included in Diagnostic Data

This reports the number of babies with a confirmed hearing loss that also have Pass, Missed, Refused or No Info as their [combined Screening outcome](#). The computed age subcategories are based on the date of the Most Conclusive Outcomes listed on the Hearing Disposition screen.

#### Additional Cases of non-permanent, transient hearing loss

This reports the number of babies with Fluctuating Conductive hearing losses.

#### EI Services for Late Onset Hearing Losses

These counts match the logic used for Not Pass hearing losses, except they account for babies that passed or missed screening.

### Demographics Report

This final report section uses the same criteria as above categorized by CDC defined demographic sub categories.

Note that for each category, missing entries in HiTrack will be included in the Unknown count, except for Mother's Education which filters by the mother's age.

Persons with no race entries or multiple race entries are counted as "Other"

## **Data Cleaning Reports**

These three reports help program coordinators to find duplicate records for the purpose of combining records in HiTrack.



These reports help answer the following questions: "How can I identify duplicate records?" and "How can I locate duplicate records sent from different hospitals when babies are transferred?"

### **Duplicate Medical ID**

This report includes babies with duplicate Medical ID numbers in the system. Duplicate Medical IDs are allowed across birthing facilities but not within a facility. This report should be used to see if a baby has been entered under more than one facility with the same Medical ID.

### **Probable Duplicates**

This report includes babies who share the same last name and date of birth and includes birth order to help distinguish records for twins and multiple births. This report should be used to look for baby records that have been entered by more than one facility, especially transferred babies.

### **Custom Duplicates**

This report allows users to build a duplicate list based on a variety of field combinations.

### **Report Specific Filtering Options**

This report uses the general filtering options.

### **Combining Duplicates**

After duplicate records have been identified with these reports they can be combined using the [Record Combine](#) feature on the Baby List.

### **Early Intervention Report**

This report lists babies who have Early Intervention entries filtered by either status or service type. Deceased babies are not included in this report

### **Report Specific Filtering Options**

In addition to the general filtering options users may filter by Enrollment Status and Type of Service.

## **Hearing Status Report**

This report may include babies with a Confirmed Loss, Confirmed Normal Hearing or Suspected Loss, depending on the selection criteria options.

This report also includes the Milestone ages and Early Intervention enrollment if it is entered. Milestone ages are the baby's age in months for Screening, Identification of Loss, and Enrollment. These ages are delimited by slashes (/) with exclamation marks (!) noting for ages that have exceeded the target milestone.

### **Report Specific Filtering Options**

In addition to the general filter options this report allows to filter by hearing disposition confirmation:

Confirmed Losses, Suspected Losses, Confirmed Normal Hearing

Any combination of these loss filters may be used to help answer the following questions: "Which infants have not had their suspected loss confirmed yet?" and "What services are being given to infants with hearing loss?"



## **Needing Evaluation Report**

This report lists babies needing evaluation grouped by their Recommended Actions.

By default the baby must meet the following criteria to be included:

1	Not Deceased
2	"In Progress" Status
3	Classified as "Evaluation in Process", "Confirmed Loss after Evaluation" or "Risk Monitoring"

Completed Actions are not included for babies on this report. If the baby has Outpatient Refer or No OP→Dx then the baby will show up under No Recommendation/None.



This report helps program coordinators answer the following questions: "Which infants are scheduled for follow-up within the next (x) days?" and "Which infants need to be scheduled for follow-up?" (The latter question can be answered by referring the "No Recommendation/None" category.)

### **Report Specific Filtering Options**

In addition to the general filter options this report can be filtered by a specific type of recommendation or the date when the recommendation is due. The date boxes at the top of the filter page allow for appointment date filtering.

Disabling the Incomplete Records filter allows HiTrack to include actions for any babies, not only ones meeting the three criteria above. This may be useful for reporting actions recorded for babies that passed their hearing screening.

## ***Needing Outpatient Report***

This report includes babies with any of the following classifications:

Inpatient	Inpatient Referred Missed Inpatient Screening Inconclusive Inpatient Results Refused Inpatient Screening Transferred Out before/during Inpatient Screening
Outpatient	Outpatient Inconclusive Results, Need to Locate for Outpatient Screening' Scheduled for Outpatient Screening' No Outpatient Information'.

### **Report Specific Filtering Options**

This report uses the general filtering options.

## ***Outpatient Screening Report***

This report lists babies known to the current Facility or Facility Group with Outpatient screening results. The final pages of the report provide summary outcome counts and nursery level counts for babies listed on the report.



This report may help answer the question: "I need to report Outpatient screening differently for billing purposes; how do I get a list of patients who came back for Outpatient screening?"

### **Report Specific Filtering Options**

This report uses the general filtering options and also allows for filtering by screening outcome classifications.

## ***Physicians Report***

This report lists babies by their Primary Care Provider and provides an alternative to printing individual letters to inform physicians of screening results. The report can be generated for all physicians or for only those who wish to receive information in this format. If no patients of a particular physician were screened during the selected date range, no report is generated for that physician.

The list of babies for each Physician starts on a new page.



This report helps answer the questions: "How can I communicate patient hearing screening outcomes to doctors?" and "Which baby records do not have their doctor specified?"

## **Report Specific Filtering Options**

In addition to the general filter options this allows for filtering by screening classification and by specific doctor.

## ***Risk Indicator Report***

On this report the JCIH risk indicators are listed as columns across the top of the report with baby records showing as rows. The risk indicators abbreviation explanations and indicator statistics are included at the end of the report. A diagnostic and Hearing status evaluation summary is also provided at the end of the report.

This report includes babies that are marked as 'Yes' for any Risk Indicator.

Starting in HiTrack version 4.4.7 several indicators have been retired to conform to the JCIH 2007 Position Statement. These indicators are provided on this report for historical reporting purposes.



The new JCIH 2007 indicators are:

"Chemotherapy" [CHEM].
------------------------

"Neonatal intensive care of more than 5 days or any of the following regardless of length of stay: ECMO, assisted ventilation, exposure to ototoxic medications (gentimycin and tobramycin) or loop diuretics (furo-semide/Lasix), and hyperbilirubinemia that requires exchange transfusion" [NN5L]
--

These indicators were added because there was not enough justifiable overlap with any existing JCIH 1994 or JCIH 2000 indicator.

## **Report Specific Filtering Options**

This report uses the general filtering options and also allows for filtering by each indicator.

## ***Screener Profile Report***

This report shows screening outcome rates broken down by screening Facility, Screener and then the stage of screening. Audiologists are included as Screeners on this report as well as other screens in HiTrack. All Most Conclusive Ear results are counted unless filtered by technology type.



This report helps hospital program coordinators answer the following question: "Which screeners at my hospital need additional training and monitoring?"

## **Report Specific Filtering Options**

In addition to the general filtering options this report provides the ability to filter by screening technology.

## ***Screening Results Report***

This report lists babies within the specified Selection Criteria range and infant screening outcomes. You can create a concise listing of babies including only name, medical record number, date of birth and screening results or a more comprehensive report including contact information. The final pages of the report provide summary outcome counts and nursery level counts for babies listed on the report.

## **Report Specific Filtering Options**

This report uses the general filtering options and also allows for filtering by screening outcome classifications.

## **Flow Chart**

The Flow Chart lets screening program coordinators see the 'flow' of babies through their program. It provides a concise one-page summary of the total number of births, the numbers and percentages of babies who passed, were referred, were not screened, the number of children undergoing specific types of diagnostic assessment, and the number found to have confirmed hearing loss. Periodic review of the HiTrack Flow Chart allows coordinators assess the quality of their program.

There are three main vertical sections on this report: Inpatient, Outpatient and Hearing Status (Evaluation). The total # of babies at the top of the chart indicates how many babies at this facility.

The Flow Chart is also clickable. Clicking on a category such as "Screened" will bring up a list of all the babies that are included in the category. These babies can then be added to the Baby Basket for further inspection or revision.

## **Inpatient Section**

The top section shows the number of babies for each IP screening outcome. A special "NICU" category is included in the Not Screening box to account for babies without screening results (Missed or No Information) because they are still in an intensive care unit. Babies with a nursery level other than "Well Baby Nursery" are counted here if they have no IP results.

## **Outpatient Section**

The OP section accounts for all babies who 'Referred' after being screened during IP, whose IP results were 'Inconclusive' or who were 'Not Screened' during IP. The OP section is divided vertically into two sections; OP outcomes for babies who were not screened during IP are shown on the right and OP outcomes for babies whose IP results were Referred or Inconclusive are on the left.



The divided OP section lets program coordinators assess OP screening efficiency. It answers the questions: "How effective are my programs' follow-up policies for babies who come back for additional screening? Are these policies working to screen babies who were missed at birth?"

## **Hearing Status**

The Hearing Status section includes a count for the Total Recommended for Evaluation. This count includes babies who were referred from OP screening, who referred directly to diagnostic evaluation ('No Out. Scrn ->Dx') and babies who received some diagnostic evaluation without being referred in a typical way. These atypical babies are reported in 'Atypical - > Dx' box in the

center of the chart. These are babies whose records indicate evaluation in process or a confirmed hearing status although they passed or were not screened during IP or OP.

Most Conclusive Outcomes for Hearing Disposition are used to categorize babies in the Hearing Status section and must be manually set for each baby.

The Total Recommended for Evaluation is divided into four sub-sections:

Normal Hearing	Babies who have 'Confirmed Normal' Hearing.
Not Evaluated	Babies for whom Tracking has been stopped because they have be marked as Discontinued, Lost/Locate, Refused or Deceased.
Evaluation in process	Babies who do not yet have a confirmed hearing status divided by their current Diagnostic Recommendation. ('In Process'). The last category in this box shows babies who have yet not received a diagnostic recommendation.
Confirmed Loss	Babies with a Confirmed Loss. Unspecified Type of Loss shows the number of babies with only Sound Field Hearing Status or whose type of loss is not yet specified.

HiTrack's Flow Chart also accommodates Single Stage protocols where babies are referred directly for evaluation when they do not pass Inpatient screening. In this case, babies who do not pass Inpatient screening will be counted under the 'No Out. Scrn ->Dx' counts in the Outpatient section chart and then represented in the Hearing Status section below.

### **Confirmed Loss - Other Category**

This count includes babies with a confirmed loss where one or both of the selected Hearing Disposition ear entries are Neural or Undetermined Loss.

### **Evaluation in Process**

This count reports babies who did not Pass Inpatient or Outpatient that have no other diagnostic Recommendations.

### **Unmonitored Infants with Risk Indicators**

This supplemental count reports babies with positive ("Yes") Risk Indicators that don't have an active Risk Monitoring Recommended Action and that are not being evaluated for hearing loss.

### **Pass Infants monitored for Risk Indicators**

This supplemental count reports babies that have an active Risk Monitoring Recommended Action that have Passed Inpatient/Outpatient Screening.

## **Report Specific Filtering Options**

This report uses the general filtering options.

## **Grouping Options**

This report can be grouped in three different ways. Grouping the data as a single report shows all babies for the selection criteria on one Flow Chart. When grouped by Birth facility or IP Screening Site the total list of babies are divided into different Flow Chart pages based on where they were born or where they were screened. Note that changing the grouping options does not change the total number of babies included; it only affects how they are shown.

## ***Hospital Summary***

The Hospital Summary report allows state and regional screening program coordinators see the 'flow' of babies through facilities in their program. It provides concise facility-oriented summaries of the total number of babies, the numbers and percentages of babies who passed, were referred, were not screened, the number of children undergoing specific types of diagnostic assessment, and the number found to have confirmed hearing loss. The numbers on this report are based on the same criteria as the HiTrack Flow Chart, except for that the OP section counts are simplified.

On this report each row of numbers and percentages represents a facility grouped as either a Birthing or Screening location followed by an 'All Facilities' total at the end of the report.



This report allows state and regional coordinators to answer the following question: "How can I compare and report each facility in my state or region in one step?" This report can also be used as a "score card" for facility coordinator training meetings!

### **Report Specific Filtering Options**

This report uses the general filtering options.

### **Grouping Options**

This report can be grouped in three different ways. Grouping the data as a single report shows all babies for the selection criteria as one row with total counts. When grouped by Birth facility or IP Screening Site the total list of babies are divided into different Facility rows based on where they were born or where they were screened. Note that changing the grouping options does not change the total number of babies included; it only affects how they are shown.

## Quarterly Report

This 12-month report is composed of four main categories: Percent Screened as Inpatients, Percent Passing Inpatient Screening, Percent Completing Outpatient Screening and Hearing Status. The first two categories are shown together in the first half of the report and the latter two categories are also shown together in the second half of the report. Each facility is shown as a row and a "Quarter Totals" summary count is included at the end of each section.



This report is helpful because it allows a state or regional coordinator to evaluate screening statistics over time. It helps answer the following questions: "I have just started a new protocol in my state/region how can I tell if overall performance has improved?" and "How can I report percents of screening completion for all facilities in my state/region?"

### Percent Screened as Inpatient

Actual: IP Pass + IP Refer + IP Inconclusive

Expected: Total Babies - IP Refused - IP Deceased

### Percent Passing Inpatient Screening

Actual: IP Pass

Expected: IP Pass + IP Refer + IP Inconclusive (Same as Actual 'Screened')

### No Info

Totals grouped by facility are shown for babies who have no information.

### Deceased

Totals grouped by facility are shown for babies who are marked as Deceased during the IP stage.

### Percent Completing Outpatient Screening

Actual: OP Refer + OP Pass + 'No Scrn ->Dx' + OP Inconclusive

Expected: IP Referred + IP Inconclusive + Missed + IP No-Info - OP Deceased.

This does not include IP Transferred, allowing individual hospitals to possibly show more than 100%.

'No Scrn -> Dx' is included to account for babies with atypical evaluation. Any sites using a [Single Stage](#) protocol will see most of their non-pass babies counted in this category.



### Referred For Diagnostics

Actual: Total Recommended for Evaluation (Includes Atypical Dx),

Expected: Total Babies - IP Deceased - OP Deceased

**Other Dx**

These numbers equal totals for each respective Hearing Status section based on the Flow Chart categories.

**No Follow-Up**

This number includes babies from the following categories: Need to Locate, Refused, Deceased, and Follow-up Discontinued.

**Report Specific Filtering Options**

This report uses the general filtering options and may be grouped by facility.

**Grouping Options**

This report can be grouped in three different ways. Grouping the data as a single report shows all babies for the selection criteria as one row with total counts. When grouped by Birth facility or IP Screening Site the total list of babies are divided into different Facility rows based on where they were born or where they were screened. Note that changing the grouping options does not change the total number of babies included; it only affects how they are shown.

## **Technology Summary**

This report shows summary information for OAE and A-ABR screening technologies. This helps answer the questions: "How can I report screening effectiveness by technology?", "Is my first stage screening technology providing a cost effective waypoint?", and "Do I need to consider changing screening technology at my facility?"



Each count represents an overall "head" screening outcome at a given stage of screening. (For example, if the Most Conclusive Ear results are both OAE then the baby is counted in the OAE column. If the Most Conclusive Ear results are both A-ABR then the baby is counted in the A-ABR column.)

The first column indicates the facility name for which these babies are reported. Inpatient and Outpatient counts are provided for each facility as separate rows of information. The total number of infants screened at the stage is given. The number of infants not screened for that stage is included in curled brackets. "{x}".

The third and fourth columns show baby counts whose head screening outcomes were either OAE or A-ABR.

The last column on the right counts babies who have a Most Conclusive Ear result with one technology for one ear and a different technology for the other ear.

A noteworthy count on this report is the {Non-Pass OAE} count under the A-ABR column. This number shows the number of babies with A-ABR Most Conclusive Results that also did not pass screening with OAE.

### **Report Specific Filtering Options**

This report uses the general filtering options.

### **Grouping Options**

When grouping the report by Birth Facility each row counts represent babies born at each facility. If grouped by IP screening site then the counts represent babies screened at each facility.

## Milestone Summary Report

The Milestone Summary Report (and the Milestone mode on the Baby List) use the following criteria for each milestone. Hearing Loss dates are based on the hearing entry specified as a Most Conclusive outcome on the Hearing Disposition screen. Deceased babies are excluded from these counts.



This report helps answer the questions: "How closely is my program meeting the recommended CDC milestones?" and "What is my programs average age of loss identification and EI enrollment?"

### "Screened by" Milestone:

This count uses the screening Most Conclusive Result date of Inpatient screening or Outpatient screening if it exists. Babies without Most Conclusive Results are considered Late or Overdue for this Milestone.

### "Hearing Identified by" Milestone:

This count reports the user-selected first date of hearing disposition for babies in the following categories:

Inpatient Referred
Inconclusive Inpatient Results
Outpatient Refer
Outpatient Inconclusive Results
Go directly to Evaluation without Outpatient Screening
<i>Not:</i> Outpatient Pass
<i>Not:</i> Follow-Up Discontinued before/during Outpatient Screening

Babies in these categories without a hearing disposition specified are considered Late or Overdue for this Milestone.

### "Enrolled by" Milestone:

This count uses the date enrolled entered from the EI Services screen or the Confirmed Appt. Date/Target date on the "Early Intervention Referral" recommendation for babies with:

Hearing Loss specified
<i>Not:</i> Refused Evaluation
<i>Not:</i> Parental Opt-Out of EI service
<i>Not:</i> Deceased before/during Eval

Babies with a hearing loss specified without EI Services or EI Referrals specified are considered Late or Overdue for this Milestone.



EI referrals are included for users with databases converted from HiTrack 3.5 because Enrollment details were not available in this older version.

### **Report Specific Filtering Options**

This report uses the general filtering options.

### **Grouping Options**

When grouping the report by Birth Facility each page represent births at each facility. If grouped by IP screening site then each page represent babies screened at each facility.

## ***User Defined Reports***

Selecting this report option will open up an "Available User Defined Reports" list. This list lets you, select an existing saved report to run, create a new report or delete an existing report.

Clicking the New button shows the main User Defined report screen that consists of three tabs:

The Column definitions tab allows you to select up to 8 columns to include on the report. Filling in the check box and using the ">" button will add the selected column to the "Selected Elements" list. Use "<" to remove a selected element.

The top of the Selection Criteria tab provides the general report filtering options. Advanced Selection Criteria can be created using the same steps as used in [Advanced Criteria](#) for other reports.

The Sorting tab provides three layers of sorting options. Uncheck the "Ascending" check box to sort the rows in descending (Z to A, Newest to Oldest) order.

Overall report options are provided at the bottom of the screen. You can name the report, save the report, Generate the report or return to the Available User Defined Reports list.

When a report is generated the results of the report, including the selected columns, are shown in a grid that can be printed or exported using buttons below the list.



When report results are saved as a CSV (comma separated values) file HiTrack creates a file that can be opened in Microsoft Excel 2003 automatically for further reporting flexibility.

## Advanced Criteria

Advanced criteria can be added to most reports in HiTrack allowing for more targeted reporting and data quality investigation.

### Adding criteria

The "Add" button below the grid shows the Advanced Selection Criteria screen. The Name field is used to describe complex criteria, such as "Low weight boys". The list of drop down boxes on the left allows you to select and combine multiple criteria elements. After an element is selected the drop down to the right provides relevant options such as "Equals", "Not Equals", "Is Blank" and other options for date and number ranges. The filter criteria for this element becomes more narrow as additional elements are selected on this screen because a Boolean "AND" action is applied.

For example if two elements are selected:

"Gender Equals Male"

"AND..."

"Birth Weight (gram) Less Than # 1200"

A baby must qualify for *both* elements to be included with this filter.

After selecting the desired elements for this criterion the "Add" button can be used to include this in the advanced criteria grid.

### Modifying criteria


Criteria rows may be edited or removed by using the "Edit" and "Delete" buttons shown below the grid. The Edit button will return you to the Advanced Selection Criteria screen to make changes to criteria elements.

### Boolean "OR" Operation

Each criteria row added to the Advanced criteria grid filters with a Boolean "OR" action. For example if you added the criteria from the example above "Low weight boys" and then added a second criteria "Girls in NICU" the following rows would be shown:

Low weight boys	Gender Equals Male AND Birth weight (gram) Less than 1200	<b>OR</b>
Girls in NICU	Gender Equals Female AND Nursery Equals NICU	

With these criteria the possible list of infants has been expanded to include babies with both criteria.



## Appendix C: Classification Combinations

This table lists the overall combined screening outcome used on some reports.

<b>Inpatient</b>	<b>Inpatient Ears</b>	<b>Outpatient</b>	<b>Outpatient Ears</b>	<b><i>CDC Outcome</i></b>
Inconclusive	0	Broken Appt.	0	Unknown
Inconclusive	0	Deceased	0	Deceased
Inconclusive	0	Discontinued	0	Unknown
Inconclusive	0	Inconclusive	0	Unknown
Inconclusive	0	Inconclusive	1	Missed
Inconclusive	0	Inconclusive	2	Not Passed
Inconclusive	0	No Info.	0	Unknown
Inconclusive	0	No Out. Scrn->Eval	0	Not Passed
Inconclusive	0	Passed	1	Missed
Inconclusive	0	Passed	2	Passed
Inconclusive	0	Referred	1	Missed
Inconclusive	0	Referred	2	Referred
Inconclusive	0	Refused	0	Refused
Inconclusive	0	Scheduled	0	Unknown
Inconclusive	1	Broken Appt.	0	Missed
Inconclusive	1	Deceased	0	Deceased
Inconclusive	1	Discontinued	0	Missed
Inconclusive	1	Inconclusive	0	Missed
Inconclusive	1	Inconclusive	1	Missed
Inconclusive	1	Inconclusive	2	Not Passed
Inconclusive	1	No Info.	0	Missed
Inconclusive	1	No Out. Scrn->Eval	0	Not Passed
Inconclusive	1	Passed	1	Not Passed
Inconclusive	1	Passed	2	Passed
Inconclusive	1	Referred	1	Not Passed
Inconclusive	1	Referred	2	Referred
Inconclusive	1	Refused	0	Refused

Inconclusive	1	Scheduled	0	Missed
Inconclusive	2	Broken Appt.	0	Not Passed
Inconclusive	2	Deceased	0	Not Passed
Inconclusive	2	Discontinued	0	Not Passed
Inconclusive	2	Inconclusive	0	Not Passed
Inconclusive	2	Inconclusive	1	Not Passed
Inconclusive	2	Inconclusive	2	Not Passed
Inconclusive	2	No Info.	0	Not Passed
Inconclusive	2	No Out. Scrn->Eval	0	Not Passed
Inconclusive	2	Passed	1	Not Passed
Inconclusive	2	Passed	2	Passed
Inconclusive	2	Referred	1	Not Passed
Inconclusive	2	Referred	2	Referred
Inconclusive	2	Refused	0	Not Passed
Inconclusive	2	Scheduled	0	Not Passed
Missed	0	Broken Appt.	0	Missed
Missed	0	Deceased	0	Deceased
Missed	0	Discontinued	0	Missed
Missed	0	Inconclusive	0	Missed
Missed	0	Inconclusive	1	Missed
Missed	0	Inconclusive	2	Not Passed
Missed	0	No Info.	0	Missed
Missed	0	No Out. Scrn->Eval	0	Not Passed
Missed	0	Passed	1	Missed
Missed	0	Passed	2	Passed
Missed	0	Referred	1	Missed
Missed	0	Referred	2	Referred
Missed	0	Refused	0	Refused
Missed	0	Scheduled	0	Missed
No Info.	0	Broken Appt.	0	Unknown
No Info.	0	Deceased	0	Deceased

No Info.	0	Discontinued	0	Unknown
No Info.	0	Inconclusive	0	Unknown
No Info.	0	Inconclusive	1	Missed
No Info.	0	Inconclusive	2	Not Passed
No Info.	0	No Info.	0	Unknown
No Info.	0	No Out. Scrn->Eval	0	Not Passed
No Info.	0	Passed	1	Missed
No Info.	0	Passed	2	Passed
No Info.	0	Referred	1	Missed
No Info.	0	Referred	2	Referred
No Info.	0	Refused	0	Refused
No Info.	0	Scheduled	0	Unknown
Passed	1	Broken Appt.	0	Missed
Passed	1	Deceased	0	Deceased
Passed	1	Discontinued	0	Missed
Passed	1	Inconclusive	0	Missed
Passed	1	Inconclusive	1	Missed
Passed	1	Inconclusive	2	Not Passed
Passed	1	No Info.	0	Missed
Passed	1	No Out. Scrn->Eval	0	Not Passed
Passed	1	Passed	1	Missed
Passed	1	Passed	2	Passed
Passed	1	Referred	1	Missed
Passed	1	Referred	2	Referred
Passed	1	Refused	0	Refused
Passed	1	Scheduled	0	Missed
Passed	2	Broken Appt.	0	Passed
Passed	2	Deceased	0	Passed
Passed	2	Discontinued	0	Passed
Passed	2	Inconclusive	0	Passed
Passed	2	Inconclusive	1	Not Passed

Passed	2	Inconclusive	2	Not Passed
Passed	2	No Info.	0	Passed
Passed	2	No Out. Scrn->Eval	0	Passed
Passed	2	Passed	1	Passed
Passed	2	Passed	2	Passed
Passed	2	Referred	1	Referred
Passed	2	Referred	2	Referred
Passed	2	Refused	0	Passed
Passed	2	Scheduled	0	Passed
Referred	1	Broken Appt.	0	Missed
Referred	1	Deceased	0	Deceased
Referred	1	Discontinued	0	Missed
Referred	1	Inconclusive	0	Missed
Referred	1	Inconclusive	1	Missed
Referred	1	Inconclusive	2	Not Passed
Referred	1	No Info.	0	Missed
Referred	1	No Out. Scrn->Eval	0	Not Passed
Referred	1	Passed	1	Missed
Referred	1	Passed	2	Passed
Referred	1	Referred	1	Missed
Referred	1	Referred	2	Referred
Referred	1	Refused	0	Refused
Referred	1	Scheduled	0	Missed
Referred	2	Broken Appt.	0	Referred
Referred	2	Deceased	0	Referred
Referred	2	Discontinued	0	Referred
Referred	2	Inconclusive	0	Referred
Referred	2	Inconclusive	1	Not Passed
Referred	2	Inconclusive	2	Not Passed
Referred	2	No Info.	0	Referred
Referred	2	No Out. Scrn->Eval	0	Referred

Referred	2	Passed	1	Referred
Referred	2	Passed	2	Passed
Referred	2	Referred	1	Referred
Referred	2	Referred	2	Referred
Referred	2	Refused	0	Referred
Referred	2	Scheduled	0	Referred
Refused	0	Broken Appt.	0	Refused
Refused	0	Deceased	0	Deceased
Refused	0	Discontinued	0	Refused
Refused	0	Inconclusive	0	Refused
Refused	0	Inconclusive	1	Missed
Refused	0	Inconclusive	2	Not Passed
Refused	0	No Info.	0	Refused
Refused	0	No Out. Scrn->Eval	0	Not Passed
Refused	0	Passed	1	Missed
Refused	0	Passed	2	Passed
Refused	0	Referred	1	Missed
Refused	0	Referred	2	Referred
Refused	0	Refused	0	Refused
Refused	0	Scheduled	0	Refused
Transferred Out	0	Broken Appt.	0	Transferred
Transferred Out	0	Deceased	0	Deceased
Transferred Out	0	Discontinued	0	Transferred
Transferred Out	0	Inconclusive	0	Unknown
Transferred Out	0	Inconclusive	1	Missed
Transferred Out	0	Inconclusive	2	Not Passed
Transferred Out	0	No Info.	0	Transferred
Transferred Out	0	No Out. Scrn->Eval	0	Not Passed
Transferred Out	0	Passed	1	Missed
Transferred Out	0	Passed	2	Passed
Transferred Out	0	Referred	1	Missed

Transferred Out	0	Referred	2	Referred
Transferred Out	0	Refused	0	Refused
Transferred Out	0	Scheduled	0	Transferred
Deceased	0	Broken Appt.	0	Deceased
Deceased	0	Deceased	0	Deceased
Deceased	0	Discontinued	0	Deceased
Deceased	0	Inconclusive	0	Deceased
Deceased	0	Inconclusive	1	Deceased
Deceased	0	Inconclusive	2	Deceased
Deceased	0	No Info.	0	Deceased
Deceased	0	No Out. Scrn->Eval	0	Deceased
Deceased	0	Passed	1	Deceased
Deceased	0	Passed	2	Deceased
Deceased	0	Referred	1	Deceased
Deceased	0	Referred	2	Deceased
Deceased	0	Refused	0	Deceased
Deceased	0	Scheduled	0	Deceased

## Appendix D: Merge Matching Criteria and Overwrite rules

HiTrack checks for existing records when merging from equipment data files or data transfer files. Existing records are found in two steps. First, the database is filtered down to a short list of possible matches. During the second step a matching “score” is determined for each record in this short list. Depending on the score the records will be merged together, added as a new record or left for user resolution.

### Step 1: Filters are used to find a list of possible matches

Based on the type of merge, HiTrack filters down the number of infants to examine for matches by performing the following filters (in order) until matching records are found.

The purpose of the filters is to narrow down a potentially huge database into a small subset that can be checked row by row.

If no records are found after checking all filters the incoming record is added new.

### Equipment Merges

<i>Search for...</i>	<i>Facility</i>
1. Babies born on the same day (excludes time of day) with an exact match on Medical ID	All
2. Babies born within 1 day with matching birth weight	Birth Facility
3. Babies born within 1 day with matching birth weight	Merging Facility

### HiTrack 3.5 ETF Merges

Search for... (across all Facilities)

1. Babies with an exact match on HT35 Unique ID and Medical ID
2. Babies with an exact match on HT35 Unique ID and Birth Date
3. Babies born the same day (excludes time of day) with an exact match on Medical ID
4. Babies with an exact match on HT35 Unique ID and birth weight
5. Babies born within 1 day with matching birth weight and with similar last name (4 letters)

### HiTrack 4 ETF4 Merges

Search for... (across all Facilities)

1. Babies with an exact match on HT4 "ID" (globally unique identifier)
2.. Uses same filters as HiTrack 3.5 ETF Merges

## Step 2: Calculating match score for possible matches

A score is calculated for each baby found with filters and HiTrack decides what to do with the record based on the following ranges of scores:

0-54	Not a match
55 - 69	This is a possible match but needs user resolution. This baby is left in Pending Data as "Needing Attention".
70+	These records will be Automatically merged together. This is a "High Confidence" match.

A score lower than the "unsure match" score will result in a new record after all the babies found by filters have been checked for matches.

Some filters from Step 1 will result in subsets of babies whose records have already been pre-matched on a given field.

### Baby Field Matching Scores

Scores are assigned for matches on the following fields:

HiTrack 4 "ID"	100
HiTrack 3.5 Unique ID	46
Baby DOB	40
Medical ID	25
Birth Weight	10
Baby lastname	5

Score can be inferred as "trust" of the field used in matching. Matches on Globally Unique IDs (GUIDS) from HiTrack 4 ETF4 files ("ID") will always result in an automatic match. (Score 100> Score 70)

### Other Scores

Demographic merges, alldata.txt and demgraph.txt have slightly different scores because they give higher score value (trust) to Medical IDs.

Merges from XML files use the same filters and scores as HiTrack 4 ETF4 rules.

### Match Score Examples

Same lastname, weight and medical ID, all else different (or blank)	Score: $25+10+5=35$	Not a Match
Same DOB and lastname, all else different (or blank)	Score: $40+5=45$	Not a match
Same HT35 Unique ID, all else different (or blank)	Score: 46	Not a match
Same HT35 Unique ID and Medical ID, all else different (or blank)	Score: $46+25=71$	High Confidence Match, Automatically merged
Same DOB, weight and last name, all else different (or blank)	Score: $40+10+5=55$	Unsure match, "Needs Attention".

## Field Overwriting

HiTrack uses rules when baby records are automatically merged from a High Confidence match to determine which fields to overwrite.

HiTrack updates and overwrites existing baby records by splitting each baby record into three areas:

<p>Baby Identifying Fields: Medical ID, Last name, first name and Birth date.</p>	<p>Other Baby demographic fields</p>	<p>Screening History, Diagnostics, Contacts and other historical lists of information for the baby</p>
---	--------------------------------------	--

Depending on the source of data the following rules are used:

<i>Source</i>	Identifying Fields	Other Fields	Historical Lists
Data Transfer from another HiTrack 4 database.	Identifying fields, such as last name and birth date, are handled individually.  Whichever value was changed most recently is retained.	If any of the incoming values for this area have changed more recently then all of the existing fields in this area are overwritten to the incoming values.	All fields for exact* entries are overwritten with incoming values if they have changed more recently.  New records are appended
Data Transfer from a HiTrack 3.5 database	If any of the incoming values have changed more recently then all of the existing fields are overwritten to the incoming values.		All fields for exact* entries are overwritten with incoming values if they have changed more recently.  New records are appended
Equipment Files	No fields are overwritten.  New values for blank fields are inserted.		Existing contacts are updated based on matching last and first names.  All other records are always appended.

### Overwriting Notes

\*Existing historical entries are matched on HiTrack's internal Primary Keys. This allows HiTrack to confidently match existing historical rows.

HiTrack's Change Log keeps logs changes to identifying fields during merge overwrites.

In prior versions of HiTrack Baby "Notes" were considered part of the Baby demographic fields. In HiTrack 4 Notes are maintained as their own separate historical list.

Customized field overwriting is allowed for baby and person matches that require user attention.

HiTrack 4 Data Transfer overwrite rules are used when merging from data in XML sources.

## Appendix E: Most Conclusive Result Calculation

### Overview

HiTrack not only allows users to save multiple screening results for an infant, it also saves data management time by automatically determining which results best reflect the child's overall screening status and need for follow-up. HiTrack uses ranking rules to: 1) select a Most Conclusive set of Left and Right ear results at each screening level (Inpatient/Outpatient), and 2) combine Left and Right ear results to produce a Most Conclusive Head result at each screening level. The Most Conclusive Head results determine how a child will be [counted in reports](#), what follow-up is needed, and which [letter will automatically be generated](#) for the child's parent and health care provider.

### Most Conclusive Result Ranking Rules

HiTrack's ranking rules have been carefully thought through and periodically reviewed to correctly determine the Most Conclusive Head results for children in the most common screening scenarios. In the small number of cases where unique circumstances would suggest that a different decision rule should be applied, a user can override HiTrack's selection and manually designate the Most Conclusive results for an infant.

HiTrack offers two ranking rules:

<i>"Results ranked independent of technology" (default)</i>
---

<i>"Results ranked with precedence given to Automated ABR"</i>
--

At the broadest level both rules use the following algorithm:

<i>MCRs...</i>	<i>Rule</i>
for Individual Ears	Best Test Result is selected. ("Pass beats Refer")
for Head Result	Worst Test Result (from Left and Right MCR) is selected.

## **Most Conclusive Result Window**

HITRCK calculates Most Conclusive Results (MCRs) using a ranking rule for test results but first filters the results by a “window of time” since the latest result was obtained. (Inpatient and Outpatient results are calculated separately.) By default this window is defined as beginning 5 days prior to the newest test and ending on the date of the newest test.

The following two examples illustrate how the Most Conclusive Result Window works.

### **Example #1**

Inpatient Day	1 (DOB)	2	3	5	6	7	8	
	Excluded results		←MCR Window→					MCR Outcome
Right Ear	Pass						Refer	Refer
Left Ear	Pass						Refer	Refer

The shaded area above shows results that will be excluded in the query for seeking Most Conclusive Results, based on this modification to HT V4. In this case, the baby’s result would be calculated as an overall “Refer” instead of an overall “Pass”

### **Example #2**

Inpatient Day	1 (DOB)	2	3	5	6	7	8	
	Excluded results		←MCR Window→					MCR Outcome
Right Ear	Pass				Pass		Refer	Pass
Left Ear	Invalid				Pass		Refer	Pass

In this example the overall results would be a “Pass” because both ear passes are within the ranking filter window.

In general, the MCR Window affects only infants in Inpatient for more than five days with developing hearing loss that results from, for example, Head Trauma, high doses of Ototoxic Medications or Meningitis. As always, clinical judgment should be used in cases where the circumstances and conditions require manual overriding.

## **“Same Day” Most Conclusive Result Option**

HiTrack also allows program coordinators to require that most conclusive results to only be selected from the most recent screenings that happened on the same day as the other ear. When this setting is enabled “Inconclusive” will be reported for any baby that has only one ear pass on the most recent day.

This setting helps programs follow the JCIH recommendation that both ears should be rescreened – even ears that passed in prior attempts.

## Appendix F: Automatic Letter Generation Criteria

### **Criteria**

The following table explains what automatic letters will be generated based on possible combinations of Inpatient and Outpatient Screening results.

<i>Letter Title</i>	<i>Inpatient Head Result</i>	<i>Outpatient Head Result</i>
Diagnostic Referral (PC)	Inconclusive	Inconclusive
Diagnostic Referral (PC)	Inconclusive	No Out. Scrn->Eval
Diagnostic Referral (PC)	Inconclusive	Referred
Diagnostic Referral (PC)	Missed	Inconclusive
Diagnostic Referral (PC)	Missed	No Out. Scrn->Eval
Diagnostic Referral (PC)	Missed	Referred
Diagnostic Referral (PC)	No Info.	Inconclusive
Diagnostic Referral (PC)	No Info.	No Out. Scrn->Eval
Diagnostic Referral (PC)	No Info.	Referred
Diagnostic Referral (PC)	Passed	Inconclusive
Diagnostic Referral (PC)	Passed	No Out. Scrn->Eval
Diagnostic Referral (PC)	Passed	Referred
Diagnostic Referral (PC)	Referred	Inconclusive
Diagnostic Referral (PC)	Referred	No Out. Scrn->Eval
Diagnostic Referral (PC)	Referred	Referred
Diagnostic Referral (PC)	Refused	Inconclusive
Diagnostic Referral (PC)	Refused	No Out. Scrn->Eval
Diagnostic Referral (PC)	Refused	Referred
Diagnostic Referral (PC)	Transferred Out	Inconclusive
Diagnostic Referral (PC)	Transferred Out	No Out. Scrn->Eval
Diagnostic Referral (PC)	Transferred Out	Referred
Diagnostic Referral (Physician)	Inconclusive	Inconclusive
Diagnostic Referral (Physician)	Inconclusive	No Out. Scrn->Eval
Diagnostic Referral (Physician)	Inconclusive	Referred
Diagnostic Referral (Physician)	Missed	Inconclusive

Diagnostic Referral (Physician)	Missed	No Out. Scrn->Eval
Diagnostic Referral (Physician)	Missed	Referred
Diagnostic Referral (Physician)	No Info.	Inconclusive
Diagnostic Referral (Physician)	No Info.	No Out. Scrn->Eval
Diagnostic Referral (Physician)	No Info.	Referred
Diagnostic Referral (Physician)	Passed	Inconclusive
Diagnostic Referral (Physician)	Passed	No Out. Scrn->Eval
Diagnostic Referral (Physician)	Passed	Referred
Diagnostic Referral (Physician)	Referred	Inconclusive
Diagnostic Referral (Physician)	Referred	No Out. Scrn->Eval
Diagnostic Referral (Physician)	Referred	Referred
Diagnostic Referral (Physician)	Refused	Inconclusive
Diagnostic Referral (Physician)	Refused	No Out. Scrn->Eval
Diagnostic Referral (Physician)	Refused	Referred
Diagnostic Referral (Physician)	Transferred Out	Inconclusive
Diagnostic Referral (Physician)	Transferred Out	No Out. Scrn->Eval
Diagnostic Referral (Physician)	Transferred Out	Referred
Inpatient Inconclusive/Missed (PC)	Inconclusive	No Info.
Inpatient Inconclusive/Missed (PC)	Missed	No Info.
Inpatient Inconclusive/Missed (Phys.)	Inconclusive	No Info.
Inpatient Inconclusive/Missed (Phys.)	Missed	No Info.
Inpatient Pass (PC)	Passed	No Info.
Inpatient Pass (Physician)	Passed	No Info.
Inpatient Refer for Outp. (Physician)	Referred	No Info.
Inpatient Refer for Outp. Scn (PC)	Referred	No Info.
Inpatient Refusal (PC)	Refused	No Info.
Inpatient Refusal (Physician)	Refused	No Info.
Outpatient Broken Appt. (PC)	Inconclusive	Broken Appt.
Outpatient Broken Appt. (PC)	Missed	Broken Appt.
Outpatient Broken Appt. (PC)	No Info.	Broken Appt.
Outpatient Broken Appt. (PC)	Passed	Broken Appt.

Outpatient Broken Appt. (PC)	Referred	Broken Appt.
Outpatient Broken Appt. (PC)	Refused	Broken Appt.
Outpatient Broken Appt. (PC)	Transferred Out	Broken Appt.
Outpatient Broken Appt. (Physician)	Inconclusive	Broken Appt.
Outpatient Broken Appt. (Physician)	Missed	Broken Appt.
Outpatient Broken Appt. (Physician)	No Info.	Broken Appt.
Outpatient Broken Appt. (Physician)	Passed	Broken Appt.
Outpatient Broken Appt. (Physician)	Referred	Broken Appt.
Outpatient Broken Appt. (Physician)	Refused	Broken Appt.
Outpatient Broken Appt. (Physician)	Transferred Out	Broken Appt.
Outpatient Pass (PC)	Inconclusive	Passed
Outpatient Pass (PC)	Missed	Passed
Outpatient Pass (PC)	No Info.	Passed
Outpatient Pass (PC)	Passed	Passed
Outpatient Pass (PC)	Referred	Passed
Outpatient Pass (PC)	Refused	Passed
Outpatient Pass (PC)	Transferred Out	Passed
Outpatient Pass (Physician)	Inconclusive	Passed
Outpatient Pass (Physician)	Missed	Passed
Outpatient Pass (Physician)	No Info.	Passed
Outpatient Pass (Physician)	Passed	Passed
Outpatient Pass (Physician)	Referred	Passed
Outpatient Pass (Physician)	Refused	Passed
Outpatient Pass (Physician)	Transferred Out	Passed
Outpatient Refusal (PC)	Inconclusive	Refused
Outpatient Refusal (PC)	Missed	Refused
Outpatient Refusal (PC)	No Info.	Refused
Outpatient Refusal (PC)	Passed	Refused
Outpatient Refusal (PC)	Referred	Refused
Outpatient Refusal (PC)	Refused	Refused
Outpatient Refusal (PC)	Transferred Out	Refused

Outpatient Refusal (Physician)	Inconclusive	Refused
Outpatient Refusal (Physician)	Missed	Refused
Outpatient Refusal (Physician)	No Info.	Refused
Outpatient Refusal (Physician)	Passed	Refused
Outpatient Refusal (Physician)	Referred	Refused
Outpatient Refusal (Physician)	Refused	Refused
Outpatient Refusal (Physician)	Transferred Out	Refused